



SURGICAL
HOME
INSTRUMENTAL Y EQUIPOS MÉDICOS


LAPAROSCOPIA

Instrumental de Acero Quirúrgico



2,8mm&5mm Forceps






Tips	Code	5mm 200mm	5mm 330mm	5mm 440mm	2.8mm 300mm	Handle A	Handle B		
 16mm	S401-251	•	•	•		•		•	
 21mm	S401-253	•	•	•	•	•		•	
 20mm	S401-254	•	•	•		•		•	
 21mm	S401-256	•	•	•	•	•		•	
 14mm	S401-258	•	•	•		•		•	
 15mm	S401-259	•	•	•		•		•	
 23mm	S401-261	•	•	•	•	•		•	
 25mm	S401-262-1	•	•	•	•	•		•	
 23mm	S401-263	•	•	•	•	•		•	
 19mm	S401-262-2	•	•	•	•	•		•	








2.8mm&5mm Forceps



Tips	Code	5mm 200mm	5mm 330mm	5mm 440mm	2.8 mm 300mm	Handle A	Handle B	Handle C	
 25mm	S401-264	Maryland forceps, 70°							
 25mm	S401-265								
 28mm	S401-262-3								
 23mm	S401-382								
 25mm	S401-266								
 18mm	S401-112								
 18mm	S401-113								
 17mm	S401-115								
 25mm	S401-271-1								
 23mm	S401-271-2								









2.8mm&5mm Forceps










Tips	Code	5mm 200mm	5mm 330mm	5mm 440mm	2.8mm 300mm	Handle A	Handle B		
	S401-271-3	Claw forceps, 2x4 greeve teeth							
	S401-272-1	☐	☐	☐	☐		☐		☐
	S401-272-2	☐	☐	☐			☐		☐
	S401-273	☐	☐	☐	☐		☐		☐
	S401-274	☐	☐	☐	☐		☐		☐
	S401-275	☐	☐	☐	☐		☐		☐
	S401-277-1	☐	☐	☐	☐		☐		☐
	S401-278-1	☐	☐	☐	☐		☐		☐
	S401-279-1	☐	☐	☐	☐		☐		☐
	S401-280	☐	☐	☐	☐		☐		☐
	S401-281-1	☐	☐	☐	☐		☐		☐
	S401-330	☐	☐	☐	☐		☐		☐
	S401-332-1	☐	☐	☐	☐		☐		☐



2.8mm&5mm Forceps

Tips	Code	5mm 200mm	5mm 330mm	5mm 440mm	2.8mm 300mm	Handle A	Handle B	Handle C	Handle D
 30mm	S401-332-2	Stomach grasper Single Cution							
 30mm	S401-333	☐	☐	☐	☐	☐	☐	☐	☐
 25mm	S401-334	☐	☐	☐	☐	☐	☐	☐	☐
 21mm	S401-335-2	☐	☐	☐	☐	☐	☐	☐	☐
 22mm	S401-335-1	☐	☐	☐	☐	☐	☐	☐	☐
 35mm	S401-336	☐	☐	☐	☐	☐	☐	☐	☐
 38mm	S401-337	☐	☐	☐	☐	☐	☐	☐	☐
 26mm	S401-338	☐	☐	☐	☐	☐	☐	☐	☐
 37mm	S401-339	☐	☐	☐	☐	☐	☐	☐	☐
 23mm	S401-278-2	☐	☐	☐	☐	☐	☐	☐	☐
 33mm	S401-281-2	☐	☐	☐	☐	☐	☐	☐	☐
 23mm	S401-385	☐	☐	☐	☐	☐	☐	☐	☐
 25mm	S401-386	☐	☐	☐	☐	☐	☐	☐	☐



Tips	Code	5mm 200mm	5mm 330mm	5mm 440mm	5mm 550mm	Handle A	Handle B		
 25mm	S401-387	Knot forceps, B-shaped	☒	☒	☒	☐	☒	☐	☒
 25mm	S401-388	Knot forceps, S-shaped	☒	☒	☒	☐	☒	☐	☒
 26mm	S401-411	Wave grasper, fenestrated	☒	☒	☒	☐	☒	☐	☒
 26mm	S401-397	Alis grasper, horizontal teeth	☒	☒	☒	☐	☒	☐	☒
 28mm	S401-398	Alis grasper, Vertical	☒	☒	☒	☐	☒	☐	☒
 38mm	S401-276	Debaquey grasper	☒	☒	☒	☐	☒	☐	☒
 25mm	S401-390	Cervix grasper	☒	☒	☒	☐	☒	☐	☒
 30mm	S401-392	Atraumatic grasper, with hollow jaws, wavy	☒	☒	☒	☐	☒	☐	☒
 10mm	S401-393	Right angle dissector	☒	☒	☒	☒	☐	☒	☐
 22mm	S401-394	Grasper, to retrieve lost gallstones	☒	☒	☒	☒	☐	☒	☐
 33mm	S401-284-1	Lung grasper	☒	☒	☒	☐	☒	☐	☒
 28mm	S401-395	Grasper Forceps	☒	☒	☒	☐	☒	☐	☒
 28mm	S401-396	Grasper Forceps	☒	☒	☒	☐	☒	☐	☒






10mm Forceps

Tips	Code	10mm 250mm	10mm 330mm	Handle E	Handle F	Handle G	Handle H	Handle I
	S401-368	Mixer						
	S401-369	Maryland forceps, small tips						
	S401-367	Maryland forceps, large tips						
	S401-334	Anvil holder forceps						
	S401-297	Grams plate grasper						
	S401-277-2	Grasper, dolphin nose						
	S401-381	Claw forceps						
	S401-283	Babcock grasper						
	S401-276-2	Broadness mouth grasper						
	S401-276-1	Atraumatic grasper						
	S401-273-1	Alligator grasper						
	S401-399	Suture grasper						
	S401-284-2	Lung grasper						



10mm Forceps

Tips	Code	10mm 250mm	10mm 330mm	Handle E	Handle F	Handle G	Handle H	Handle I
	S401-282	Hemostatic forceps						
		☐	☐		☐			
	S401-105	Biopsy forceps				☐		
		☐	☐			☐		
	S401-334-1	Target type forceps						☐
		☐	☐					☐
	S301-111	Lymph gland forceps			☐			
		☐	☐		☐			
	S301-112	Pulmonary lobe forceps			☐			
		☐	☐		☐			
	S401-376-2	Maryland forceps, articulating					☐	
			☐				☐	
	S302-293	Babcock forceps, articulating					☐	
			☐				☐	
	S401-298	Deba key grasper, articulating					☐	
			☐				☐	
	S401-276-3	Deba key grasper			☐			
		☐	☐		☐			
	S401-279-2	Atraumatic forceps, fenestrated			☐			
		☐	☐		☐			
	S401-383	Scissors, straight		☐				
		☐	☐	☐				
	S401-384	Scissors, curved		☐				
		☐	☐	☐				
	S401-382	Intestine grasper			☐			
		☐	☐		☐			



Trocar, magnetic



Code	3.5mm 95mm	5.5mm 95mm	10.5mm 95mm	12.5mm 95mm	 Magnetic	 With screw	 With safety	 Screw&safety
S401-619.1	•				•			
S401-619.2		•			•			
S401-619.3			•		•			
S401-619.4				•	•			
S401-605.1		•				•		
S401-605.2			•			•		
S401-621.1		•					•	
S401-621.2			•				•	
S401-621.3				•			•	
S401-617.1		•						•
S401-617.2			•					•

Trocar, cross



Code	3.5mm 95mm	5.5mm 95mm	10.5mm 95mm	12.5mm 95mm	 Cross	 With screw	 With safety	 Screw&safety
S401-611.1	•				•			
S401-611.2		•			•			
S401-611.3			•		•			
S401-611.4				•	•			
S401-614.1		•				•		
S401-614.2			•			•		
S401-603.1		•					•	
S401-603.2			•				•	
S401-603.3				•			•	
S401-618.1		•						•
S401-618.2			•					•



Trocar

Trocar, spring



Code	3.5mm 95mm	5.5mm 95mm	10.5mm 95mm	12.5mm 95mm	Code	With screw	With safety	Screw&safety
S401-615.2		•			•			
S401-615.3			•		•			
S401-623.1		•				•		
S401-623.2			•			•		
S401-616.1		•					•	
S401-616.2			•				•	
S401-624.1		•						•
S401-624.2			•					•



Conical Tips, Blunt



Conical Tips, Sharp

TSCS also offer you Trocars with Conical Tips, Blunt and Conical Tips, Sharp as you needed.

9

Reducer



S401-491 Reducer ϕ 10mm- ϕ 5mm



S401-492 Reducer ϕ 10mm- ϕ 5mm



S401-493



ϕ 5.5mm S401-802
 ϕ 10.5mm S401-801



ϕ 5.5mm S401-804
 ϕ 10.5mm S401-803



ϕ 5.5mm S401-806
 ϕ 10.5mm S401-805




+56 9 2015 0890



surgicalhome.cl



Surgical Home Chile



Mode V1 $\phi 5 \times 330 \text{ mm}$



S401-341
Straight

S401-342
Curved

Mode Pistol $\phi 5 \times 330 \text{ mm}$



S401-341
Straight

S401-342
Curved

Mode V3 $\phi 5 \times 330 \text{ mm}$



S401-353
Straight

S401-354
Curved

Mode A $\phi 5 \times 330 \text{ mm}$



S401-341
Straight

S401-342
Curved

S401-345
Automatic resetting

S401-346
Multi-function

Mode V2 $\phi 5 \times 330 \text{ mm}$



S401-343
Straight

S401-344
Curved

S401-347
Automatic resetting


S401-348
Multi-function



Clip Applicator



Hem-o-lok



11

Clip Remover



+56 9 2015 0890



surgicalhome.cl



Surgical Home Chile



Clip Applicator



S401-104 With tyco clips
φ 10x330mm



S401-421
With Johnson & Johnson clips (L)
φ 10x330mm



S401-422
With Johnson & Johnson clips (S)
φ 10x330mm



Bulldog Forceps



S602-121 Bulldog Forceps
φ 10x330mm



S602-120 Bulldog Forceps
S602-120.1
S602-120.2
φ 10x330mm

S602-120.1 Tips

S602-120.2 Tips



Veress Needle

- S401-481 φ 3x120mm
- S401-482 φ 3x140mm
- S401-483 φ 2.5x120mm
- S401-484 φ 2x120mm
- S401-485 φ 2.5x140mm
- S401-486 φ 2x140mm



Spring forceps



S401-389-1
Spring Forceps




S401-389-2
Spring Forceps



S401-389-3
Spring Forceps





Monopolar forceps&scissors


Tips	Code	5mm 350mm	5mm 400mm	Single 360°	Double 360°	Handle A	Handle B		
	S407-256	☐	☐	☐	☐	☐		☐	
	S407-262	☐	☐	☐	☐	☐		☐	
	S407-264	☐	☐	☐	☐	☐		☐	
	S407-281	☐	☐	☐	☐		☐		☐
	S407-280	☐	☐	☐	☐		☐		☐
	S407-275	☐	☐	☐	☐		☐		☐
	S407-265	☐	☐	☐	☐	☐			☐
	S407-271	☐	☐	☐	☐		☐		☐
	S407-398	☐	☐	☐	☐		☐		☐
	S407-336	☐	☐	☐	☐		☐		☐

Bipolar coagulation forceps

Tips	Code	5mm 350mm	5mm 400mm	Single 360°	Double 360°	Handle A	Handle B	Handle C	Handle D
	S407-538	☐	☐	☐	☐	☐		☐	
	S407-536	☐	☐	☐	☐	☐		☐	
	S407-537	☐	☐	☐	☐	☐		☐	
	S407-535	☐	☐	☐	☐	☐		☐	



Articulating Electrodes



Bipolar Coagulation Forceps



14



+56 9 2015 0890




surgicalhome.cl



Surgical Home Chile




Monopolar Electrode



Monopolar cable



Cholangiography Forceps




15

Stone Forceps For bile Duct



Bile Duct Knife



+56 9 2015 0890



surgicalhome.cl



Surgical Home Chile




Hernia Repair Instrument



Golden Finger




Suction&irrigation Tube




Ligature Probe



Single Port Laparoscopy




Suspensory Laparoscopic Instrument



Pediatric Laparoscopy Instrument



Small Incision Surgical Instruments



Fan-Shaped Forcep



Reusable

Disposable shaft with reusable handle
Perfect compatible
Very economical

Disposable shaft

Handle A
Without ratchet

Handle B
With ratchet

Reusable handle



S408-256R
Scissors, double action



S408-262R
Maryland dissector, 30°



S408-264R
Maryland dissector, 70°



S408-281R
Dorsey grasper



S408-280R
Fenestrated grasper



S408-275R
Babcock grasper



S408-265R
Mixer



S408-271.1R
Claw forceps, 2x3 teeth



S408-398R
Allis grasper



S408-336R
Bowel grasper

Articulating Disposable Instruments

Handle
Without Ratchet

Handle A
With Ratchet

Handle B
With Ratchet

Rotating Knob

20

Rotating Knob Color Indication

Tips

Code

5mm
350mm



S408-256F

Scissors, double action



S408-262F

Maryland dissector, 30°



S408-264F

Maryland dissector, 70°



S408-281F

Dorsey grasper



S408-280F

Fenestrated grasper



S408-275F

Babcock grasper



S408-265F

Mixer



S408-271.1F

Claw forceps, 2x3 teeth



S408-398F

Allis grasper



S408-336F

Bowel grasper



+56 9 2015 0890



surgicalhome.cl



Surgical Home Chile



Disposable Instruments




Tips	Code	Rotating Knob Color Indication											
		5mm 350mm											
	S408-256	Scissors, double action	*	*									
	S408-262	Maryland dissector, 30°	*	*									
	S408-264	Maryland dissector, 70°	*		*								
	S408-281	Dorsey grasper	*		*								
	S408-280	Fenestrated grasper	*			*							
	S408-276	Babcock grasper	*				*						
	S408-257	Scissors, high insulation	*					*					
	S408-265	Mixer	*						*				
	S408-271.1	Claw forceps, 2x3 teeth	*							*			
	S408-398	Allis grasper	*								*		
	S408-336	Bowel grasper	*									*	

Disposable Retrieval Bag

TSCS retrieval bag is used in all type of laparoscopic surgery for removing excised tissues, specimen through a small surgical incision, it can be introduced through 10mm or larger size trocar canula.



Ergonomic handle design, push-pull operation comfort, easy in, easy out, easy open, easy close, save time.



TPU high strength medical bags, transparent soft, Ni-Ti open device, fully open, when removed, automatically tighten the bag, to prevent overflow.



Disposable Trocars

TSCS offer you flexibility across surgical specialities, so you can choose the right trocar with the right level of performance when and where it's needed.

Unique puncture needle design, less trauma, faster healing

Scalable valve with Freee reducer.

Casing in accordance with serrated ergonomic design, both fixation device to prevent slippage, and reduce damage to the skin.

Direct visualisation eliminates blind entry by enabling visualisation of tissue layers during insertion and multiple needle tips for choice.

S408-207-1 ϕ 12mm

S408-204-1 ϕ 10mm

S408-205-1 ϕ 5mm

Disposable suction&irrigation tube

S408-501-1 With Hook

S408-501-2 Without Hook

Disposable veress needle

S408-210 ϕ 2x120mm
S408-211 ϕ 2x140mm

S408-212

S408-213

Disposable Retractor Rings

S408-601 Retractor Rings

S408-601-2 Hook, sharp ϕ 5mm


S408-601-3 Hook, blunt ϕ 10mm



Disposable Bipolar Coagulation Forceps




Instruments For Thoracoscopic Surgery



23

Sternum retractor



+56 9 2015 0890



surgicalhome.cl



Surgical Home Chile



Instruments for Thoracoscopic Surgery



S301-S-1
Scissors Curved
340mm



S301-S-7
Needle Holder
340mm



S301-S-5
Dissector Forceps Straight
340mm



S301-S-11-1
Tissue Forceps (Oval type)
340mm



S301-S-11-2
Tissue Forceps (Without teeth)
340mm



S301-S-3
Lung Forceps
340mm



S301-S-3-1
Lung Forceps (L)
340mm



S301-S-3-2
Lung Forceps (M)
340mm



S301-S-3-2
Lung Forceps (S)
340mm



S301-S-4
Dissector Forceps Curved
340mm



S301-S-6
Grasper
340mm



S301-S-10
Dissector Forceps Right Angled
340mm




S301-S-9
Lymph Gland Forceps
340mm




S301-S-8
Tissue Forceps
340mm



S301-S-2
Tissue Forceps (Without teeth)
340mm



T140 1902
Three Joint
340mm



Single Joint



S301-D-1
Scissors Curved
340mm



S301-D-3
Lung Forceps
340mm



S301-D-3-1
Lung Forceps (L)
340mm



S301-D-3-2
Lung Forceps (M)
340mm



S301-D-6
Grasper
340mm



S301-D-7
Needle Holder
340mm



S301-D-3-3
Lung Forceps (S)
340mm



S301-D-4
Dissector Forceps Curved
340mm



S301-D-5
Dissector Forceps Straight
340mm



S301-D-11-1
Tissue Forceps (Oval type)
340mm



S301-D-11-21
Tissue Forceps (Without teeth)
340mm



S301-D-8
Tissue Forceps
340mm



S301-D-9
Lymph Gland Forceps
340mm



S301-D-10
Dissector Forceps Right Angled
340mm



Multifunctional Uterus Apparatus



Combined Multifunctional Uterus Apparatus



NEW Uterus Apparatus




Uterus Apparatus



Uterus Apparatus A




Uterus Apparatus B





Gynecological Forceps




Forceps




Special Vaginal Retractor(Suture Needle)




Instruments for Hysterectomy surgery



Hysteroscopy(Inspection)




Hysteroscopy



Hysteroresectoscope




Hysteroscopy(work)



29

Hysteroscopy



+56 9 2015 0890



surgicalhome.cl



Surgical Home Chile




Morcellator



S504-125
Control Unit

Motor Cable



S504-106
Handwheel

S401-258
Scissors
φ5x330mm

S501-109
Uterus Grasping Forceps
φ10x420mm

S504-110
Uterus Grasping Forceps
φ5x420mm

30



S504-101 Hysterectomy Apparatus φ15x330mm

S504-103 Hysterectomy Apparatus φ18x330mm

S504-102 Hysterectomy Apparatus φ15x330mm

S504-104 Hysterectomy Apparatus φ18x330mm

S504-121 Aliming Stick

S504-105 Uterus Measuring Stick φ5x500mm

S504-119 Knot-Pushing Stick

S504-107 Large Bodkin φ10x420mm

S401-461 Irrigating Tube



S504-112 Dilator φ10-15mm

S504-113 Dilator φ10-18mm

S504-114 Reducer φ18-5mm

S504-115 Reducer φ15-5mm

S504-116 Puncture Cannula φ15x95mm

S504-117 Puncture Cannula φ18x95mm

S504-123 Cervix Forceps

S504-122 Cervix Forceps

S504-118A
Electric Motor



+56 9 2015 0890



surgicalhome.cl



Surgical Home Chile



Morcellator(continuous resection)






S504-101C
Big Uterus Grasping Forceps
φ10x375mm

S504-102C
Small Uterus Grasping Forceps
φ5x375mm




Optional instrument






Pedestal Laryngoscope



Otoscope Instrument



Bronchoscopy

S207-101 Bronchoscope ϕ 7x280mm
 S207-102 Bronchoscope ϕ 6x280mm
 S207-103 Bronchoscope ϕ 5x260mm
 S207-104 Bronchoscope ϕ 4x250mm
 S207-105 Bronchoscope ϕ 3.5x250mm

S207-110 ϕ 2x360mm
 Foreign body forceps, for trachea, peanut-shaped

S207-111 ϕ 2x360mm
 Foreign body forceps, for trachea, alligator

S207-112 ϕ 2x360mm
 Foreign body forceps, for trachea, wheat-shaped

S207-113 Suction Tube ϕ 2.5x320mm
 S207-114 Suction Tube ϕ 3x320mm
 S207-115 Suction Tube ϕ 4x320mm

Shaver (For ENT)




S210-101 Handpiece



S210-104 Serration blade

S210-106 Spherical Blade



S210-105 Smooth Blade

S210-107 Serration Blade

S210-108 Blade, Curved



S210-103 Foot Switch



S210-102 Control Unit

Lithotripsy

S603-102 Strong-force Lithoclast, Curved



Obturator



S603-106
 Evacuator (Elric)

S603-107
 Suction Connector (With Tube)



S602-116-1
 In-flow Valve



S602-116-2
 Out-flow Valve

S602-103
 Endoscope 70° ϕ 4x302mm

Lithotripsy

S602-102 Strong-force Lithoclast, straight



S606-103 Sheath



Obturator



S602-107 Suction Connector (With Tube)




S603-116 Evacuator (Elric)


S602-104
 Endoscope 30°
 ϕ 4x302mm




Pelvis Bladder Endoscope



Resectoscopy(worked with Green Laser)



Cystoscope



Urethrotomy



Resectoscope




Electrosurgical Unit



S601-102
Electrosurgical Unit

- ★ 400W
- ★ Working Modes and output in every Mode
Cut 5–400W
Falgurate 5–150W
- ★ Bipolar 7–70W
- ★ Voltage: 100–220V 50/60Hz
- ★ Used For Laparoscopy any Resectoscopy

Ureterorenoscope



S604-101 Ureterorenoscope Fr7.5

Thyroid Instruments



S208-103 Dissecting Forceps ϕ 3x330mm




S208-101 Trocar ϕ 5x95mm
S208-102 Trocar ϕ 3x95mm



S208-105 Needle Holder Forceps ϕ 5x330mm



S208-106 Suction & Irrigation tube ϕ 3x330mm



S208-108 Dissecting ϕ 2mm

S208-112 Injection Needle ϕ 2.8mm

S208-109 Dissecting ϕ 6mm

S208-110 Thyroid Dissecting ϕ 10mm




S208-111 Hook



Arthroscope Instrument



S205-101 Endoscope 0° .4mm
S205-102 Endoscope 30° .4mm
S205-103 Endoscope 70° .4mm



Shaver For Arthroscopy



Diskoscope



S802-102 Sheath



S802-101
Diskoscope Endoscope 30°
Ø 4mm



S802-112
Bayoneted Forceps
90° 3mm



S802-113
Bayoneted Forceps
90° 5mm



S802-114
Bayoneted Forceps
40° 3mm



S802-115
Bayoneted Forceps
40° 5mm



S802-103 Dilator



S802-104 Dissector



S802-105
Short Ball Probe



S802-106
Guide Pen upward length 2.5mm
Guide Pen upward length 3.5mm



S802-107
Long Ball Probe



S802-108
Dissector

Sigmoidoscope



S701-106
Inflating Ball



S701-102
Observing Endoscope



S701-103 Operators



S602-102
Endoscope 30°



S701-104 Obturator



S602-110
Flexible Biopsy Forceps




S602-111
Flexible Grasper



S602-112
Flexible Foreign Body Forceps



S602-113
Flexible Scissors



S701-110 Suction



S701-107 Cleaning Brush



S701-105 Cotton Forceps



Esophagoscope

S206-101 Esophagoscope $\phi 9.5 \times 13 \times 180 \text{ mm}$
 S206-102 Esophagoscope $\phi 10 \times 14 \times 250 \text{ mm}$
 S206-103 Esophagoscope $\phi 11 \times 15 \times 300 \text{ mm}$
 S206-104 Esophagoscope $\phi 11 \times 15 \times 350 \text{ mm}$
 S206-105 Esophagoscope $\phi 12 \times 16 \times 400 \text{ mm}$
 S206-106 Esophagoscope $\phi 13 \times 17 \times 450 \text{ mm}$



S206-109 Suction Tube $\phi 2.5 \times 450 \text{ mm}$

S206-110 Suction Tube $\phi 3 \times 450 \text{ mm}$

S206-111 Suction Tube $\phi 4 \times 450 \text{ mm}$

S206-112 Handle



S206-107
 $\phi 4 \times 500 \text{ mm}$

Foreign Body Forceps, Duck-Mouth



S206-108
 $\phi 4 \times 500 \text{ mm}$

Foreign Body Forceps, Backward teeth



S206-109
 $\phi 3.5 \times 400 \text{ mm}$

Polypus forceps, Brace-shaped

Beauty Surgery Instruments



S105-001 Endoscope Sheath

S105-008 Mercedes Casing $\phi 5 \times 300 \text{ mm}$
 S105-009 Mercedes Casing $\phi 4 \times 300 \text{ mm}$
 S105-010 Mercedes Casing $\phi 3 \times 250 \text{ mm}$
 S105-011 Mercedes Casing $\phi 3 \times 150 \text{ mm}$
 S105-012 Mercedes Casing $\phi 2 \times 150 \text{ mm}$
 S105-013 Mercedes Casing $\phi 2 \times 100 \text{ mm}$

S105-002 Mercedes $\phi 5 \times 300 \text{ mm}$
 S105-003 Mercedes $\phi 4 \times 300 \text{ mm}$
 S105-004 Mercedes $\phi 3 \times 250 \text{ mm}$
 S105-005 Mercedes $\phi 3 \times 150 \text{ mm}$
 S105-006 Mercedes $\phi 2 \times 150 \text{ mm}$
 S105-007 Mercedes $\phi 2 \times 100 \text{ mm}$



S105-016 Stripper



S105-014
 Stripper Down



S105-015
 Stripper Up



S105-016
 Stripper Flat



S105-017
 Stripper Rhombus



S105-018
 Stripper Arc




S105-019
 Stripper



S105-020
 Stripper



S105-021
 Stripper



S105-022
 Stripper



S105-023
 Stripper



S105-024
 Scissors (Straight)



S105-025
 Scissors (Curved)



S105-026
 Scissors (Hooked)



S105-027
 Scissors (Straight)



S105-028
 Scissors (Curved)



40



+56 9 2015 0890



surgicalhome.cl



Surgical Home Chile





Otoscope

ANGLE	CODE	SPECIFICATION	
	S216-107 S216-101 S212-136	$\varnothing$ 2.7x60mm 0° $\varnothing$ 2.7x95mm 0° $\varnothing$ 4x50mm/60mm 0°	
	S216-102 S212-137	$\varnothing$ 2.7x95mm 30° $\varnothing$ 4x50mm/60mm 30°	
	S212-138	$\varnothing$ 4x50mm/60mm 70°	

Sinuscope


ANGLE	CODE	SPECIFICATION	
	S212-101 S212-104	$\varnothing$ 2.7x175mm 0° $\varnothing$ 4x175mm 0° (110° Wide Angle)	
	S203-102 S205-105	$\varnothing$ 2.7x175mm 30° $\varnothing$ 4x175mm 30° (110° Wide Angle)	
	S212-107	$\varnothing$ 4x175mm 45°	
	S212-108	$\varnothing$ 4x175mm 70°	

Laryngoscope

ANGLE	CODE	SPECIFICATION	Laryngoscope
	S212-222	$\varnothing$ 8x185mm 70°	
	S212-221	$\varnothing$ 8x185mm 90°	
	S212-220	$\varnothing$ 5x175mm 70°	
	S212-219	$\varnothing$ 4x185mm 12°	



Cystoscope




ANGLE	CODE	SPECIFICATION
0°	S603-101	Φ 2.7x302mm 0°
	S602-105	Φ 3x302mm 0°
	S602-101	Φ 4x302mm 0°
12°	S602-104	Φ 4x302mm 12°
30°	S603-102	Φ 2.7x302mm 30°
	S602-106	Φ 3x302mm 30°
	S602-102	Φ 4x302mm 30°
45°	S602-107	Φ 4x302mm 45°
70°	S602-103	Φ 4x302mm 70°

Cystoscope



Laparoscope




ANGLE	CODE	SPECIFICATION
0°	S401-011	Φ 5x3.05 mm 0°
	S401-001	Φ 10x3.30mm 0° (110° Wide Angle)
30°	S401-012	Φ 5x3.05 mm 30°
	S401-002	Φ 10x3.30mm 30° (110° Wide Angle)
0°	S401-003	Φ 10x2.90mm 0°

Laparoscope

Laparoscope



Hysteroscope




ANGLE	CODE	SPECIFICATION
0°	S501-100	Φ 5.6x202mm 0°
	S501-111	Φ 6.5x197mm 0°
22°	S501-099	Φ 5.6x182mm 22°
	S501-101	Φ 6.8x188mm 22°
30°	S501-098	Φ 4x302 mm 30°

Hysteroscope

Hysteroscope




Arthroscope



ANGLE	CODE	SPECIFICATION
0°	S801-101	Φ 4x175mm 0°
30°	S801-102	Φ 4x175mm 30°
70°	S801x103	Φ 4x175mm 70°

Arthroscope





V230 V-3

HD V230-HD Waterproof camera

- ▶ Resolution:1920x1080
- ▶ Display:9 inch LCD screen
- ▶ TV System:PAL&NTSC
- ▶ SNR:54db
- ▶ Video output:S-video,RGB,HDMI
- ▶ White Balance:AWC,Manu,Auto Function:Auto Freeze,Zoom,record



S904-100

LED light source

- ▶ 100W LED
- ▶ Direct current fan Quiet operation
- ▶ Lamp life ≥ 50,000 hour
- ▶ Low temperature no infrared and ultraviolet rays
- ▶ 6500k color temp. approx. = Natural sunlight
- ▶ No high voltage components very Safe and efficient environment



S905-106

Suction irrigation pump

- ▶ Fully automatic, easy intuitive pressure suction and irrigation can be adjusted according to surgery setting
- ▶ Maximum liquid volume ≥ 5000ml/min
- ▶ Low noise with superior performance, special filters ensure optimally sterile conditions for minimally invasive surgery
- ▶ Range of pos-pressure:10-50kPa
- ▶ Range of neg-pressure:-50-10kPa



S401-480

Inflator

- ▶ Fully automatic, quiet operation
- ▶ range of flow:10-30L/3 level flow rate
- ▶ Automatic safety alarm function
- ▶ Gas supply design for great safety
- ▶ Dual gas filtering and decompression design
- ▶ Range of press:5-30mmHg




S905-103

Auto pump

- ▶ used for Arthroscopy, Hysteroscopy and Urology systems
- ▶ Sterile tubing sets with sensor, a peristaltic pump, with two channels
- ▶ Emulsion tube and shim are replaceable(5pcs, L: approx.16cm)
- ▶ High sensitivity pressure sensor
- ▶ Range of press:10-400mmHg(2-53.3kPa)
- ▶ Range of flow:0.1-1.0 L/min



Portable LED Light Source



Value



Accessories





SURGICAL
HOME
INSTRUMENTAL Y EQUIPOS MÉDICOS

CONTÁCTANOS



+56 9 2015 0890



[surgicalhome.cl](https://www.instagram.com/surgicalhome.cl)



Surgical Home Chile