

**SURGINEST
GLOBAL**
INSTRUMENTAL Y EQUIPOS MÉDICOS





AS101 (Arthroscope 4mm/0°)

A1011 Arthroscope sheath
ø5.5mm×140mm×0°

A1012 Blunt obturator
ø4.3mm×175mm

A1013 Puncture needle
ø4.3mm×175mm



4175.00 ø4mm×175mm×0°

4175.00W1 ø4mm×175mm×0° (90° wide angle)

4175.00W2 ø4mm×175mm×0° (110° wide angle)

AS102 (Arthroscope 4mm/30°)

A1021 Arthroscope sheath
ø5.5mm×140mm×30°

A1022 Blunt obturator
ø4.3mm×175mm

A1023 Puncture needle
ø4.3mm×175mm



4175.30 ø4mm×175mm×30°

4175.30W1 ø4mm×175mm×30° (90° wide angle)

4175.30W2 ø4mm×175mm×30° (110° wide angle)

AS103 (Arthroscope 4mm/70°)

A1031 Arthroscope sheath
ø5.5mm×140mm×70°

A1032 Puncture needle
ø4.3mm×175mm

A1033 Blunt obturator
ø4.3mm×175mm



4175.70 ø4mm×175mm×70°

4175.70W1 ø4mm×175mm×70° (90° wide angle)


4175.70W2 ø4mm×175mm×70° (110° wide angle)

Arthroscope set

AS201/AS202/AS203 (click locked)




AS201 (Arthroscope 4mm/0°)




4175.00 ø4mm×175mm×0°
4175.00W1 ø4mm×175mm×0° (90° wide angle)
4175.00W2 ø4mm×175mm×0° (110° wide angle)

AS202 (Arthroscope 4mm/30°)



4175.30 ø4mm×175mm×30°
4175.30W1 ø4mm×175mm×30° (90° wide angle)
4175.30W2 ø4mm×175mm×30° (110° wide angle)

AS203 (Arthroscope 4mm/70°)




4175.70 ø4mm×175mm×70°
4175.70W1 ø4mm×175mm×70° (90° wide angle)
4175.70W2 ø4mm×175mm×70° (110° wide angle)

Arthroscope set

AS301/AS302/AS303 (click locked)




AS301 (Arthroscope 4mm/0°)




3175.00 ø3mm×175mm×0°

AS302 (Arthroscope 4mm/30°)



3175.30 ø3mm×175mm×30°

AS203 (Arthroscope 4mm/70°)



3175.70 ø3mm×175mm×70°


Arthroscope set

AS401/AS402/AS403 (click locked)



AS401 (Arthroscope 2.7mm/0°)

- A4011 Arthroscope sheath
ø4mm×78mm×0°
- A4012 Blunt obturator
ø3mm×110mm
- A4013 Puncture needle
ø3mm×110mm



7110.00

ø2.7mm×110mm×0°

AS402 (Arthroscope 2.7mm/30°)

- A4021 Arthroscope sheath
ø4mm×78mm×30°
- A4022 Blunt obturator
ø3mm×110mm
- A4023 Puncture needle
ø3mm×110mm



7110.30

ø2.7mm×110mm×30°

AS403 (Arthroscope 2.7mm/70°)

- A4031 Arthroscope sheath
ø4mm×78mm×70°
- A4032 Puncture needle
ø3mm×110mm
- A4033 Blunt obturator
ø3mm×110mm




7110.70

ø2.7mm×110mm×70°

Arthroscope set


AS501/AS502/AS503 (rotating locked, Arthrex type)



AS501 (Arthroscope 4mm/0°)

A5011 Arthroscope sheath
ø6mm×155mm×0°

A5012 Blunt obturator
ø4mm×160mm




4155.00 ø4mm×155mm×0°

AS502 (Arthroscope 4mm/30°)

A5021 Arthroscope sheath
ø6mm×155mm×30°

A5022 Blunt obturator
ø4mm×160mm




4155.30 ø4mm×155mm×30°

AS503 (Arthroscope 4mm/70°)

A5031 Arthroscope sheath
ø6mm×155mm×70°

A5032 Blunt obturator
ø4mm×160mm




4155.70 ø4mm×155mm×70°

Arthroscope set


AS601/AS602/AS603




AS601 (Arthroscope 2.7mm/30°)



AS602 (Arthroscope 1.9mm/30°)



AS603 (Arthroscope 4mm/30°)



Arthroscope instruments



A3001 Trocar cannula
ø5.8mm×74mm



A3002 Trocar (needle head)
ø5mm×93mm



A3003 Trocar (blunt head)
ø5mm×93mm



A3004 Foreign body forceps
ø4.5mm×135mm



A3005 Meniscus forceps (straight)
ø4.5mm×135mm×0°



A3006 Meniscus forceps (left)
ø4.5mm×135mm×45°



A3007 Meniscus forceps (right)
ø4.5mm×135mm×45°



A3008 Synovium forceps (straight)
ø3.4mm×140mm



A3009 Needle location biopsy forceps
ø3.4mm×145mm



A3010 Hook scissors (straight)
ø3.4mm×140mm



A3011 Curved scissors (left)
ø3.4mm×140mm



A3012 Curved scissors (right)
ø3.4mm×140mm



A3013 Punch forceps
ø5mm×150mm



Arthroscope instruments



A4001 Hock knife
ø3mm×105mm



A4002 Backward knife
ø3mm×105mm



A4003 Meniscus knife
ø3mm×105mm



A4004 Rosette knife
ø3mm×105mm



A4005 Half-round knife
ø3mm×105mm



A4006 Spade knife (straight)
ø3mm×105mm



A4007 Curette
ø3mm×105mm



A4008 Curette
ø3mm×105mm



A4009 Bone file
ø3mm×105mm



A4010 Banana knife
ø3mm×105mm



A4011 Right angle probe (3mm)

A4012 Right angle probe (5mm)



A4013 Punch forceps (right)
ø4.5mm×175mm×45°



A4014 Punch forceps (left)
ø4.5mm×175mm×45°



A4015 Punch forceps (right)
ø4.5mm×175mm×90°



A4016 Punch forceps (left)
ø4.5mm×175mm×90°



A4017 Joint pliers forceps
ø4.5mm×175mm×0°



Arthroscope Instruments



1. Ergonomic handles & pinless joint boost power/control
2. Precise force distribution for uniform tissue handling
3. 15+ punch/jaw options adapt to diverse procedures
4. Medical-grade steel/polymer ensure extended longevity
5. Optimized gaps + Luer-lock port simplify precision/cleaning

Comparison Chart



1 Joint



2 Knee



3 Shoulder



4 Hip

Working length

90mm

130mm

160mm

185mm

Handle A, Classic



Handle B, Ergo



Handle C, Positioning



Handle D, Slide-lock



Arthroscope Instruments



Straight



Upbiter



Upswept

Bite width



Tip width



Curved left



Curved right



Joint Instruments Set (handle A, D)

J1001 Basket Punch
(straight)
ø1.2x90mm



J1002 Basket Punch
(left)
ø1.2x90mm



J1003 Basket Punch
(right)
ø1.2x90mm



J1004 Basket Punch
(straight)
ø1.7x90mm



J1005 Basket Punch
(left)
ø1.7x90mm



J1006 Basket Punch
(right)
ø1.7x90mm



J1007 Basket Punch
(upward)
ø1.7x90mm



J1008 Basket Punch
(straight)
ø2.2x90mm



J1009 Basket Punch
(left)
ø2.2x90mm



J1010 Basket Punch
(right)
ø2.2x90mm



J1011 Basket Punch
(upward)
ø2.2x90mm



J1012 Tissue Grasper
(allis)
ø3.0x90mm



J1013 Tissue Grasper
(loose)
ø3.0x90mm



J1014 Tissue Grasper
(biopsy)
ø3.0x90mm



J1015 Tissue Grasper
(saw-tooth)
ø3.0x90mm



Arthroscope Instruments

Joint Instruments Set



J2001 Probe
ø3x120mm



J2002 Hook Knife
ø3x120mm



J2003 Knock Bar
ø5x60mm



J2004 Curette
ø3x160mm



J2005 Chondro Pick
ø4x135mm, 15°



J2006 Chondro Pick
ø4x135mm, 45°



J2007 Chondro Pick
ø4x135mm, 60°



J2008 Curette(Annular)




J2009 Trocar Sheath
ø5.2x80mm




J2010 Trocar Sheath
ø4.3x95mm



J2011 Knife Handle



J2012 Knife, 15°
(compatible with J1011)



J2013 Knife, 45°
(compatible with J1011)



J2014 Knife, 90°
(compatible with J1011)



J2015 Elevator
(blunt head)
ø3x120mm



J2016 Elevator
(pointed head)
ø3x120mm

Arthroscope Instruments

Knee Instruments Set (handle A, B, D)

K1001 Basket Punch
(straight, round)
ø4.4x120mm



K1002 Basket Punch
(upward, round)
ø4.4x120mm



K1003 Basket Punch
(right, round)
ø4.4x120mm



K1004 Basket Punch
(left, round)
ø4.4x120mm



K1005 Basket Punch
(straight)
ø3.2x120mm



K1006 Basket Punch
(upward)
ø3.2x120mm



K1007 Basket Punch
(rod upward)
ø3.2x120mm



K1008 Basket Punch
(right)
ø3.2x120mm



K1009 Basket Punch
(left)
ø3.2x120mm



K1010 Basket Punch
(right, upward)
ø3.2x120mm



K1011 Basket Punch
(left, upward)
ø3.2x120mm



K1012 Basket Punch
(right, 90°)
ø3.4x120mm



K1013 Basket Punch
(left, 90°)
ø3.4x120mm



K1014 Basket Punch
(left backbiter)
ø2.4x120mm



K1015 Basket Punch
(right backbiter)
ø2.4x120mm



K1016 Basket Punch
(backbiter)
ø2.6x120mm



K1017 Basket Punch
(straight)
ø2.5x120mm



K1018 Basket Punch
(upward)
ø2.5x120mm



K1019 Basket Punch
(straight)
ø2.2x120mm



K1020 Basket Punch
(upward)
ø2.2x120mm



K1021 Basket Punch
(straight)
ø1.8x120mm



K1022 Basket Punch
(upward)
ø1.8x120mm



Arthroscope Instruments

Knee Instruments Set (handle A, B, D)

K1023 Basket Punch (right) ø1.8x120mm	
K1024 Basket Punch (left) ø1.8x120mm	
K1025 Basket Punch (straight) ø4.5x120mm	
K1026 Basket Punch (upward) ø4.5x120mm	
K1027 Basket Punch (right) ø3.2x120mm	
K1028 Basket Punch (left) ø3.2x120mm	
K1029 Tissue Grasper (allis) ø4.0x120mm	
K1030 Tissue Grasper (allis) ø4.0x160mm	
K1031 Tissue Grasper (saw-tooth) ø4.0x120mm	
K1032 Tissue Grasper (saw-tooth) ø4.0x160mm	
K1033 Tissue Grasper (loose) ø4.0x120mm	
K1034 Tissue Grasper (loose) ø4.0x160mm	
K1035 Tissue Grasper (canie-tooth) ø4.0x120mm	
K1036 Tissue Grasper (canie-tooth) ø4.0x160mm	
K1037 Suture Grasper (straight) ø4.0x120mm	
K1038 Suture Grasper (straight) ø4.0x160mm	

K1039 Suture Grasper (right) ø4.0x120mm	
K1040 Suture Grasper (right) ø4.0x160mm	
K1041 Suture Grasper (left) ø4.0x120mm	
K1042 Suture Grasper (left) ø4.0x160mm	
K1043 Suture Grasper (round) ø4.0x120mm	
K1044 Suture Grasper (round) ø4.0x160mm	
K1045 Suture Grasper (tip tooth) ø4.0x160mm	
K1046 Suture Grasper (lock, close) ø3.2x160mm	
K1047 Suture Grasper (lock, open) ø3.2x160mm	
K1048 Tissue Scissors (straight) ø4.0x120mm	
K1049 Tissue Scissors (straight) ø4.0x160mm	
K1050 Tissue Scissors (right) ø4.0x120mm	
K1051 Tissue Scissors (right) ø4.0x160mm	
K1052 Tissue Scissors (left) ø4.0x120mm	
K1053 Tissue Scissors (left) ø4.0x160mm	

Arthroscope Instruments

Knee Instruments Set



K2001 Suture Cutter
ø4x160mm



K2002 Suture Passer
ø4x140mm



K2003 Suture Passer Needle
ø1.5x210mm
(matched with K2002)

Nucleus Pulposus Forceps



- R10160 ø1x160mm, 0deg
- R20160 ø2x160mm, 0deg
- R30160 ø3x160mm, 0deg
- R40160 ø4x160mm, 0deg
- R50160 ø5x160mm, 0deg



- R10164 ø1x160mm, 45deg
- R20164 ø2x160mm, 45deg
- R30164 ø3x160mm, 45deg
- R40164 ø4x160mm, 45deg
- R50164 ø5x160mm, 45deg



K2011 Puncture Needle
ø0.8x105mm



K2012 Puncture Needle
ø1.3x105mm



K2013 Puncture Needle (curved)
ø1.3x105mm



K2014 Suture Hook
ø2.5x100mm



K2015 Suture Hook
ø2.5x100mm



K2016 Suture Hook
ø2.5x100mm



K2017 Suture Hook
ø2.5x100mm



K2018 Suture Hook
ø2.5x100mm



K2019 Suture Hook
ø2.5x100mm



K2020 Guide Needle
ø0.8x200mm

Arthroscope Instruments

Knee Instruments Set



K2021 Trocar Sheath
ø5.8x105mm



K2022 Trocar (sharp)
ø4x110mm



K2023 Trocar (blunt)
ø4x110mm



K2024 Hook Knife
ø3x120mm



K2025 Banana Shaped Knife
ø3x120mm



K2026 Semicircle Shaped Knife
ø3x120mm



K2027 Rose Shaped Knife
ø3x120mm



K2028 Cutting Knife
ø3x120mm



K2029 Contra Angled Knife
ø3x120mm

Meniscus Heel Instruments Set



M1001 Drill Guide Arm
M1002 Drill Guide Sleeve

M1003 Pin Aimer (left)



M1004 Pin Aimer (right)



M1005 Measuring Cylinder



M2001 Guide Wire
ø2.8x200mm



M2002 Guide Wire
ø2x220mm



M3001 Curette
ø5x120mm



Arthroscope Instruments

Cruciate Ligament Reconstruction Instruments Set



Femoral Cannulated Drill Bit

Model number	Diameter	Working length	Picture
F400190	4.0mm	190mm	
F400190	4.5mm	190mm	
F400190	5.0mm	190mm	
F400190	5.5mm	190mm	
F400190	6.0mm	190mm	
F400190	6.5mm	190mm	
F400190	7.0mm	190mm	
F400190	7.5mm	190mm	
F400190	8.0mm	190mm	
F400190	8.5mm	190mm	
F400190	9.0mm	190mm	
F400190	9.5mm	190mm	
F400190	10.0mm	190mm	
F400190	10.5mm	190mm	
F400190	11.0mm	190mm	
F400190	11.5mm	190mm	
F400190	12.0mm	190mm	



Tibia Cannulated Drill Bit

Model number	Diameter	Working length	Picture
T400190	4.0mm	190mm	
T400190	4.5mm	190mm	
T400190	5.0mm	190mm	
T400190	5.5mm	190mm	
T400190	6.0mm	190mm	
T400190	6.5mm	190mm	
T400190	7.0mm	190mm	
T400190	7.5mm	190mm	
T400190	8.0mm	190mm	
T400190	8.5mm	190mm	
T400190	9.0mm	190mm	
T400190	9.5mm	190mm	
T400190	10.0mm	190mm	
T400190	10.5mm	190mm	
T400190	11.0mm	190mm	
T400190	11.5mm	190mm	
T400190	12.0mm	190mm	



C1001 Drill Guide Arm



C1002 Drill Guide Sleeve



C1003 ACL, Tip Aimer



C1004 ACL, Elbow Aimer



C1005 PCL, Tibial Aimer



C1006 PCL, Femoral Aimer

Arthroscope Instruments

Cruciate Ligament Reconstruction Instruments Set



C2001 PCL Elevator



C2002 Needle Drill, Trocar-Tip
ø2.3x452mm



C2003 Needle Drill, Trocar-Tip
ø2.3x400mm



C2004 Needle Drill, Drill-Tip
ø2.3x440mm



C2005 Needle Drill, Drill-Tip
ø2.3x400mm



C2006 Depth Probe
340mm



C3015 Tendon Strippers, Closed, ø5mm



C3016 Tendon Strippers, Closed, ø6mm



C3017 Tendon Strippers, Closed, ø7mm



C3018 Tendon Strippers, Closed, ø8mm



C3019 Tendon Strippers, Closed, ø9mm



C3025 Tendon Strippers, Slotted, ø5mm




C3026 Tendon Strippers, Slotted, ø6mm



C3027 Tendon Strippers, Slotted, ø7mm



C3028 Tendon Strippers, Slotted, ø8mm



C3029 Tendon Strippers, Slotted, ø9mm



C4001 Calibrated Probe



C4011 Aiming Device, ø5mm

C4012 Aiming Device, ø6mm

C4013 Aiming Device, ø7mm



C4021 Detacher, Upbiter



C4022 Detacher, Downbiter



C4031 Rasp, Downswept 20°/20°



C4032 Rasp, Upswept 20°/20°



C4041 Screwdriver

Arthroscope Instruments

Cruciate Ligament Reconstruction Instruments Set



- C5050 Tendon Cannula, $\varnothing$ 5.0mm
- C5055 Tendon Cannula, $\varnothing$ 5.5mm
- C5060 Tendon Cannula, $\varnothing$ 6.0mm
- C5065 Tendon Cannula, $\varnothing$ 6.5mm
- C5070 Tendon Cannula, $\varnothing$ 7.0mm
- C5075 Tendon Cannula, $\varnothing$ 7.5mm
- C5080 Tendon Cannula, $\varnothing$ 8.0mm
- C5085 Tendon Cannula, $\varnothing$ 8.5mm
- C5090 Tendon Cannula, $\varnothing$ 9.0mm
- C5095 Tendon Cannula, $\varnothing$ 9.5mm




C5101 Graft Preparation Working Station



C5111 Size Block, Closed
4.5-12mm



C5112 Size Block, Slotted
4.5-12mm



C6001 Muscle Tendon Scraper



C6002 Tie Tensioner



C6003 Bone Hammer

C7001 Suture Grasper
(Straight)
 $\varnothing$ 4.0x160mm



Arthroscope Instruments

Shoulder Instruments Set (handle A, C, D)

S1001 Suture Cutter
ø4.0x160mm



S1002 Loose Body Grasper
ø4.0x160mm



S1003 Suture Grasper
(straight)
ø4.0x160mm



S1004 Suture Grasper
(round)
ø4.0x160mm



S1005 Suture Grasper
(with tooth)
ø4.0x160mm



S1006 Tissue Forceps
(serration)
ø4.4x160mm



S1007 Penetrator
(straight)
ø2.7x160mm



S1008 Penetrator
(straight, upward)
ø2.7x160mm



S1009 Penetrator
(left)
ø2.7x160mm



S1010 Penetrator
(right)
ø2.7x160mm



S1011 Penetrator
(straight, upward)
ø3.4x160mm



S1012 Penetrator
(left)
ø3.4x160mm



S1013 Penetrator
(right)
ø3.4x160mm



S1014 Penetrator
(left, upward)
ø3.4x160mm



S1015 Penetrator
(right, upward)
ø3.4x160mm



S1016 Suture Cutter
(open)
ø3.2x160mm



S1017 Suture Cutter
(close)
ø3.2x160mm



Arthroscope Instruments

Shoulder Instruments Set



S2001 Suture Passer
ø4.0x170mm



S2002 Suture Passer Needle
ø1.5x237mm



S2011 Knot Pusher
ø3.2x160mm



S2012 Switching Stick
ø4.0x330mm



S2013 Suture Hook
(straight, 80°)
ø4.0x150mm



S2014 Suture Hook
(left, 45°)
ø4.0x150mm



S2015 Suture Hook
(right, 45°)
ø4.0x150mm



S2016 Elevator, 15°
ø5.0x160mm



S2017 Elevator, 35°
ø5.0x160mm



S2018 Curette
(round)
ø4.2x160mm



S2019 Glenoid Rasp
(flat)
ø6.5x160mm



S2020 Probe
ø3.0x160mm

S3001 Trocar, Thread
ø7.0x70mm

S3002 Trocar, Thread
ø8.0x70mm

S3003 Trocar, Thread
ø9.0x70mm



Arthroscope Instruments

Shoulder Instruments Set



S4001 Suture Hook
(straight, 15°)
ø4.0x150mm



S4002 Suture Hook
(straight, 45°)
ø4.0x150mm



S4003 Suture Hook
(straight, 80°)
ø4.0x150mm



S4004 Suture Hook
(left, 25°)
ø4.0x150mm



S4005 Suture Hook
(right, 25°)
ø4.0x150mm



S4006 Suture Hook
(left, 45°)
ø4.0x150mm



S4007 Suture Hook
(right, 45°)
ø4.0x150mm



S4008 Suture Hook
(left, 65°)
ø4.0x150mm



S4009 Suture Hook
(right, 65°)
ø4.0x150mm



S4010 Suture Hook
(left, 90°)
ø4.0x150mm



S4011 Suture Hook
(right, 90°)
ø4.0x150mm



Roll Type Set



S5011 Suture Hook
(right, 45°)
ø4.0x150mm



S5012 Suture Hook
(left, 45°)
ø4.0x150mm



S5013 Suture Hook
(straight, 45°)
ø4.0x150mm



S5014 Suture Hook
(left, 90°)
ø4.0x150mm



S5015 Suture Hook
(right, 90°)
ø4.0x150mm



S6001 Loop Wire
ø0.8mm

Arthroscope Instruments

Shoulder Instruments Set



S7001 Glenoid Rasp
(flat)
ø6.5x160mm



S7002 Glenoid Rasp
(semicircle)
ø5.0x150mm



S7003 Tissue Elevator
(upward)
ø5.0x150mm



S7004 Tissue Elevator
(downward)
ø5.0x150mm



S7005 Wire Hook
ø3.2x210mm



S7011 Chondro Pick 25°, ø5.0x170mm

S7012 Chondro Pick 40°, ø5.0x170mm

S7013 Chondro Pick 60°, ø5.0x170mm

S7014 Chondro Pick 80°, ø5.0x170mm



S7021 Curette
(round)
ø4.2x160mm



S7022 Curette
(closed)
ø3.2x125mm



S7031 Hole Puncher, ø2.9x130mm

S7032 Hole Puncher, ø3.5x130mm

S7033 Hole Puncher, ø4.5x130mm

S7034 Hole Puncher, ø5.1x130mm



S7041 Starter Drill, ø2.7x140mm

S7042 Starter Drill, ø3.8x140mm

S7043 Starter Drill, ø4.5x140mm



S7051 Reamer Drill, ø4.6x140mm

S7052 Reamer Drill, ø4.8x140mm

S7053 Reamer Drill, ø5.2x140mm



S7061 Guide Drill, ø2.8x190mm

S7062 Guide Drill, ø4.8x190mm



S7071 Thread Pin, ø1.5x250mm

S7072 Thread Pin, ø1.8x250mm

S7073 Thread Pin, ø2.4x250mm

S7074 Thread Pin, ø2.9x250mm



S7081 Guide Pin (round), ø2.1x235mm

S7082 Guide Pin (round), ø2.4x235mm

S7083 Guide Pin (round), ø3.0x235mm




S7091 Guide Pin (sharp), ø2.1x235mm


S7092 Guide Pin (sharp), ø2.4x235mm

S7093 Guide Pin (sharp), ø3.0x235mm

S7094 Guide Pin (sharp), ø3.7x235mm



C8001 Guide Cannula, ø1.7x210mm



C8011 Guide Cannula, ø5.5x6.5mm

C8012 Guide Cannula, ø3.8x5.0mm

C8013 Guide Cannula, ø3.0x4.5mm

Arthroscope Instruments

Hip Instruments Set (handle A, B, D)

H1001 Basket Punch
(straight)
ø3.2x185mm



H1002 Basket Punch
(straight, upward)
ø3.2x185mm



H1003 Basket Punch
(right)
ø3.2x185mm



H1004 Basket Punch
(left)
ø3.2x185mm



H1005 Basket Punch
(straight)
ø1.8x185mm



H1006 Basket Punch
(straight, upward)
ø1.8x185mm



H1007 Basket Punch
(right)
ø1.8x185mm



H1008 Basket Punch
(left)
ø1.8x185mm



H1009 Basket Punch
(right, backbiter)
ø2.4x185mm



H1010 Basket Punch
(left, backbiter)
ø2.4x185mm



H1011 Basket Punch
(backbiter)
ø2.6x185mm



H1012 Loose Body Grasper
ø4.0x185mm



H1013 Tissue Catcher
(straight)
ø4.0x185mm



H1014 Tissue Catcher
(right)
ø4.0x185mm



H1015 Tissue Catcher
(left)
ø4.0x185mm



H1016 Penetrator
(straight)
ø2.7x185mm



H1017 Penetrator
(straight, upward)
ø2.7x185mm



H1018 Penetrator
(straight, downward)
ø3.4x185mm



H1019 Penetrator
(right, downward)
ø3.4x185mm



H1020 Penetrator
(left, downward)
ø3.4x185mm



H1021 Suture Cutter
ø4.0x185mm



H1022 Suture Grasper
(round)
ø2.6x185mm



H1023 Suture Grasper
(small)
ø4.0x185mm



H1024 Suture Grasper
(big)
ø4.0x185mm



Arthroscope Instruments

Hip Instruments Set



H2001 Bone Rasp, $\varnothing 6.5 \times 220$ mm



H2002 Elevator, $\varnothing 5.0 \times 220$ mm



H2003 Curette, $\varnothing 5.0 \times 220$ mm



H2004 Curette (straight, upward), $\varnothing 5.0 \times 220$ mm



H2011 Chondro Pick 25°, $\varnothing 5.0 \times 220$ mm

H2012 Chondro Pick 40°, $\varnothing 5.0 \times 220$ mm

H2013 Chondro Pick 60°, $\varnothing 5.0 \times 220$ mm

H2014 Chondro Pick 80°, $\varnothing 5.0 \times 220$ mm



H2021 Cutting Tip (right), $\varnothing 3.8 \times 250$ mm



H2022 Cutting Tip (left), $\varnothing 3.8 \times 250$ mm



H2023 Puncture Needle, $\varnothing 1.6 \times 190$ mm




H2024 Puncture Needle, $\varnothing 2.5 \times 215$ mm

H2031 Guide Needle, $\varnothing 0.8$ mm

H2031 Guide Needle, $\varnothing 1.0$ mm

H2031 Guide Needle, $\varnothing 1.2$ mm


H2031 Guide Needle, $\varnothing 1.5$ mm



H2041 Guide Needle, $\varnothing 4.3 \times 140$ mm

H2042 Guide Needle, $\varnothing 4.9 \times 140$ mm

H2043 Guide Needle, $\varnothing 5.5 \times 140$ mm



H2051 Sem-Pipe (upper slotting)
 $\varnothing 8.0 \times 135$ mm




H2052 Sem-Pipe (upper slotting)
 $\varnothing 8.0 \times 120$ mm



H2061 Trocar Sheath, $\varnothing 5.0 \times 155$ mm

H2052 Trocar, $\varnothing 4.3 \times 160$ mm



H2053 Trocar Sheath, $\varnothing 5.6 \times 155$ mm

H2054 Trocar, $\varnothing 4.9 \times 160$ mm



H2055 Trocar Sheath, $\varnothing 6.4 \times 155$ mm

H2056 Trocar, $\varnothing 5.5 \times 160$ mm



H2057 Water Valve Bridge



H2071 Knot Pusher, $\varnothing 3.2 \times 225$ mm



H2072 Guide Bar (hollow)
 $\varnothing 4.5 \times 330$ mm



**SURGINEST
GLOBAL**

INSTRUMENTAL Y EQUIPOS MÉDICOS