


# LOCKING PLATE SYSTEM

Instrumental de Acero Quirúrgico

CATÁLOGO



**SURGINEST  
GLOBAL**  
INSTRUMENTAL Y EQUIPOS MÉDICOS



VA Proximal Lateral Tibia Locking Plate



| REF       | Holes      | Thickness | Width  | Length |
|-----------|------------|-----------|--------|--------|
| 5200-0701 | 4 holes L  | 3.8mm     | 14.2mm | 87mm   |
| 5200-0702 | 6 holes L  | 3.8mm     | 14.2mm | 117mm  |
| 5200-0703 | 8 holes L  | 3.8mm     | 14.2mm | 147mm  |
| 5200-0704 | 10 holes L | 3.8mm     | 14.2mm | 177mm  |
| 5200-0707 | 4 holes R  | 3.8mm     | 14.2mm | 87mm   |
| 5200-0708 | 6 holes R  | 3.8mm     | 14.2mm | 117mm  |
| 5200-0709 | 8 holes R  | 3.8mm     | 14.2mm | 147mm  |
| 5200-0710 | 10 holes R | 3.8mm     | 14.2mm | 177mm  |

Use 3.5 VA Locking Screw/3.5 Cortical Screw

VA Distal Tibia Locking L-Plate



| REF       | Holes     | Thickness | Width  | Length |
|-----------|-----------|-----------|--------|--------|
| 5200-0501 | 4 holes L | 2.2mm     | 10.8mm | 72mm   |
| 5200-0502 | 6 holes L | 2.2mm     | 10.8mm | 90mm   |
| 5200-0503 | 4 holes R | 2.2mm     | 10.8mm | 72mm   |
| 5200-0504 | 6 holes R | 2.2mm     | 10.8mm | 90mm   |

Use 2.7 VA Locking Screw/2.7 Cortical Screw

VA Distal Lateral Tibia Locking Plate



| REF       | Holes      | Thickness | Width  | Length |
|-----------|------------|-----------|--------|--------|
| 5100-5801 | 4 holes L  | 3.8mm     | 13.5mm | 82mm   |
| 5100-5802 | 6 holes L  | 3.8mm     | 13.5mm | 112mm  |
| 5100-5803 | 8 holes L  | 3.8mm     | 13.5mm | 142mm  |
| 5100-5804 | 10 holes L | 3.8mm     | 13.5mm | 172mm  |
| 5100-5805 | 12 holes L | 3.8mm     | 13.5mm | 202mm  |
| 5100-5809 | 4 holes R  | 3.8mm     | 13.5mm | 82mm   |
| 5100-5810 | 6 holes R  | 3.8mm     | 13.5mm | 112mm  |
| 5100-5811 | 8 holes R  | 3.8mm     | 13.5mm | 142mm  |
| 5100-5812 | 10 holes R | 3.8mm     | 13.5mm | 172mm  |
| 5100-5813 | 12 holes R | 3.8mm     | 13.5mm | 202mm  |

Use 2.7 VA Locking Screw/3.5 VA Locking Screw/3.5 Cortical Screw

VA Distal Tibia Locking T-Plate



| REF       | Holes   | Thickness | Width  | Length |
|-----------|---------|-----------|--------|--------|
| 5200-0601 | 4 holes | 2.2mm     | 10.8mm | 72mm   |
| 5200-0602 | 6 holes | 2.2mm     | 10.8mm | 90mm   |

Use 2.7 VA Locking Screw/2.7 Cortical Screw

VA Distal Lateral Tibia Locking Plate



| REF       | Holes      | Thickness | Width  | Length |
|-----------|------------|-----------|--------|--------|
| 5100-5701 | 4 holes L  | 4mm       | 13.5mm | 112mm  |
| 5100-5702 | 6 holes L  | 4mm       | 13.5mm | 142mm  |
| 5100-5703 | 8 holes L  | 4mm       | 13.5mm | 172mm  |
| 5100-5704 | 10 holes L | 4mm       | 13.5mm | 202mm  |
| 5100-5705 | 12 holes L | 4mm       | 13.5mm | 232mm  |
| 5100-5708 | 4 holes R  | 4mm       | 13.5mm | 112mm  |
| 5100-5709 | 6 holes R  | 4mm       | 13.5mm | 142mm  |
| 5100-5710 | 8 holes R  | 4mm       | 13.5mm | 172mm  |
| 5100-5711 | 10 holes R | 4mm       | 13.5mm | 202mm  |
| 5100-5712 | 12 holes R | 4mm       | 13.5mm | 232mm  |

Use 2.7 VA Locking Screw/3.5 VA Locking Screw/3.5 Cortical Screw

Distal Lateral Fibular Locking Plate



| REF       | Holes     | Thickness | Width | Length |
|-----------|-----------|-----------|-------|--------|
| 5200-0401 | 3 holes L | 2.5mm     | 10mm  | 79mm   |
| 5200-0402 | 4 holes L | 2.5mm     | 10mm  | 92mm   |
| 5200-0403 | 5 holes L | 2.5mm     | 10mm  | 105mm  |
| 5200-0404 | 6 holes L | 2.5mm     | 10mm  | 118mm  |
| 5200-0405 | 7 holes L | 2.5mm     | 10mm  | 131mm  |
| 5200-0406 | 9 holes L | 2.5mm     | 10mm  | 157mm  |
| 5200-0410 | 3 holes R | 2.5mm     | 10mm  | 79mm   |
| 5200-0411 | 4 holes R | 2.5mm     | 10mm  | 92mm   |
| 5200-0412 | 5 holes R | 2.5mm     | 10mm  | 105mm  |
| 5200-0413 | 6 holes R | 2.5mm     | 10mm  | 118mm  |
| 5200-0414 | 7 holes R | 2.5mm     | 10mm  | 131mm  |
| 5200-0415 | 9 holes R | 2.5mm     | 10mm  | 157mm  |

Use 2.7 VA Locking Screw/2.7 Cortical Screw / 3.5 Cortical Screw



3.5mm Pediatric Eight-plate



| REF        | Holes        | Length |
|------------|--------------|--------|
| 01.1261-01 | 2 holes*10mm | 10mm   |
| 01.1261-02 | 2 holes*12mm | 12mm   |
| 01.1261-03 | 2 holes*16mm | 16mm   |

Use 3.5mm Cannulated Screw Full-threaded

3.5mm Cannulated Screw Full-threaded



| REF        | Specification | REF        | Specification |
|------------|---------------|------------|---------------|
| 01.1261-04 | 3.5*10mm      | 01.1261-12 | 3.5*26mm      |
| 01.1261-05 | 3.5*12mm      | 01.1261-13 | 3.5*28mm      |
| 01.1261-06 | 3.5*14mm      | 01.1261-14 | 3.5*30mm      |
| 01.1261-07 | 3.5*16mm      | 01.1261-15 | 3.5*32mm      |
| 01.1261-08 | 3.5*18mm      | 01.1261-16 | 3.5*35mm      |
| 01.1261-09 | 3.5*20mm      | 01.1261-17 | 3.5*40mm      |
| 01.1261-10 | 3.5*22mm      | 01.1261-18 | 3.5*45mm      |
| 01.1261-11 | 3.5*24mm      | 01.1261-19 | 3.5*50mm      |

4.5 Pediatric Hinge Locking plate



| REF        | Holes        | Length |
|------------|--------------|--------|
| 01.1264-01 | 2 holes*12mm | 12mm   |
| 01.1264-02 | 2 holes*16mm | 16mm   |

Use 4.5mm Cannulated Screw Full-threaded

4.5mm Pediatric Eight-plate



| REF        | Holes        | Length |
|------------|--------------|--------|
| 01.1262-01 | 2 holes*12mm | 12mm   |
| 01.1262-02 | 2 holes*16mm | 16mm   |
| 01.1262-03 | 2 holes*20mm | 20mm   |
| 01.1262-04 | 2 holes*25mm | 25mm   |

Use 4.5mm Cannulated Screw Full-threaded

4.5mm Pediatric Plate



| REF        | Holes        | Length |
|------------|--------------|--------|
| 01.1263-01 | 4 holes*12mm | 12mm   |
| 01.1263-02 | 4 holes*16mm | 16mm   |
| 01.1263-03 | 4 holes*20mm | 20mm   |

Use 4.5mm Cannulated Screw Full-threaded

4.5mm Cannulated Screw Full-threaded



| REF        | Specification | REF        | Specification |
|------------|---------------|------------|---------------|
| 01.1262-05 | 4.5*10        | 01.1262-14 | 4.5*28        |
| 01.1262-06 | 4.5*12        | 01.1262-15 | 4.5*30        |
| 01.1262-07 | 4.5*14        | 01.1262-16 | 4.5*32        |
| 01.1262-08 | 4.5*16        | 01.1262-21 | 4.5*34        |
| 01.1262-09 | 4.5*18        | 01.1262-22 | 4.5*36        |
| 01.1262-10 | 4.5*20        | 01.1262-23 | 4.5*38        |
| 01.1262-11 | 4.5*22        | 01.1262-18 | 4.5*40        |
| 01.1262-12 | 4.5*24        | 01.1262-19 | 4.5*45        |
| 01.1262-13 | 4.5*26        | 01.1262-20 | 4.5*50        |

3.5mm Pediatric Sliding Barrel Locking Plate



| REF        | Holes     | Length |
|------------|-----------|--------|
| 01.1273-01 | 4 holes L | 80mm   |
| 01.1273-02 | 4 holes R | 80mm   |

Use 3.5mm Locking Screw

3.5mm Pediatric Sliding Straight Locking Plate



| REF        | Holes    | Length |
|------------|----------|--------|
| 01.1274-01 | 5 holes  | 60mm   |
| 01.1274-02 | 7 holes  | 80mm   |
| 01.1274-03 | 9 holes  | 100mm  |
| 01.1274-04 | 11 holes | 120mm  |

Use 3.5mm Locking Screw

2.7mm LCP Pediatric Hip Plate



| REF        | Holes         | Length |
|------------|---------------|--------|
| 01.7003-01 | 3 holes*100°  | 48.2mm |
| 01.7003-02 | 3 holes*110°  | 48.2mm |
| 01.7003-03 | 4 holes*120°  | 58.2mm |
| 01.7003-04 | 6 holes*120°  | 78.2mm |
| 01.7003-05 | 3 holes*130°  | 48.2mm |
| 01.7003-06 | 10 holes*130° | 118.2m |

Use 2.7mm Locking Screw, 2.7mm Cortical Screw

3.5mm LCP Pediatric Hip Plate



| REF        | Holes         | Length |
|------------|---------------|--------|
| 01.7004-01 | 3 holes*90°   | 75mm   |
| 01.7004-02 | 4 holes*90°   | 87mm   |
| 01.7004-03 | 5 holes*90°   | 99mm   |
| 01.7004-04 | 6 holes*90°   | 111mm  |
| 01.7004-05 | 7 holes*90°   | 123mm  |
| 01.7004-06 | 8 holes*90°   | 135mm  |
| 01.7004-07 | 9 holes*90°   | 147mm  |
| 01.7004-08 | 3 holes*100°  | 75mm   |
| 01.7004-09 | 3 holes*110°  | 75mm   |
| 01.7004-10 | 3 holes*120°  | 62mm   |
| 01.7004-11 | 4 holes*120°  | 75mm   |
| 01.7004-12 | 6 holes*120°  | 101mm  |
| 01.7004-13 | 8 holes*120°  | 127mm  |
| 01.7004-14 | 10 holes*120° | 153mm  |
| 01.7004-15 | 3 holes*130°  | 75mm   |
| 01.7004-16 | 4 holes*130°  | 87mm   |
| 01.7004-17 | 5 holes*130°  | 99mm   |
| 01.7004-18 | 7 holes*130°  | 123mm  |
| 01.7004-19 | 9 holes*130°  | 147mm  |
| 01.7004-20 | 3 holes*150°  | 64mm   |

Use 3.5mm Locking Screw, 3.5mm Cortical Screw



### 5.0mm LCP Pediatric Hip Plate



| REF        | Holes         | Length |
|------------|---------------|--------|
| 01.7005-01 | 3 holes*90°   | 96mm   |
| 01.7005-02 | 4 holes*90°   | 112mm  |
| 01.7005-03 | 5 holes*90°   | 122mm  |
| 01.7005-04 | 7 holes*90°   | 220mm  |
| 01.7005-05 | 9 holes*90°   | 252mm  |
| 01.7005-06 | 3 holes*100°  | 101mm  |
| 01.7005-07 | 3 holes*110°  | 100mm  |
| 01.7005-08 | 3 holes*120°  | 99mm   |
| 01.7005-09 | 4 holes*120°  | 115mm  |
| 01.7005-10 | 6 holes*120°  | 147mm  |
| 01.7005-11 | 8 holes*120°  | 179mm  |
| 01.7005-12 | 10 holes*120° | 211mm  |
| 01.7005-13 | 3 holes*130°  | 96mm   |
| 01.7005-14 | 4 holes*130°  | 112mm  |
| 01.7005-15 | 5 holes*130°  | 122mm  |
| 01.7005-16 | 7 holes*130°  | 220mm  |
| 01.7005-17 | 9 holes*130°  | 252mm  |
| 01.7005-18 | 3 holes*140°  | 99mm   |
| 01.7005-19 | 4 holes*140°  | 115mm  |
| 01.7005-20 | 6 holes*140°  | 147mm  |
| 01.7005-21 | 3 holes*150°  | 77.5mm |

Use 5.0mm Locking Screw, 4.5mm Cortical Screw

### 2.7mm LCP Pediatric Straight Plate



| REF        | Holes    | Length |
|------------|----------|--------|
| 01.7010-01 | 5 holes  | 49.5mm |
| 01.7010-02 | 6 holes  | 58.5mm |
| 01.7010-03 | 7 holes  | 67.5mm |
| 01.7010-04 | 8 holes  | 76.5mm |
| 01.7010-05 | 9 holes  | 85.5mm |
| 01.7010-06 | 10 holes | 94.5mm |

Use 2.7mm Locking Screw, 2.7mm Cortical Screw

### 2.7mm Cortex Screw, T8 Stardrive Self-tapping



| REF       | Specification | REF       | Specification |
|-----------|---------------|-----------|---------------|
| 030390010 | 2.7*10mm      | 030390022 | 2.7*22mm      |
| 030390012 | 2.7*12mm      | 030390024 | 2.7*24mm      |
| 030390014 | 2.7*14mm      | 030390026 | 2.7*26mm      |
| 030390016 | 2.7*16mm      | 030390028 | 2.7*28mm      |
| 030390018 | 2.7*18mm      | 030390030 | 2.7*30mm      |
| 030390020 | 2.7*20mm      |           |               |

### 2.7mm Locking Screw, T8 Stardrive Self-tapping



| REF       | Specification | REF       | Specification |
|-----------|---------------|-----------|---------------|
| 030340010 | 2.7*10mm      | 030340026 | 2.7*26mm      |
| 030340012 | 2.7*12mm      | 030340028 | 2.7*28mm      |
| 030340014 | 2.7*14mm      | 030340030 | 2.7*30mm      |
| 030340016 | 2.7*16mm      | 030340032 | 2.7*32mm      |
| 030340018 | 2.7*18mm      | 030340034 | 2.7*34mm      |
| 030340020 | 2.7*20mm      | 030340036 | 2.7*36mm      |
| 030340022 | 2.7*22mm      | 030340038 | 2.7*38mm      |
| 030340024 | 2.7*24mm      | 030340040 | 2.7*40mm      |

### 3.5mm Locking Screw (Stardrive)



| REF       | Specification | REF       | Specification |
|-----------|---------------|-----------|---------------|
| 5100-3901 | 3.5*12mm      | 5100-3917 | 3.5*44mm      |
| 5100-3902 | 3.5*14mm      | 5100-3918 | 3.5*46mm      |
| 5100-3903 | 3.5*16mm      | 5100-3919 | 3.5*48mm      |
| 5100-3904 | 3.5*18mm      | 5100-3920 | 3.5*50mm      |
| 5100-3905 | 3.5*20mm      | 5100-3921 | 3.5*52mm      |
| 5100-3906 | 3.5*22mm      | 5100-3922 | 3.5*54mm      |
| 5100-3907 | 3.5*24mm      | 5100-3923 | 3.5*56mm      |
| 5100-3908 | 3.5*26mm      | 5100-3924 | 3.5*58mm      |
| 5100-3909 | 3.5*28mm      | 5100-3925 | 3.5*60mm      |
| 5100-3910 | 3.5*30mm      | 5100-3926 | 3.5*65mm      |
| 5100-3911 | 3.5*32mm      | 5100-3927 | 3.5*70mm      |
| 5100-3912 | 3.5*34mm      | 5100-3928 | 3.5*75mm      |
| 5100-3913 | 3.5*36mm      | 5100-3929 | 3.5*80mm      |
| 5100-3914 | 3.5*38mm      | 5100-3930 | 3.5*85mm      |
| 5100-3915 | 3.5*40mm      | 5100-3931 | 3.5*90mm      |
| 5100-3916 | 3.5*42mm      | 5100-3932 | 3.5*95mm      |

### 5.0mm Locking Screw (Stardrive)



| REF       | Specification | REF       | Specification |
|-----------|---------------|-----------|---------------|
| 5100-4001 | 5.0*22mm      | 5100-4015 | 5.0*50mm      |
| 5100-4002 | 5.0*24mm      | 5100-4016 | 5.0*52mm      |
| 5100-4003 | 5.0*26mm      | 5100-4017 | 5.0*54mm      |
| 5100-4004 | 5.0*28mm      | 5100-4018 | 5.0*56mm      |
| 5100-4005 | 5.0*30mm      | 5100-4019 | 5.0*58mm      |
| 5100-4006 | 5.0*32mm      | 5100-4020 | 5.0*60mm      |
| 5100-4007 | 5.0*34mm      | 5100-4021 | 5.0*65mm      |
| 5100-4008 | 5.0*36mm      | 5100-4022 | 5.0*70mm      |
| 5100-4009 | 5.0*38mm      | 5100-4023 | 5.0*75mm      |
| 5100-4010 | 5.0*40mm      | 5100-4024 | 5.0*80mm      |
| 5100-4011 | 5.0*42mm      | 5100-4025 | 5.0*85mm      |
| 5100-4012 | 5.0*44mm      | 5100-4026 | 5.0*90mm      |
| 5100-4013 | 5.0*46mm      | 5100-4027 | 5.0*95mm      |
| 5100-4014 | 5.0*48mm      |           |               |

### 3.5mm Cortical Screw (Stardrive)



| REF       | Specification | REF       | Specification |
|-----------|---------------|-----------|---------------|
| 5100-4101 | 3.5*12        | 5100-4112 | 3.5*34        |
| 5100-4102 | 3.5*14        | 5100-4113 | 3.5*36        |
| 5100-4103 | 3.5*16        | 5100-4114 | 3.5*38        |
| 5100-4104 | 3.5*18        | 5100-4115 | 3.5*40        |
| 5100-4105 | 3.5*20        | 5100-4116 | 3.5*42        |
| 5100-4106 | 3.5*22        | 5100-4117 | 3.5*44        |
| 5100-4107 | 3.5*24        | 5100-4118 | 3.5*46        |
| 5100-4108 | 3.5*26        | 5100-4119 | 3.5*48        |
| 5100-4109 | 3.5*28        | 5100-4120 | 3.5*50        |
| 5100-4110 | 3.5*30        | 5100-4121 | 3.5*55        |
| 5100-4111 | 3.5*32        | 5100-4122 | 3.5*60        |

### 4.5mm Cortical Screw (Stardrive)



| REF       | Specification | REF       | Specification |
|-----------|---------------|-----------|---------------|
| 5100-4201 | 4.5*22mm      | 5100-4211 | 4.5*42mm      |
| 5100-4202 | 4.5*24mm      | 5100-4212 | 4.5*44mm      |
| 5100-4203 | 4.5*26mm      | 5100-4213 | 4.5*46mm      |
| 5100-4204 | 4.5*28mm      | 5100-4214 | 4.5*48mm      |
| 5100-4205 | 4.5*30mm      | 5100-4215 | 4.5*50mm      |
| 5100-4206 | 4.5*32mm      | 5100-4216 | 4.5*52mm      |
| 5100-4207 | 4.5*34mm      | 5100-4217 | 4.5*54mm      |
| 5100-4208 | 4.5*36mm      | 5100-4218 | 4.5*56mm      |
| 5100-4209 | 4.5*38mm      | 5100-4219 | 4.5*58mm      |
| 5100-4210 | 4.5*40mm      | 5100-4220 | 4.5*60mm      |



### 1.5S Mini Locking Plate



| REF       | Holes   | Thickness |
|-----------|---------|-----------|
| 021260004 | 4 holes | 23mm      |
| 021260006 | 6 holes | 36mm      |

Thickness:1.0mm, Width: 4.5mm

### 1.5S Mini Reconstruction Locking Plate



| REF       | Holes    | Thickness |
|-----------|----------|-----------|
| 021270004 | 4 holes  | 20mm      |
| 021270006 | 6 holes  | 30mm      |
| 021270008 | 8 holes  | 40mm      |
| 021270010 | 10 holes | 50mm      |

Thickness:1.0mm, Width: 4.5mm

### 1.5S Mini T Locking Plate (3 holes)



| REF       | Holes   | Thickness |
|-----------|---------|-----------|
| 021280004 | 4 holes | 25mm      |
| 021280006 | 6 holes | 35mm      |
| 021280008 | 8 holes | 45mm      |

Thickness:1.0mm, Width: 4.5mm

### 1.5S Mini T Locking Plate (4 holes)



| REF       | Holes   | Thickness |
|-----------|---------|-----------|
| 020980004 | 4 holes | 25mm      |
| 020980006 | 6 holes | 35mm      |
| 020980008 | 8 holes | 45mm      |

Thickness:1.0mm, Width: 4.5mm

### 1.5S Mini Y Locking Plate



| REF       | Holes   | Thickness |
|-----------|---------|-----------|
| 020990004 | 4 holes | 26mm      |
| 020990006 | 6 holes | 36mm      |
| 020990008 | 8 holes | 46mm      |

Thickness:1.0mm, Width: 4.5mm

### 1.5S Mini Condylar Locking Plate



| REF       | Holes   | Thickness |
|-----------|---------|-----------|
| 021000004 | 4 holes | 26mm      |
| 021000006 | 6 holes | 36mm      |
| 021000008 | 8 holes | 46mm      |

Thickness:1.0mm, Width: 4.5mm

### 1.5S Mini Strut Locking Plate, oblique-angled



| REF       | Holes     |
|-----------|-----------|
| 021011000 | 8 holes L |
| 021012000 | 8 holes R |

Thickness:1.0mm

### 2.0S Mini Locking Plate



| REF       | Holes   | Thickness |
|-----------|---------|-----------|
| 021020004 | 4 holes | 31mm      |
| 021020005 | 5 holes | 38mm      |
| 021020006 | 6 holes | 45mm      |
| 021020008 | 8 holes | 59mm      |

Thickness:1.5mm, Width: 5mm

### 2.0S Mini Reconstruction Locking Plate



| REF       | Holes    | Thickness |
|-----------|----------|-----------|
| 021030004 | 4 holes  | 27mm      |
| 021030006 | 6 holes  | 40mm      |
| 021030008 | 8 holes  | 54mm      |
| 021030010 | 10 holes | 67mm      |

Thickness:1.5mm, Width: 5mm

### 2.0S Mini T Locking Plate



| REF       | Holes   | Thickness |
|-----------|---------|-----------|
| 021040003 | 3 holes | 26mm      |
| 021040005 | 5 holes | 39mm      |
| 021040007 | 7 holes | 53mm      |

Thickness:1.5mm, Width: 5mm

### 2.0S Mini Y Locking Plate



| REF       | Holes   | Thickness |
|-----------|---------|-----------|
| 021050003 | 3 holes | 27mm      |
| 021050005 | 5 holes | 41mm      |
| 021050007 | 7 holes | 55mm      |

Thickness:1.5mm, Width: 5mm

### 2.0S Mini T Reconstruction Locking Plate



| REF       | Holes   | Thickness |
|-----------|---------|-----------|
| 021060003 | 3 holes | 26mm      |
| 021060005 | 5 holes | 39mm      |
| 021060007 | 7 holes | 53mm      |

Thickness:1.5mm, Width: 5mm

### 2.0S Mini Condylar Locking Plate



| REF       | Holes   | Thickness |
|-----------|---------|-----------|
| 021070003 | 3 holes | 27mm      |
| 021070005 | 5 holes | 41mm      |
| 021070007 | 7 holes | 55mm      |

Thickness:1.5mm, Width: 5mm



### 2.4S Mini Locking Plate



| REF       | Holes    | Thickness |
|-----------|----------|-----------|
| 021080004 | 4 holes  | 36mm      |
| 021080005 | 5 holes  | 44mm      |
| 021080006 | 6 holes  | 52mm      |
| 021080008 | 8 holes  | 68mm      |
| 021080010 | 10 holes | 84mm      |

Thickness:1.5mm, Width: 6.5mm


### 2.4S Mini Reconstruction Locking Plate



| REF       | Holes    | Thickness |
|-----------|----------|-----------|
| 021090004 | 4 holes  | 30mm      |
| 021090006 | 6 holes  | 44mm      |
| 021090008 | 8 holes  | 58mm      |
| 021090010 | 10 holes | 72mm      |

Thickness:1.3mm, Width: 6.5mm

### 2.4S Mini T Locking Plate



| REF       | Holes   | Thickness |
|-----------|---------|-----------|
| 021100003 | 3 holes | 28mm      |
| 021100005 | 5 holes | 43mm      |
| 021100007 | 7 holes | 57mm      |

Thickness:1.5mm, Width: 6.5mm


### 2.4S Mini Y Locking Plate



| REF       | Holes   | Thickness |
|-----------|---------|-----------|
| 021110003 | 3 holes | 31mm      |
| 021110005 | 5 holes | 46mm      |
| 021110007 | 7 holes | 60mm      |

Thickness:1.5mm, Width: 6.5mm

### 2.4S Mini T Locking Plate



| REF       | Holes   | Thickness |
|-----------|---------|-----------|
| 021120003 | 3 holes | 29mm      |
| 021120005 | 5 holes | 44mm      |
| 021120007 | 7 holes | 58mm      |

Thickness:1.5mm, Width: 6.5mm

### 2.4S Mini T Locking Plate



| REF       | Holes   | Thickness |
|-----------|---------|-----------|
| 021130003 | 3 holes | 30mm      |
| 021130005 | 5 holes | 45mm      |
| 021130007 | 7 holes | 59mm      |

Thickness:1.5mm, Width: 6.5mm

### 2.7S Mini Locking Plate



| REF       | Holes    | Thickness |
|-----------|----------|-----------|
| 021140004 | 4 holes  | 40mm      |
| 021140005 | 5 holes  | 49mm      |
| 021140006 | 6 holes  | 58mm      |
| 021140008 | 8 holes  | 76mm      |
| 021140010 | 10 holes | 94mm      |

Thickness:2.0mm, Width: 7.5mm

### 2.7S Mini Reconstruction Locking Plate



| REF       | Holes    | Thickness |
|-----------|----------|-----------|
| 021150004 | 4 holes  | 33mm      |
| 021150006 | 6 holes  | 49mm      |
| 021150008 | 8 holes  | 65mm      |
| 021150010 | 10 holes | 81mm      |

Thickness:1.5mm, Width: 6.5mm

### 2.7S Mini T Locking Plate



| REF       | Holes   | Thickness |
|-----------|---------|-----------|
| 021160003 | 3 holes | 32mm      |
| 021160004 | 4 holes | 40mm      |

Thickness:1.5mm, Width: 7mm

### 2.7S Mini L Locking Plate



| REF       | Holes     | Thickness |
|-----------|-----------|-----------|
| 021181003 | 3 holes L | 32mm      |
| 021181004 | 4 holes L | 40mm      |
| 021181005 | 3 holes R | 32mm      |
| 021181006 | 4 holes R | 40mm      |

Thickness:1.5mm, Width: 7.5mm


### 2.7S Mini L Locking Plate



| REF       | Holes     | Thickness |
|-----------|-----------|-----------|
| 021171003 | 3 holes L | 35mm      |
| 021171004 | 4 holes L | 42mm      |
| 021171005 | 3 holes R | 35mm      |
| 021171006 | 4 holes R | 42mm      |

Thickness:1.5mm, Width: 7.5mm

### 2.7S Mini Condylar Locking Plate



| REF       | Holes   | Thickness |
|-----------|---------|-----------|
| 021190003 | 3 holes | 34mm      |
| 021190005 | 5 holes | 50mm      |
| 021190007 | 7 holes | 66mm      |

Thickness:1.5mm, Width: 7.0mm



### 1/3 Tubular Locking Plate



| REF       | Holes    | Thickness | Width | Length |
|-----------|----------|-----------|-------|--------|
| 5100-0201 | 5 holes  | 2mm       | 10mm  | 71mm   |
| 5100-0202 | 6 holes  | 2mm       | 10mm  | 84mm   |
| 5100-0203 | 7 holes  | 2mm       | 10mm  | 97mm   |
| 5100-0204 | 8 holes  | 2mm       | 10mm  | 110mm  |
| 5100-0205 | 9 holes  | 2mm       | 10mm  | 123mm  |
| 5100-0206 | 10 holes | 2mm       | 10mm  | 136mm  |
| 5100-0207 | 12 holes | 2mm       | 10mm  | 162mm  |

Use 3.5 Locking Screw/3.5 Cortical Screw

### Ulna & Radius Locking Plate



| REF       | Holes    | Thickness | Width | Length |
|-----------|----------|-----------|-------|--------|
| 5100-0301 | 6 holes  | 3.2mm     | 11mm  | 92mm   |
| 5100-0302 | 7 holes  | 3.2mm     | 11mm  | 105mm  |
| 5100-0303 | 8 holes  | 3.2mm     | 11mm  | 118mm  |
| 5100-0304 | 9 holes  | 3.2mm     | 11mm  | 131mm  |
| 5100-0305 | 10 holes | 3.2mm     | 11mm  | 144mm  |
| 5100-0306 | 12 holes | 3.2mm     | 11mm  | 170mm  |

Use 3.5 Locking Screw/3.5 Cortical Screw

### Acromion of Scapula Locking Plate



| REF        | Holes     | Thickness | Width  | Length |
|------------|-----------|-----------|--------|--------|
| 01.1154-01 | 9 holes L | 1.5mm     | 34.2mm | 77.3mm |
| 01.1154-11 | 9 holes R | 1.5mm     | 34.2mm | 77.3mm |

Use 2.7 Locking Screw

### Scapula Medial Locking Plate



| REF        | Holes      | Thickness | Width | Length |
|------------|------------|-----------|-------|--------|
| 01.1157-01 | 10 holes L | 1.5mm     | 82mm  | 82mm   |
| 01.1157-11 | 10 holes R | 1.5mm     | 82mm  | 82mm   |

Use 2.7 Locking Screw

### Scapula Lateral Locking Plate



| REF        | Holes     | Thickness | Width | Length |
|------------|-----------|-----------|-------|--------|
| 01.1155-01 | 6 holes L | 1.5mm     | 8mm   | 89mm   |
| 01.1155-11 | 6 holes R | 1.5mm     | 8mm   | 89mm   |

Use 2.7 Locking Screw

### Glenoid Cavity Locking Plate



### S-clavicle Locking Plate



### S-clavicle Locking Plate-III



### S-clavicle Locking Plate-II



### Distal Clavicular Locking Plate-I



| REF        | Holes     | Thickness | Width  | Length |
|------------|-----------|-----------|--------|--------|
| 01.1156-01 | 4 holes L | 1.5mm     | 20.1mm | 42mm   |
| 01.1156-11 | 4 holes R | 1.5mm     | 20.1mm | 42mm   |

Use 2.7 Locking Screw

| REF       | Holes      | Thickness | Width  | Length |
|-----------|------------|-----------|--------|--------|
| 5100-0401 | 5 holes L  | 3mm       | 10.5mm | 54mm   |
| 5100-0402 | 6 holes L  | 3mm       | 10.5mm | 67mm   |
| 5100-0403 | 7 holes L  | 3mm       | 10.5mm | 81mm   |
| 5100-0404 | 8 holes L  | 3mm       | 10.5mm | 94mm   |
| 5100-0405 | 9 holes L  | 3mm       | 10.5mm | 107mm  |
| 5100-0406 | 10 holes L | 3mm       | 10.5mm | 131mm  |
| 5100-0407 | 5 holes R  | 3mm       | 10.5mm | 54mm   |
| 5100-0408 | 6 holes R  | 3mm       | 10.5mm | 67mm   |
| 5100-0409 | 7 holes R  | 3mm       | 10.5mm | 81mm   |
| 5100-0410 | 8 holes R  | 3mm       | 10.5mm | 94mm   |
| 5100-0411 | 9 holes R  | 3mm       | 10.5mm | 107mm  |
| 5100-0412 | 10 holes R | 3mm       | 10.5mm | 131mm  |

Use 3.5 Locking Screw/3.5 Cortical Screw/4.0 Cancellous Screw

| REF       | Holes   | Thickness | Width | Length |
|-----------|---------|-----------|-------|--------|
| 5100-8701 | 6 holes | 3.2mm     | 10mm  | 92mm   |
| 5100-8702 | 7 holes | 3.2mm     | 10mm  | 105mm  |
| 5100-8703 | 8 holes | 3.2mm     | 10mm  | 119mm  |
| 5100-8704 | 6 holes | 3.2mm     | 10mm  | 92mm   |
| 5100-8705 | 7 holes | 3.2mm     | 10mm  | 105mm  |
| 5100-8706 | 8 holes | 3.2mm     | 10mm  | 119mm  |

Use 3.5 Locking Screw/3.5 Cortical Screw

| REF       | Holes      | Thickness | Width  | Length  |
|-----------|------------|-----------|--------|---------|
| 5200-0301 | 5 holes L  | 3.2mm     | 10.2mm | 73.4mm  |
| 5200-0302 | 6 holes L  | 3.2mm     | 10.2mm | 87mm    |
| 5200-0303 | 7 holes L  | 3.2mm     | 10.2mm | 100.9mm |
| 5200-0304 | 8 holes L  | 3.2mm     | 10.2mm | 112.6mm |
| 5200-0305 | 9 holes L  | 3.2mm     | 10.2mm | 126.6mm |
| 5200-0306 | 10 holes L | 3.2mm     | 10.2mm | 140.4mm |
| 5200-0307 | 5 holes R  | 3.2mm     | 10.2mm | 73.4mm  |
| 5200-0308 | 6 holes R  | 3.2mm     | 10.2mm | 87mm    |
| 5200-0309 | 7 holes R  | 3.2mm     | 10.2mm | 100.9mm |
| 5200-0310 | 8 holes R  | 3.2mm     | 10.2mm | 112.6mm |
| 5200-0311 | 9 holes R  | 3.2mm     | 10.2mm | 126.6mm |
| 5200-0312 | 10 holes R | 3.2mm     | 10.2mm | 140.4mm |

Use 3.5 Locking Screw/3.5 Cortical Screw/4.0 Cancellous Screw

| REF       | Holes     | Thickness | Width | Length |
|-----------|-----------|-----------|-------|--------|
| 5100-0501 | 3 holes L | 3mm       | 10mm  | 55mm   |
| 5100-0502 | 4 holes L | 3mm       | 10mm  | 68mm   |
| 5100-0503 | 5 holes L | 3mm       | 10mm  | 81mm   |
| 5100-0504 | 6 holes L | 3mm       | 10mm  | 94mm   |
| 5100-0505 | 7 holes L | 3mm       | 10mm  | 107mm  |
| 5100-0506 | 3 holes R | 3mm       | 10mm  | 55mm   |
| 5100-0507 | 4 holes R | 3mm       | 10mm  | 68mm   |
| 5100-0508 | 5 holes R | 3mm       | 10mm  | 81mm   |
| 5100-0509 | 6 holes R | 3mm       | 10mm  | 94mm   |
| 5100-0510 | 7 holes R | 3mm       | 10mm  | 107mm  |

Use 2.7/3.5 Locking Screw/3.5 Cortical Screw/4.0 Cancellous Screw



Distal Clavicular Locking Plate-I



| REF       | Holes     | Thickness | Width | Length |
|-----------|-----------|-----------|-------|--------|
| 5100-8801 | 3holes L  | 3mm       | 10mm  | 67mm   |
| 5100-8802 | 4 holes L | 3mm       | 10mm  | 81mm   |
| 5100-8803 | 5 holes L | 3mm       | 10mm  | 95mm   |
| 5100-8804 | 6 holes L | 3mm       | 10mm  | 108mm  |
| 5100-8805 | 7 holes L | 3mm       | 10mm  | 122mm  |
| 5100-8806 | 8 holes L | 3mm       | 10mm  | 136mm  |
| 5100-8807 | 3 holes R | 3mm       | 10mm  | 67mm   |
| 5100-8808 | 4 holes R | 3mm       | 10mm  | 81mm   |
| 5100-8809 | 5 holes R | 3mm       | 10mm  | 95mm   |
| 5100-8810 | 6 holes R | 3mm       | 10mm  | 108mm  |
| 5100-8811 | 7 holes R | 3mm       | 10mm  | 122mm  |
| 5100-8812 | 8 holes R | 3mm       | 10mm  | 136mm  |

Use 2.7/3.5 Locking Screw/3.5 Cortical Screw/4.0 Cancellous Screw

Clavicle Hook Locking Plate



| REF       | Holes     | Thickness | Width  | Length |
|-----------|-----------|-----------|--------|--------|
| 5100-0601 | 2 holes L | 3mm       | 11.5mm | 89mm   |
| 5100-0602 | 3 holes L | 3mm       | 11.5mm | 103mm  |
| 5100-0603 | 4 holes L | 3mm       | 11.5mm | 119mm  |
| 5100-0604 | 5 holes L | 3mm       | 11.5mm | 134mm  |
| 5100-0605 | 6 holes L | 3mm       | 11.5mm | 148mm  |
| 5100-0606 | 2 holes R | 3mm       | 11.5mm | 89mm   |
| 5100-0607 | 3 holes R | 3mm       | 11.5mm | 103mm  |
| 5100-0608 | 4 holes R | 3mm       | 11.5mm | 119mm  |
| 5100-0609 | 5 holes R | 3mm       | 11.5mm | 134mm  |
| 5100-0610 | 6 holes R | 3mm       | 11.5mm | 148mm  |

Use 3.5 Locking Screw/3.5 Cortical Screw/4.0 Cancellous Screw

Distal Clavicular Locking Plate



| REF        | Holes     | Thickness | Width | Length |
|------------|-----------|-----------|-------|--------|
| 01.1158-01 | 3 holes L | 4.2mm     | 11mm  | 61mm   |
| 01.1158-02 | 4 holes L | 4.2mm     | 11mm  | 71mm   |
| 01.1158-03 | 5 holes L | 4.2mm     | 11mm  | 81mm   |
| 01.1158-11 | 3 holes R | 4.2mm     | 11mm  | 61mm   |
| 01.1158-12 | 4 holes R | 4.2mm     | 11mm  | 71mm   |
| 01.1158-13 | 5 holes R | 4.2mm     | 11mm  | 81mm   |

Use 2.7/3.5 Locking Screw/3.5 Cortical Screw

Olecranon Locking Plate



| REF       | Holes      | Thickness | Width | Length |
|-----------|------------|-----------|-------|--------|
| 5100-0701 | 3 holes L  | 2.5mm     | 11mm  | 107mm  |
| 5100-0702 | 4 holes L  | 2.5mm     | 11mm  | 120mm  |
| 5100-0703 | 6 holes L  | 2.5mm     | 11mm  | 146mm  |
| 5100-0704 | 8 holes L  | 2.5mm     | 11mm  | 172mm  |
| 5100-0705 | 10 holes L | 2.5mm     | 11mm  | 198mm  |
| 5100-0706 | 3 holes R  | 2.5mm     | 11mm  | 107mm  |
| 5100-0707 | 4 holes R  | 2.5mm     | 11mm  | 120mm  |
| 5100-0708 | 6 holes R  | 2.5mm     | 11mm  | 146mm  |
| 5100-0709 | 8 holes R  | 2.5mm     | 11mm  | 172mm  |
| 5100-0710 | 10 holes R | 2.5mm     | 11mm  | 198mm  |

Use 3.5 Locking Screw/3.5 Cortical Screw/4.0 Cancellous Screw

Olecranon Locking Plate II



| REF        | Holes   | Thickness | Width  | Length |
|------------|---------|-----------|--------|--------|
| 01.1062-01 | 3 holes | 3mm       | 10.2mm | 49mm   |
| 01.1062-02 | 4 holes | 3mm       | 10.2mm | 63mm   |
| 01.1062-03 | 5 holes | 3mm       | 10.2mm | 77mm   |
| 01.1062-04 | 6 holes | 3mm       | 10.2mm | 91mm   |
| 01.1062-05 | 7 holes | 3mm       | 10.2mm | 105mm  |
| 01.1062-06 | 8 holes | 3mm       | 10.2mm | 119mm  |

Use 3.5 Locking Screw/3.5 Cortical Screw

Distal Ulnar Locking Plate



| REF       | Holes    | Thickness | Width | Length |
|-----------|----------|-----------|-------|--------|
| 5100-0801 | 6 hoels  | 1.3mm     | 5.1mm | 39mm   |
| 5100-0802 | 8 hoels  | 1.3mm     | 5.1mm | 49mm   |
| 5100-0803 | 10 holes | 1.3mm     | 5.1mm | 59mm   |

Use 2.4 Locking Screw/2.4 Cortical Screw

Distal Ulnar Locking Plate



| REF       | Holes   | Thickness | Width | Length |
|-----------|---------|-----------|-------|--------|
| 5100-9501 | 3 holes | 1.8mm     | 7mm   | 46mm   |
| 5100-9502 | 4 holes | 1.8mm     | 7mm   | 57mm   |
| 5100-9503 | 5 holes | 1.8mm     | 7mm   | 68mm   |

(Use 2.7 Locking Screw/2.7 Cortical Screw)

VA Distal Lateral Radius Locking Plate



| REF       | Holes   | Thickness | Width | Length |
|-----------|---------|-----------|-------|--------|
| 5100-0901 | 5 holes | 2mm       | 6.7mm | 47mm   |
| 5100-0902 | 6 holes | 2mm       | 6.7mm | 55mm   |

(Use 2.4/2.7 Universal Locking Screw/2.7 Locking Screw/2.7 Cortical Screw)

VA Distal Medial Radius Locking Plate



| REF       | Holes     | Thickness | Width | Length |
|-----------|-----------|-----------|-------|--------|
| 5100-1001 | 4 holes L | 2mm       | 7.2mm | 41mm   |
| 5100-1002 | 5 holes L | 2mm       | 7.2mm | 48mm   |
| 5100-1003 | 6 holes L | 2mm       | 7.2mm | 55mm   |
| 5100-1004 | 4 holes R | 2mm       | 7.2mm | 41mm   |
| 5100-1005 | 5 holes R | 2mm       | 7.2mm | 48mm   |
| 5100-1006 | 6 holes R | 2mm       | 7.2mm | 55mm   |

(Use 2.4/2.7 Universal Locking Screw/2.7 Locking Screw/2.7 Cortical Screw)

Distal Dorsal Radial Delta Locking Plate



| REF       | Holes     | Thickness | Width | Length |
|-----------|-----------|-----------|-------|--------|
| 5100-1101 | 3 holes L | 1.8mm     | 8.6mm | 65mm   |
| 5100-1102 | 4 holes L | 1.8mm     | 8.6mm | 73mm   |
| 5100-1103 | 5 holes L | 1.8mm     | 8.6mm | 82mm   |
| 5100-1104 | 7 holes L | 1.8mm     | 8.6mm | 99mm   |
| 5100-1105 | 9 holes L | 1.8mm     | 8.6mm | 116mm  |
| 5100-1106 | 3 holes R | 1.8mm     | 8.6mm | 65mm   |
| 5100-1107 | 4 holes R | 1.8mm     | 8.6mm | 73mm   |
| 5100-1108 | 5 holes R | 1.8mm     | 8.6mm | 82mm   |
| 5100-1109 | 7 holes R | 1.8mm     | 8.6mm | 99mm   |
| 5100-1110 | 9 holes R | 1.8mm     | 8.6mm | 116mm  |

Use 2.7 Locking Screw/2.7 Cortical screw



VA Distal Radius Locking Plate



| REF       | Holes     | Thickness | Width | Length |
|-----------|-----------|-----------|-------|--------|
| 5100-1212 | 2 holes L | 2.2mm     | 8mm   | 52mm   |
| 5100-1201 | 3 holes L | 2.2mm     | 8mm   | 60mm   |
| 5100-1202 | 4 holes L | 2.2mm     | 8mm   | 69mm   |
| 5100-1203 | 5 holes L | 2.2mm     | 8mm   | 77mm   |
| 5100-1204 | 6 holes L | 2.2mm     | 8mm   | 86mm   |
| 5100-1205 | 7 holes L | 2.2mm     | 8mm   | 94mm   |
| 5100-1211 | 2 holes R | 2.2mm     | 8mm   | 52mm   |
| 5100-1206 | 3 holes R | 2.2mm     | 8mm   | 60mm   |
| 5100-1207 | 4 holes R | 2.2mm     | 8mm   | 69mm   |
| 5100-1208 | 5 holes R | 2.2mm     | 8mm   | 77mm   |
| 5100-1209 | 6 holes R | 2.2mm     | 8mm   | 86mm   |
| 5100-1210 | 7 holes R | 2.2mm     | 8mm   | 94mm   |

(Use 2.4/2.7 Universal Locking Screw/2.7 Locking Screw/2.7 Cortical Screw)

VA Distal Radius Locking Plate



| REF       | Holes     | Thickness | Width | Length |
|-----------|-----------|-----------|-------|--------|
| 5100-6601 | 2 holes L | 2.2       | 8mm   | 52mm   |
| 5100-6602 | 3 holes L | 2.2       | 8mm   | 60mm   |
| 5100-6603 | 4 holes L | 2.2       | 8mm   | 69mm   |
| 5100-6604 | 5 holes L | 2.2       | 8mm   | 77mm   |
| 5100-6605 | 6 holes L | 2.2       | 8mm   | 86mm   |
| 5100-6606 | 7 holes L | 2.2       | 8mm   | 94mm   |
| 5100-6607 | 2 holes R | 2.2       | 8mm   | 52mm   |
| 5100-6608 | 3 holes R | 2.2       | 8mm   | 60mm   |
| 5100-6609 | 4 holes R | 2.2       | 8mm   | 69mm   |
| 5100-6610 | 5 holes R | 2.2       | 8mm   | 77mm   |
| 5100-6611 | 6 holes R | 2.2       | 8mm   | 86mm   |
| 5100-6612 | 7 holes R | 2.2       | 8mm   | 94mm   |

(Use 2.4/2.7 Universal Locking Screw/2.7 Locking Screw/2.7 Cortical Screw)

Distal Volar Radial Locking Plate With Drill Guide



| REF       | Holes     | Thickness | Width | Length |
|-----------|-----------|-----------|-------|--------|
| 5100-1301 | 3 holes L | 2.5mm     | 9mm   | 49mm   |
| 5100-1302 | 4 holes L | 2.5mm     | 9mm   | 58mm   |
| 5100-1303 | 5 holes L | 2.5mm     | 9mm   | 66mm   |
| 5100-1304 | 7 holes L | 2.5mm     | 9mm   | 83mm   |
| 5100-1305 | 9 holes L | 2.5mm     | 9mm   | 99mm   |
| 5100-1306 | 3 holes R | 2.5mm     | 9mm   | 49mm   |
| 5100-1307 | 4 holes R | 2.5mm     | 9mm   | 58mm   |
| 5100-1308 | 5 holes R | 2.5mm     | 9mm   | 66mm   |
| 5100-1309 | 7 holes R | 2.5mm     | 9mm   | 83mm   |
| 5100-1310 | 9 holes R | 2.5mm     | 9mm   | 99mm   |

Use 2.7 Locking Screw/2.7 Cortical screw

Proximal Radius Locking Plate



| REF       | Holes     | Thickness | Width | Length |
|-----------|-----------|-----------|-------|--------|
| 5100-1401 | 3 holes L | 1.8mm     | 8.7mm | 53mm   |
| 5100-1402 | 4 holes L | 1.8mm     | 8.7mm | 63mm   |
| 5100-1403 | 5 holes L | 1.8mm     | 8.7mm | 72mm   |
| 5100-1404 | 3 holes R | 1.8mm     | 8.7mm | 53mm   |
| 5100-1405 | 4 holes R | 1.8mm     | 8.7mm | 63mm   |
| 5100-1406 | 5 holes R | 1.8mm     | 8.7mm | 72mm   |

Use 2.4 Locking Screw/2.4 Cortical Screw

Proximal Radius Locking Plate-II



| REF       | Holes     | Thickness | Width | Length |
|-----------|-----------|-----------|-------|--------|
| 5100-9201 | 2 holes L | 1.7mm     | 6.9mm | 35mm   |
| 5100-9202 | 3 holes L | 1.7mm     | 6.9mm | 45mm   |
| 5100-9203 | 4 holes L | 1.7mm     | 6.9mm | 54mm   |
| 5100-9204 | 2 holes R | 1.7mm     | 6.9mm | 35mm   |
| 5100-9205 | 3 holes R | 1.7mm     | 6.9mm | 45mm   |
| 5100-9206 | 4 holes R | 1.7mm     | 6.9mm | 54mm   |

Use 2.4 Locking Screw/2.4 Cortical Screw

Proximal Radius Locking Plate-III



| REF       | Holes     | Thickness | Width | Length |
|-----------|-----------|-----------|-------|--------|
| 5100-9301 | 2 holes L | 1.5mm     | 6.2mm | 26mm   |
| 5100-9302 | 2 holes R | 1.5mm     | 6.2mm | 26mm   |

Use 2.4 Locking Screw/2.4 Cortical Screw

Coronoid Locking Plate



| REF       | Holes     | Thickness | Width  | Length |
|-----------|-----------|-----------|--------|--------|
| 5200-0201 | 6 holes L | 1.4mm     | 22.6mm | 39mm   |
| 5200-0202 | 6 holes R | 1.4mm     | 22.6mm | 39mm   |

Use 2.4 Locking Screw

Coronoid Locking Plate-II



| REF       | Holes | Thickness | Width | Length |
|-----------|-------|-----------|-------|--------|
| 5100-9001 | L     | 1.4mm     | 7.3mm | 29.3mm |
| 5100-9002 | R     | 1.4mm     | 7.3mm | 29.3mm |

Use 2.0/2.4 Locking Screw

Coronoid Locking Plate-III



| REF       | Holes | Thickness | Width | Length |
|-----------|-------|-----------|-------|--------|
| 5100-9101 | L     | 1.5mm     | 8.3mm | 17.6mm |
| 5100-9102 | R     | 1.5mm     | 8.3mm | 17.6mm |

Use 2.0/2.4 Locking Screw



Distal Diaphyseal Radius Locking Plate-II



| REF        | Holes      | Thickness | Width | Length  |
|------------|------------|-----------|-------|---------|
| 5100-9401  | 6 holes L  | 2.5mm     | 10mm  | 95.2mm  |
| 5100-9402  | 8 holes L  | 2.5mm     | 10mm  | 121.9mm |
| 5100-9403  | 10 holes L | 2.5mm     | 10mm  | 148.7mm |
| 5100-9404  | 12 holes L | 2.5mm     | 10mm  | 175.3mm |
| 5100-9405  | 13 holes L | 2.5mm     | 10mm  | 188.5mm |
| 5100-9406  | 6 holes R  | 2.5mm     | 10mm  | 95.2mm  |
| 5100-9407  | 8 holes R  | 2.5mm     | 10mm  | 121.9mm |
| 5100-9408  | 10 holes R | 2.5mm     | 10mm  | 148.7mm |
| 5100-94019 | 12 holes R | 2.5mm     | 10mm  | 175.3mm |
| 5100-9410  | 13 holes R | 2.5mm     | 10mm  | 188.5mm |

Use 3.5 Locking Screw/3.5 Cortical Screw

Distal Dorsal Radial Delta Locking Plate



| REF        | Holes     | Thickness | Width | Length |
|------------|-----------|-----------|-------|--------|
| 01.1065-01 | 3 holes L | 2.5mm     | 8.5mm | 46mm   |
| 01.1065-02 | 4 holes L | 2.5mm     | 8.5mm | 56mm   |
| 01.1065-03 | 5 holes L | 2.5mm     | 8.5mm | 66mm   |
| 01.1065-04 | 6 holes L | 2.5mm     | 8.5mm | 76mm   |
| 01.1065-11 | 3 holes R | 2.5mm     | 8.5mm | 46mm   |
| 01.1065-12 | 4 holes R | 2.5mm     | 8.5mm | 56mm   |
| 01.1065-13 | 5 holes R | 2.5mm     | 8.5mm | 66mm   |
| 01.1065-14 | 6 holes R | 2.5mm     | 8.5mm | 76mm   |

Use 2.7 Locking Screw/2.7 Cortical Screw

Wrist Fusion Locking Plate



| REF        | Holes   | Thickness | Width  | Length |
|------------|---------|-----------|--------|--------|
| 01.1067-01 | 9 holes | 3.2mm     | 11.2mm | 117mm  |
| 01.1067-11 | 8 holes | 3.2mm     | 11.2mm | 112mm  |

Use 2.7/3.5 Locking Screw/2.7/3.5 Cortical Screw

Humeral Shaft Locking Plate



| REF       | Holes    | Thickness | Width | Length |
|-----------|----------|-----------|-------|--------|
| 5100-0101 | 6 holes  | 3.6mm     | 13mm  | 92mm   |
| 5100-0102 | 7 holes  | 3.6mm     | 13mm  | 105mm  |
| 5100-0103 | 8 holes  | 3.6mm     | 13mm  | 118mm  |
| 5100-0104 | 9 holes  | 3.6mm     | 13mm  | 131mm  |
| 5100-0105 | 10 holes | 3.6mm     | 13mm  | 144mm  |
| 5100-0106 | 12 holes | 3.6mm     | 13mm  | 170mm  |
| 5100-0107 | 14 holes | 3.6mm     | 13mm  | 196mm  |

Use 3.5 Locking Screw/3.5 Cortical Screw

Proximal Humeral Locking Plate



| REF       | Holes    | Thickness | Width | Length |
|-----------|----------|-----------|-------|--------|
| 5100-1501 | 3 holes  | 4mm       | 12mm  | 90mm   |
| 5100-1502 | 4 holes  | 4mm       | 12mm  | 102mm  |
| 5100-1503 | 5 holes  | 4mm       | 12mm  | 114mm  |
| 5100-1504 | 6 holes  | 4mm       | 12mm  | 126mm  |
| 5100-1505 | 7 holes  | 4mm       | 12mm  | 138mm  |
| 5100-1506 | 8 holes  | 4mm       | 12mm  | 150mm  |
| 5100-1507 | 10 holes | 4mm       | 12mm  | 174mm  |
| 5100-1508 | 12 holes | 4mm       | 12mm  | 198mm  |

Use 3.5 Locking Screw/3.5 Cortical Screw/4.0 Cancellous Screw

Proximal Humeral Locking Plate-III



| REF       | Holes      | Thickness | Width  | Length |
|-----------|------------|-----------|--------|--------|
| 5100-8901 | 2 holes L  | 4.2mm     | 12.5mm | 86mm   |
| 5100-8902 | 3 holes L  | 4.2mm     | 12.5mm | 99mm   |
| 5100-8903 | 4 holes L  | 4.2mm     | 12.5mm | 112mm  |
| 5100-8904 | 6 holes L  | 4.2mm     | 12.5mm | 138mm  |
| 5100-8905 | 8 holes L  | 4.2mm     | 12.5mm | 164mm  |
| 5100-8906 | 10 holes L | 4.2mm     | 12.5mm | 190mm  |
| 5100-8907 | 12 holes L | 4.2mm     | 12.5mm | 216mm  |
| 5100-8908 | 14holes L  | 4.2mm     | 12.5mm | 242mm  |
| 5100-899  | 2 holes R  | 4.2mm     | 12.5mm | 86mm   |
| 5100-8910 | 3 holes R  | 4.2mm     | 12.5mm | 99mm   |
| 5100-8911 | 4 holes R  | 4.2mm     | 12.5mm | 112mm  |
| 5100-8912 | 6 holes R  | 4.2mm     | 12.5mm | 138mm  |
| 5100-8913 | 8 holes R  | 4.2mm     | 12.5mm | 164mm  |
| 5100-8914 | 10 holes R | 4.2mm     | 12.5mm | 190mm  |
| 5100-8915 | 12 holes R | 4.2mm     | 12.5mm | 216mm  |
| 5100-8916 | 14holes R  | 4.2mm     | 12.5mm | 242mm  |

Use 3.5 Locking Screw/3.5 Cortical Screw

Proximal Humeral Greater Tuberosity Locking Plate



| REF       | Holes     | Thickness | Width | Length |
|-----------|-----------|-----------|-------|--------|
| 5100-1601 | 5 holes L | 1.5mm     | 13mm  | 44mm   |
| 5100-1602 | 5 holes R | 1.5mm     | 13mm  | 44mm   |

Use 2.7/3.5 Locking Screw/2.7/3.5 Cortical Screw

Distal Humeral Extra-articular Locking Plate



| REF        | Holes     | Thickness | Width  | Length |
|------------|-----------|-----------|--------|--------|
| 01.1107-01 | 4 holes L | 3.2mm     | 12.3mm | 114mm  |
| 01.1107-02 | 5 holes L | 3.2mm     | 12.3mm | 128mm  |
| 01.1107-03 | 6 holes L | 3.2mm     | 12.3mm | 142mm  |
| 01.1107-04 | 7 holes L | 3.2mm     | 12.3mm | 156mm  |
| 01.1107-05 | 8 holes L | 3.2mm     | 12.3mm | 170mm  |
| 01.1107-06 | 9 holes L | 3.2mm     | 12.3mm | 184mm  |
| 01.1107-11 | 4 holes R | 3.2mm     | 12.3mm | 114mm  |
| 01.1107-12 | 5 holes R | 3.2mm     | 12.3mm | 128mm  |
| 01.1107-13 | 6 holes R | 3.2mm     | 12.3mm | 142mm  |
| 01.1107-14 | 7 holes R | 3.2mm     | 12.3mm | 156mm  |
| 01.1107-15 | 8 holes R | 3.2mm     | 12.3mm | 170mm  |
| 01.1107-16 | 9 holes R | 3.2mm     | 12.3mm | 184mm  |

Use 3.5 Locking Screw/3.5 Cortical Screw

Distal Lateral Humeral Locking Plate



| REF       | Holes      | Thickness | Width | Length |
|-----------|------------|-----------|-------|--------|
| 5100-1701 | 4 holes L  | 3.2mm     | 12mm  | 86mm   |
| 5100-1702 | 6 holes L  | 3.2mm     | 12mm  | 112mm  |
| 5100-1703 | 8 holes L  | 3.2mm     | 12mm  | 138mm  |
| 5100-1704 | 10 holes L | 3.2mm     | 12mm  | 164mm  |
| 5100-1705 | 12 holes L | 3.2mm     | 12mm  | 190mm  |
| 5100-1706 | 4 holes R  | 3.2mm     | 12mm  | 86mm   |
| 5100-1707 | 6 holes R  | 3.2mm     | 12mm  | 112mm  |
| 5100-1708 | 8 holes R  | 3.2mm     | 12mm  | 138mm  |
| 5100-1709 | 10 holes R | 3.2mm     | 12mm  | 164mm  |
| 5100-1710 | 12 holes R | 3.2mm     | 12mm  | 190mm  |

Use 2.7/3.5 Locking Screw/3.5 Cortical Screw



Distal Lateral Humeral Locking Plate



| REF       | Holes      | Thickness | Width | Length |
|-----------|------------|-----------|-------|--------|
| 5100-0101 | 4 holes L  | 2.5mm     | 11mm  | 81mm   |
| 5100-0102 | 6 holes L  | 2.5mm     | 11mm  | 107mm  |
| 5100-0103 | 8 holes L  | 2.5mm     | 11mm  | 133mm  |
| 5100-0104 | 10 holes L | 2.5mm     | 11mm  | 159mm  |
| 5100-0105 | 12 holes L | 2.5mm     | 11mm  | 185mm  |
| 5100-0106 | 4 holes R  | 2.5mm     | 11mm  | 81mm   |
| 5100-0107 | 6 holes R  | 2.5mm     | 11mm  | 107mm  |
| 5100-0108 | 8 holes R  | 2.5mm     | 11mm  | 133mm  |
| 5100-0109 | 10 holes R | 2.5mm     | 11mm  | 159mm  |
| 5100-0110 | 12 holes R | 2.5mm     | 11mm  | 185mm  |

Use 3.5 Locking Screw/3.5 Cortical Screw

Distal Medial Humeral Locking Plate



| REF       | Holes      | Thickness | Width  | Length |
|-----------|------------|-----------|--------|--------|
| 5100-1801 | 4 holes L  | 3mm       | 11.5mm | 69mm   |
| 5100-1802 | 6 holes L  | 3mm       | 11.5mm | 95mm   |
| 5100-1803 | 8 holes L  | 3mm       | 11.5mm | 121mm  |
| 5100-1804 | 10 holes L | 3mm       | 11.5mm | 147mm  |
| 5100-1805 | 12 holes L | 3mm       | 11.5mm | 173mm  |
| 5100-1806 | 4 holes R  | 3mm       | 11.5mm | 69mm   |
| 5100-1807 | 6 holes R  | 3mm       | 11.5mm | 95mm   |
| 5100-1808 | 8 holes R  | 3mm       | 11.5mm | 121mm  |
| 5100-1809 | 10 holes R | 3mm       | 11.5mm | 147mm  |
| 5100-1810 | 12 holes R | 3mm       | 11.5mm | 173mm  |

Use 2.7/3.5 Locking Screw/3.5 Cortical Screw

Distal Fibular Locking Plate



| REF       | Holes     | Thickness | Width | Length |
|-----------|-----------|-----------|-------|--------|
| 5100-1901 | 3 holes L | 2.5mm     | 10mm  | 81mm   |
| 5100-1902 | 4 holes L | 2.5mm     | 10mm  | 93mm   |
| 5100-1903 | 5 holes L | 2.5mm     | 10mm  | 105mm  |
| 5100-1904 | 6 holes L | 2.5mm     | 10mm  | 117mm  |
| 5100-1905 | 7 holes L | 2.5mm     | 10mm  | 129mm  |
| 5100-1906 | 8 holes L | 2.5mm     | 10mm  | 141mm  |
| 5100-1907 | 3 holes R | 2.5mm     | 10mm  | 81mm   |
| 5100-1908 | 4 holes R | 2.5mm     | 10mm  | 93mm   |
| 5100-1909 | 5 holes R | 2.5mm     | 10mm  | 105mm  |
| 5100-1910 | 6 holes R | 2.5mm     | 10mm  | 117mm  |
| 5100-1911 | 7 holes R | 2.5mm     | 10mm  | 129mm  |
| 5100-1912 | 8 holes R | 2.5mm     | 10mm  | 141mm  |

Use 3.5 Locking Screw /3.5 Cortical Screw

Distal Fibular Locking Plate-II



| REF        | Holes      | Thickness | Width  | Length |
|------------|------------|-----------|--------|--------|
| 01.1222-01 | 3 holes L  | 2.5mm     | 10.4mm | 84mm   |
| 01.1222-02 | 4 holes L  | 2.5mm     | 10.4mm | 97mm   |
| 01.1222-03 | 5 holes L  | 2.5mm     | 10.4mm | 110mm  |
| 01.1222-04 | 6 holes L  | 2.5mm     | 10.4mm | 123mm  |
| 01.1222-05 | 7 holes L  | 2.5mm     | 10.4mm | 136mm  |
| 01.1222-06 | 9 holes L  | 2.5mm     | 10.4mm | 162mm  |
| 01.1222-07 | 11 holes L | 2.5mm     | 10.4mm | 188mm  |
| 01.1222-11 | 3 holes R  | 2.5mm     | 10.4mm | 84mm   |
| 01.1222-12 | 4 holes R  | 2.5mm     | 10.4mm | 97mm   |
| 01.1222-13 | 5 holes R  | 2.5mm     | 10.4mm | 110mm  |
| 01.1222-14 | 6 holes R  | 2.5mm     | 10.4mm | 123mm  |
| 01.1222-15 | 7 holes R  | 2.5mm     | 10.4mm | 136mm  |
| 01.1222-16 | 9 holes R  | 2.5mm     | 10.4mm | 162mm  |
| 01.1222-17 | 11 holes R | 2.5mm     | 10.4mm | 188mm  |

Use 2.7/3.5 Locking Screw/2.7/3.5 Cortical Screw

Distal Fibular Locking Plate-III



| REF       | Holes      | Thickness | Width | Length |
|-----------|------------|-----------|-------|--------|
| 5100-6901 | 3 holes L  | 2.2mm     | 10mm  | 72mm   |
| 5100-6902 | 4 holes L  | 2.2mm     | 10mm  | 85mm   |
| 5100-6903 | 5 holes L  | 2.2mm     | 10mm  | 98mm   |
| 5100-6904 | 7 holes L  | 2.2mm     | 10mm  | 124mm  |
| 5100-6905 | 9 holes L  | 2.2mm     | 10mm  | 150mm  |
| 5100-6906 | 11 holes L | 2.2mm     | 10mm  | 176mm  |
| 5100-6907 | 3 holes R  | 2.2mm     | 10mm  | 72mm   |
| 5100-6908 | 4 holes R  | 2.2mm     | 10mm  | 85mm   |
| 5100-6909 | 5 holes R  | 2.2mm     | 10mm  | 98mm   |
| 5100-6910 | 7 holes R  | 2.2mm     | 10mm  | 124mm  |
| 5100-6911 | 9 holes R  | 2.2mm     | 10mm  | 150mm  |
| 5100-6912 | 11 holes R | 2.2mm     | 10mm  | 176mm  |

Use 2.7/3.5 Locking Screw/3.5 Cortical Screw

Distal Fibular Locking Plate-IV



| REF       | Holes      | Thickness | Width | Length |
|-----------|------------|-----------|-------|--------|
| 5100-7001 | 3 holes L  | 2.2mm     | 11mm  | 82mm   |
| 5100-7002 | 4 holes L  | 2.2mm     | 11mm  | 94mm   |
| 5100-7003 | 5 holes L  | 2.2mm     | 11mm  | 106mm  |
| 5100-7004 | 6 holes L  | 2.2mm     | 11mm  | 118mm  |
| 5100-7005 | 7 holes L  | 2.2mm     | 11mm  | 130mm  |
| 5100-7006 | 8 holes L  | 2.2mm     | 11mm  | 142mm  |
| 5100-7007 | 10 holes L | 2.2mm     | 11mm  | 166mm  |
| 5100-7008 | 3 holes R  | 2.2mm     | 11mm  | 82mm   |
| 5100-7009 | 4 holes R  | 2.2mm     | 11mm  | 94mm   |
| 5100-7010 | 5 holes R  | 2.2mm     | 11mm  | 106mm  |
| 5100-7011 | 6 holes R  | 2.2mm     | 11mm  | 118mm  |
| 5100-7012 | 7 holes R  | 2.2mm     | 11mm  | 130mm  |
| 5100-7013 | 8 holes R  | 2.2mm     | 11mm  | 142mm  |
| 5100-7014 | 10 holes R | 2.2mm     | 11mm  | 166mm  |

Use 2.7/3.5 Locking Screw/3.5 Cortical Screw

Tibia Locking Plate Narrow



| REF       | Holes    | Thickness | Width | Length |
|-----------|----------|-----------|-------|--------|
| 5100-2201 | 7 holes  | 5mm       | 15mm  | 139mm  |
| 5100-2202 | 8 holes  | 5mm       | 15mm  | 157mm  |
| 5100-2203 | 9 holes  | 5mm       | 15mm  | 175mm  |
| 5100-2204 | 10 holes | 5mm       | 15mm  | 193mm  |
| 5100-2205 | 12 holes | 5mm       | 15mm  | 229mm  |
| 5100-2206 | 14 holes | 5mm       | 15mm  | 265mm  |

Use 5.0 Locking Screw/4.5 Cortical Screw



Proximal Medial Tibial Osteotomy Locking Plate



| REF                                      | Holes   | Thickness | Width | Length |
|------------------------------------------|---------|-----------|-------|--------|
| 5100-2301                                | 5 holes | 2.8mm     | 16mm  | 115mm  |
| Use 5.0 Locking Screw/4.5 Cortical Screw |         |           |       |        |

Proximal Lateral Tibial Osteotomy Locking Plate



| REF                                      | Holes     | Thickness | Width | Length |
|------------------------------------------|-----------|-----------|-------|--------|
| 01.1223-01                               | 3 holes L | 4.6mm     | 16mm  | 104mm  |
| 01.1223-11                               | 3 holes R | 4.6mm     | 16mm  | 104mm  |
| Use 5.0 Locking Screw/4.5 Cortical Screw |           |           |       |        |

Distal Medial Femoral Osteotomy Locking Plate



| REF                                      | Holes     | Thickness | Width  | Length |
|------------------------------------------|-----------|-----------|--------|--------|
| 01.1275-01                               | 3 holes L | 4.2mm     | 18.2mm | 121mm  |
| 01.1275-02                               | 4 holes L | 4.2mm     | 18.2mm | 101mm  |
| 01.1275-11                               | 3 holes R | 4.2mm     | 18.2mm | 121mm  |
| 01.1275-12                               | 4 hoels R | 4.2mm     | 18.2mm | 101mm  |
| Use 5.0 Locking Screw/4.5 Cortical Screw |           |           |        |        |

Distal Lateral Femoral Osteotomy Locking Plate



| REF                                      | Holes     | Thickness | Width | Length |
|------------------------------------------|-----------|-----------|-------|--------|
| 5100-5403                                | 5 holes L | 5.6mm     | 17mm  | 161mm  |
| 5100-5404                                | 5 holes R | 5.6mm     | 17mm  | 161mm  |
| Use 5.0 Locking Screw/4.5 Cortical Screw |           |           |       |        |

Lateral Tibial Head Buttress Locking Plate



| REF                                      | Holes      | Thickness | Width | Length |
|------------------------------------------|------------|-----------|-------|--------|
| 5100-2401                                | 5 holes L  | 4.6mm     | 15mm  | 144mm  |
| 5100-2402                                | 7 holes L  | 4.6mm     | 15mm  | 182mm  |
| 5100-2403                                | 9 holes L  | 4.6mm     | 15mm  | 220mm  |
| 5100-2404                                | 11 holes L | 4.6mm     | 15mm  | 258mm  |
| 5100-2405                                | 13 holes L | 4.6mm     | 15mm  | 296mm  |
| 5100-2406                                | 5 holes R  | 4.6mm     | 15mm  | 144mm  |
| 5100-2407                                | 7 holes R  | 4.6mm     | 15mm  | 182mm  |
| 5100-2408                                | 9 holes R  | 4.6mm     | 15mm  | 220mm  |
| 5100-2409                                | 11 holes R | 4.6mm     | 15mm  | 258mm  |
| 5100-2410                                | 13 holes R | 4.6mm     | 15mm  | 296mm  |
| Use 5.0 Locking Screw/4.5 Cortical Screw |            |           |       |        |

Proximal Lateral Tibial Locking Plate-I



| REF                                           | Holes      | Thickness | Width | Length |
|-----------------------------------------------|------------|-----------|-------|--------|
| 5100-2501                                     | 3 holes L  | 4.6mm     | 14mm  | 117mm  |
| 5100-2502                                     | 5 holes L  | 4.6mm     | 14mm  | 155mm  |
| 5100-2503                                     | 7 holes L  | 4.6mm     | 14mm  | 193mm  |
| 5100-2504                                     | 9 holes L  | 4.6mm     | 14mm  | 231mm  |
| 5100-2505                                     | 11 holes L | 4.6mm     | 14mm  | 269mm  |
| 5100-2506                                     | 3 holes R  | 4.6mm     | 14mm  | 117mm  |
| 5100-2507                                     | 5 holes R  | 4.6mm     | 14mm  | 155mm  |
| 5100-2508                                     | 7 holes R  | 4.6mm     | 14mm  | 193mm  |
| 5100-2509                                     | 9 holes R  | 4.6mm     | 14mm  | 231mm  |
| 5100-2510                                     | 11 holes R | 4.6mm     | 14mm  | 269mm  |
| Use 3.5/5.0 Locking Screw/ 3.5 Cortical Screw |            |           |       |        |

Proximal Lateral Tibial Locking Plate-II



| REF                                      | Holes      | Thickness | Width  | Length |
|------------------------------------------|------------|-----------|--------|--------|
| 5100-2601                                | 5 holes L  | 4.2mm     | 12.5mm | 97mm   |
| 5100-2602                                | 7 holes L  | 4.2mm     | 12.5mm | 124mm  |
| 5100-2603                                | 9 holes L  | 4.2mm     | 12.5mm | 151mm  |
| 5100-2604                                | 11 holes L | 4.2mm     | 12.5mm | 178mm  |
| 5100-2605                                | 13 holes L | 4.2mm     | 12.5mm | 205mm  |
| 5100-2606                                | 5 holes R  | 4.2mm     | 12.5mm | 97mm   |
| 5100-2607                                | 7 holes R  | 4.2mm     | 12.5mm | 124mm  |
| 5100-2608                                | 9 holes R  | 4.2mm     | 12.5mm | 151mm  |
| 5100-2609                                | 11 holes R | 4.2mm     | 12.5mm | 178mm  |
| 5100-2610                                | 13 holes R | 4.2mm     | 12.5mm | 205mm  |
| Use 3.5 Locking Screw/3.5 Cortical Screw |            |           |        |        |

Proximal Medial Tibial Locking Plate



| REF                                      | Holes      | Thickness | Width | Length |
|------------------------------------------|------------|-----------|-------|--------|
| 5100-2701                                | 4 holes L  | 4.2mm     | 13mm  | 83mm   |
| 5100-2702                                | 6 holes L  | 4.2mm     | 13mm  | 109mm  |
| 5100-2703                                | 8 holes L  | 4.2mm     | 13mm  | 135mm  |
| 5100-2704                                | 10 holes L | 4.2mm     | 13mm  | 161mm  |
| 5100-2705                                | 12 holes L | 4.2mm     | 13mm  | 187mm  |
| 5100-2706                                | 4 holes R  | 4.2mm     | 13mm  | 83mm   |
| 5100-2707                                | 6 holes R  | 4.2mm     | 13mm  | 109mm  |
| 5100-2708                                | 8 holes R  | 4.2mm     | 13mm  | 135mm  |
| 5100-2709                                | 10 holes R | 4.2mm     | 13mm  | 161mm  |
| 5100-2710                                | 12 holes R | 4.2mm     | 13mm  | 187mm  |
| Use 3.5 Locking Screw/3.5 Cortical Screw |            |           |       |        |

Distal Lateral Tibial Locking Plate-I



| REF                                      | Holes      | Thickness | Width  | Length |
|------------------------------------------|------------|-----------|--------|--------|
| 5100-2801                                | 5 holes L  | 3.6mm     | 16.5mm | 122mm  |
| 5100-2802                                | 7 holes L  | 3.6mm     | 16.5mm | 154mm  |
| 5100-2803                                | 9 holes L  | 3.6mm     | 16.5mm | 186mm  |
| 5100-2804                                | 11 holes L | 3.6mm     | 16.5mm | 218mm  |
| 5100-2805                                | 13 holes L | 3.6mm     | 16.5mm | 250mm  |
| 5100-2806                                | 5 holes R  | 3.6mm     | 16.5mm | 122mm  |
| 5100-2807                                | 7 holes R  | 3.6mm     | 16.5mm | 154mm  |
| 5100-2808                                | 9 holes R  | 3.6mm     | 16.5mm | 186mm  |
| 5100-2809                                | 11 holes R | 3.6mm     | 16.5mm | 218mm  |
| 5100-2810                                | 13 holes R | 3.6mm     | 16.5mm | 250mm  |
| Use 5.0 Locking Screw/4.5 Cortical Screw |            |           |        |        |



Distal Lateral Tibial Locking Plate-II



| REF       | Holes      | Thickness | Width | Length |
|-----------|------------|-----------|-------|--------|
| 5100-2901 | 5 holes L  | 4.2mm     | 13mm  | 136mm  |
| 5100-2902 | 7 holes L  | 4.2mm     | 13mm  | 168mm  |
| 5100-2903 | 9 holes L  | 4.2mm     | 13mm  | 200mm  |
| 5100-2904 | 11 holes L | 4.2mm     | 13mm  | 232mm  |
| 5100-2905 | 5 holes R  | 4.2mm     | 13mm  | 136mm  |
| 5100-2906 | 7 holes R  | 4.2mm     | 13mm  | 168mm  |
| 5100-2907 | 9 holes R  | 4.2mm     | 13mm  | 200mm  |
| 5100-2908 | 11 holes R | 4.2mm     | 13mm  | 232mm  |

Use 3.5 Locking Screw/3.5 Cortical Screw

Distal Medial Tibial Locking Plate-I



| REF       | Holes      | Thickness | Width | Length |
|-----------|------------|-----------|-------|--------|
| 5100-3001 | 5 holes L  | 4.2mm     | 14mm  | 147mm  |
| 5100-3002 | 7 holes L  | 4.2mm     | 14mm  | 179mm  |
| 5100-3003 | 9 holes L  | 4.2mm     | 14mm  | 211mm  |
| 5100-3004 | 11 holes L | 4.2mm     | 14mm  | 243mm  |
| 5100-3005 | 13 holes L | 4.2mm     | 14mm  | 275mm  |
| 5100-3006 | 5 holes R  | 4.2mm     | 14mm  | 147mm  |
| 5100-3007 | 7 holes R  | 4.2mm     | 14mm  | 179mm  |
| 5100-3008 | 9 holes R  | 4.2mm     | 14mm  | 211mm  |
| 5100-3009 | 11 holes R | 4.2mm     | 14mm  | 243mm  |
| 5100-3010 | 13 holes R | 4.2mm     | 14mm  | 275mm  |

Use 3.5 Locking Screw/3.5 Cortical Screw

Distal Medial Tibial Locking Plate-II



| REF       | Holes      | Thickness | Width  | Length |
|-----------|------------|-----------|--------|--------|
| 5100-3101 | 5 holes L  | 4.2mm     | 15.5mm | 125mm  |
| 5100-3102 | 7 holes L  | 4.2mm     | 15.5mm | 151mm  |
| 5100-3103 | 9 holes L  | 4.2mm     | 15.5mm | 177mm  |
| 5100-3104 | 11 holes L | 4.2mm     | 15.5mm | 203mm  |
| 5100-3105 | 5 holes R  | 4.2mm     | 15.5mm | 125mm  |
| 5100-3106 | 7 holes R  | 4.2mm     | 15.5mm | 151mm  |
| 5100-3107 | 9 holes R  | 4.2mm     | 15.5mm | 177mm  |
| 5100-3108 | 11 holes R | 4.2mm     | 15.5mm | 203mm  |

Use 3.5 Locking Screw/3.5 Cortical Screw

Distal Tibial T-Locking Plate



| REF       | Holes   | Thickness | Width  | Length  |
|-----------|---------|-----------|--------|---------|
| 5100-6303 | 3 holes | 3mm       | 12.3mm | 56.5mm  |
| 5100-6304 | 4 holes | 3mm       | 12.3mm | 73.5mm  |
| 5100-6305 | 6 holes | 3mm       | 12.3mm | 93.5mm  |
| 5100-6306 | 8 holes | 3mm       | 12.3mm | 113.5mm |

Use 3.5 Locking Screw/3.5 Cortical Screw

MD Distal Tibial Locking Plate



| REF       | Holes   | Thickness | Width  | Length |
|-----------|---------|-----------|--------|--------|
| 5100-6307 | 3 holes | 2mm       | 11.2mm | 53mm   |
| 5100-6308 | 4 holes | 2mm       | 11.2mm | 65mm   |

Use 3.5 Locking Screw/3.5 Cortical Screw

Distal Tibial Medial Locking Plate



| REF       | Holes   | Thickness | Width  | Length |
|-----------|---------|-----------|--------|--------|
| 5100-6301 | 4 holes | 3mm       | 10.3mm | 88mm   |
| 5100-6302 | 6 holes | 3mm       | 10.3mm | 112mm  |

Use 3.5 Locking Screw/3.5 Cortical Screw

Distal Tibial Medial Locking Plate



| REF        | Holes     | Thickness | Width | Length |
|------------|-----------|-----------|-------|--------|
| 01.1220-01 | 3 holes L | 2.1mm     | 10mm  | 70mm   |
| 01.1220-02 | 4 holes L | 2.1mm     | 10mm  | 70mm   |
| 01.1220-11 | 3 holes R | 2.1mm     | 10mm  | 93mm   |
| 01.1220-12 | 4 holes R | 2.1mm     | 10mm  | 93mm   |

Use 3.5 Locking Screw/3.5 Cortical Screw

Medium Anterolateral Ankle Fusion Plate



| REF       | Holes     | Thickness | Width | Length |
|-----------|-----------|-----------|-------|--------|
| 5100-8401 | 3 holes L | 3mm       | 13mm  | 60mm   |
| 5100-8402 | 3 holes R | 3mm       | 13mm  | 60mm   |

Use 3.5 Locking Screw/3.5 Cortical Screw

Anterolateral Ankle Fusion Locking Plate



| REF       | Holes     | Thickness | Width | Length |
|-----------|-----------|-----------|-------|--------|
| 5100-8501 | 4 holes L | 4mm       | 14mm  | 96mm   |
| 5100-8502 | 4 holes R | 4mm       | 14mm  | 96mm   |

Use 3.5 Locking Screw/3.5 Cortical Screw

Medium Anterolateral Ankle Fusion Plate



| REF        | Holes     | Thickness | Width  | Length |
|------------|-----------|-----------|--------|--------|
| 01.1224-01 | 4 holes L | 4mm       | 12.4mm | 60mm   |
| 01.1224-11 | 4 holes R | 4mm       | 12.4mm | 60mm   |

Use 3.5 Locking Screw/3.5 Cortical Screw

Medium Anterolateral Ankle Fusion Plate



| REF       | Holes   | Thickness | Width | Length |
|-----------|---------|-----------|-------|--------|
| 5100-8601 | 4 holes | 3mm       | 10mm  | 110mm  |

Use 3.5 Locking Screw/3.5 Cortical Screw



Curved Femur Locking Compression Plate  
Broad



| REF       | Holes    | Thickness | Width | Length |
|-----------|----------|-----------|-------|--------|
| 5100-3201 | 8 holes  | 6mm       | 18mm  | 158mm  |
| 5100-3202 | 9 holes  | 6mm       | 18mm  | 176mm  |
| 5100-3203 | 10 holes | 6mm       | 18mm  | 194mm  |
| 5100-3204 | 12 holes | 6mm       | 18mm  | 230mm  |
| 5100-3205 | 14 holes | 6mm       | 18mm  | 266mm  |
| 5100-3206 | 16 holes | 6mm       | 18mm  | 302mm  |

Use 5.0 Locking Screw/4.5 Cortical Screw

Straight Reconstruction Locking Plate



| REF       | Holes    | Thickness | Width | Length |
|-----------|----------|-----------|-------|--------|
| 5100-3301 | 5 holes  | 3.2mm     | 11mm  | 66mm   |
| 5100-3302 | 6 holes  | 3.2mm     | 11mm  | 79mm   |
| 5100-3303 | 7 holes  | 3.2mm     | 11mm  | 92mm   |
| 5100-3304 | 8 holes  | 3.2mm     | 11mm  | 105mm  |
| 5100-3305 | 9 holes  | 3.2mm     | 11mm  | 118mm  |
| 5100-3306 | 10 holes | 3.2mm     | 11mm  | 131mm  |
| 5100-3307 | 12 holes | 3.2mm     | 11mm  | 157mm  |

Use 3.5 Locking Screw/3.5 Cortical Screw

Anterior Sacroilia c Joint Plate



| REF       | Specification | Thickness | Width | Length |
|-----------|---------------|-----------|-------|--------|
| 5100-7101 | 4 holes L     | 3.1mm     | 10mm  | 78mm   |
| 5100-7102 | 4 holes R     | 3.1mm     | 10mm  | 78mm   |

Use 3.5 Locking Screw/3.5 Cortical Screw/4.0 Cancellous Screw

Anterior Pelvic Plate



| REF       | Specification | Thickness | Width | Length |
|-----------|---------------|-----------|-------|--------|
| 5100-7201 | 7 holes L     | 3.1mm     | 10mm  | 71mm   |
| 5100-7202 | 9 holes L     | 3.1mm     | 10mm  | 95mm   |
| 5100-7203 | 11 holes L    | 3.1mm     | 10mm  | 119mm  |
| 5100-7204 | 7 holes R     | 3.1mm     | 10mm  | 71mm   |
| 5100-7205 | 9 holes R     | 3.1mm     | 10mm  | 95mm   |
| 5100-7206 | 11 holes R    | 3.1mm     | 10mm  | 119mm  |

Use 3.5 Locking Screw/3.5 Cortical Screw/4.0 Cancellous Screw

Posterior Pelvic Plate



| REF       | Specification | Thickness | Width | Length |
|-----------|---------------|-----------|-------|--------|
| 5100-7301 | 6 holes L     | 3.1mm     | 10mm  | 72mm   |
| 5100-7302 | 8 holes L     | 3.1mm     | 10mm  | 96mm   |
| 5100-7303 | 10 holes L    | 3.1mm     | 10mm  | 120mm  |
| 5100-7304 | 6 holes R     | 3.1mm     | 10mm  | 72mm   |
| 5100-7305 | 8 holes R     | 3.1mm     | 10mm  | 96mm   |
| 5100-7306 | 10 holes R    | 3.1mm     | 10mm  | 120mm  |

Use 3.5 Locking Screw/3.5 Cortical Screw/4.0 Cancellous Screw


Anatomical Acetabular Plate



| REF       | Specification | Thickness | Width | Length |
|-----------|---------------|-----------|-------|--------|
| 5100-7401 | 6 holes L     | 3.1mm     | 10mm  | 63mm   |
| 5100-7402 | 8 holes L     | 3.1mm     | 10mm  | 87mm   |
| 5100-7403 | 10 holes L    | 3.1mm     | 10mm  | 111mm  |
| 5100-7404 | 6 holes R     | 3.1mm     | 10mm  | 63mm   |
| 5100-7405 | 8 holes R     | 3.1mm     | 10mm  | 87mm   |
| 5100-7406 | 10 holes R    | 3.1mm     | 10mm  | 111mm  |

Use 3.5 Locking Screw/3.5 Cortical Screw/4.0 Cancellous Screw

Patella Mesh Locking Plate



| REF       | Holes           | Thickness | Width | Length |
|-----------|-----------------|-----------|-------|--------|
| 5100-3401 | 16 holes Small  | 1mm       | 30mm  | 38mm   |
| 5100-3402 | 16 holes Medium | 1mm       | 33mm  | 42mm   |
| 5100-3403 | 16 holes Large  | 1mm       | 36mm  | 46mm   |

Use 2.7 Locking Screw

Distal Femoral Locking Plate



| REF       | Holes      | Thickness | Width  | Length |
|-----------|------------|-----------|--------|--------|
| 5100-3501 | 5 holes L  | 6mm       | 17.5mm | 164mm  |
| 5100-3502 | 7 holes L  | 6mm       | 17.5mm | 204mm  |
| 5100-3503 | 9 holes L  | 6mm       | 17.5mm | 244mm  |
| 5100-3504 | 11 holes L | 6mm       | 17.5mm | 284mm  |
| 5100-3505 | 13 holes L | 6mm       | 17.5mm | 324mm  |
| 5100-3506 | 5 holes R  | 6mm       | 17.5mm | 164mm  |
| 5100-3507 | 7 holes R  | 6mm       | 17.5mm | 204mm  |
| 5100-3508 | 9 holes R  | 6mm       | 17.5mm | 244mm  |
| 5100-3509 | 11 holes R | 6mm       | 17.5mm | 284mm  |
| 5100-3510 | 13 holes R | 6mm       | 17.5mm | 324mm  |

Use 5.0 Locking Screw/4.5 Cortical Screw

Proximal Femoral Locking Plate-I



| REF       | Holes      | Thickness | Width | Length |
|-----------|------------|-----------|-------|--------|
| 5100-3601 | 5 holes L  | 6mm       | 18mm  | 153mm  |
| 5100-3602 | 7 holes L  | 6mm       | 18mm  | 189mm  |
| 5100-3603 | 9 holes L  | 6mm       | 18mm  | 225mm  |
| 5100-3604 | 11 holes L | 6mm       | 18mm  | 261mm  |
| 5100-3605 | 13 holes L | 6mm       | 18mm  | 297mm  |
| 5100-3606 | 5 holes R  | 6mm       | 18mm  | 153mm  |
| 5100-3607 | 7 holes R  | 6mm       | 18mm  | 189mm  |
| 5100-3608 | 9 holes R  | 6mm       | 18mm  | 225mm  |
| 5100-3609 | 11 holes R | 6mm       | 18mm  | 261mm  |
| 5100-3610 | 13 holes R | 6mm       | 18mm  | 297mm  |

Use 5.0 Locking Screw/4.5 Cortical Screw



Proximal Femur Locking Plate-II



| REF       | Holes      | Thickness | Width | Length |
|-----------|------------|-----------|-------|--------|
| 5100-3701 | 3 holes L  | 6.5mm     | 19mm  | 133mm  |
| 5100-3702 | 4 holes L  | 6.5mm     | 19mm  | 151mm  |
| 5100-3703 | 5 holes L  | 6.5mm     | 19mm  | 169mm  |
| 5100-3704 | 7 holes L  | 6.5mm     | 19mm  | 205mm  |
| 5100-3705 | 9 holes L  | 6.5mm     | 19mm  | 241mm  |
| 5100-3706 | 11 holes L | 6.5mm     | 19mm  | 277mm  |
| 5100-3707 | 3 holes R  | 6.5mm     | 19mm  | 133mm  |
| 5100-3708 | 4 holes R  | 6.5mm     | 19mm  | 151mm  |
| 5100-3709 | 5 holes R  | 6.5mm     | 19mm  | 169mm  |
| 5100-3710 | 7 holes R  | 6.5mm     | 19mm  | 205mm  |
| 5100-3711 | 9 holes R  | 6.5mm     | 19mm  | 241mm  |
| 5100-3712 | 11 holes R | 6.5mm     | 19mm  | 277mm  |

Use 5.0 Locking Screw/4.5 Cortical Screw,6.5 Cannulated Full-threaded Locking Screw

Proximal Femoral Hook Locking Plate with Cable



| REF       | Holes      | Thickness | Width  | Length |
|-----------|------------|-----------|--------|--------|
| 5100-4701 | 4 holes L  | 6.2mm     | 17.5mm | 233mm  |
| 5100-4702 | 6 holes L  | 6.2mm     | 17.5mm | 258mm  |
| 5100-4703 | 8 holes L  | 6.2mm     | 17.5mm | 283mm  |
| 5100-4704 | 10 holes L | 6.2mm     | 17.5mm | 308mm  |
| 5100-4705 | 12 holes L | 6.2mm     | 17.5mm | 333mm  |
| 5100-4706 | 4 holes R  | 6.2mm     | 17.5mm | 233mm  |
| 5100-4707 | 6 holes R  | 6.2mm     | 17.5mm | 258mm  |
| 5100-4708 | 8 holes R  | 6.2mm     | 17.5mm | 283mm  |
| 5100-4709 | 10 holes R | 6.2mm     | 17.5mm | 308mm  |
| 5100-4710 | 12 holes R | 6.2mm     | 17.5mm | 333mm  |

Use 5.0 Locking Screw/ 4.5 Cortical Screw/1.7mm Cable

Proximal Femoral Hook Locking Plate



| REF       | Holes      | Thickness | Width  | Length  |
|-----------|------------|-----------|--------|---------|
| 5100-5301 | 3 holes L  | 5.6mm     | 17.7mm | 157.5mm |
| 5100-5302 | 4 holes L  | 5.6mm     | 17.7mm | 175.5mm |
| 5100-5303 | 5 holes L  | 5.6mm     | 17.7mm | 193.5mm |
| 5100-5304 | 7 holes L  | 5.6mm     | 17.7mm | 229.5mm |
| 5100-5305 | 9 holes L  | 5.6mm     | 17.7mm | 265.5mm |
| 5100-5306 | 11 holes L | 5.6mm     | 17.7mm | 301.5mm |
| 5100-5307 | 13 holes L | 5.6mm     | 17.7mm | 337.5mm |
| 5100-5308 | 3 holes R  | 5.6mm     | 17.7mm | 157.5mm |
| 5100-5309 | 4 holes R  | 5.6mm     | 17.7mm | 175.5mm |
| 5100-5310 | 5 holes R  | 5.6mm     | 17.7mm | 193.5mm |
| 5100-5311 | 7 holes R  | 5.6mm     | 17.7mm | 229.5mm |
| 5100-5312 | 9 holes R  | 5.6mm     | 17.7mm | 265.5mm |
| 5100-5313 | 11 holes R | 5.6mm     | 17.7mm | 301.5mm |
| 5100-5314 | 13 holes R | 5.6mm     | 17.7mm | 337.5mm |

Use 5.0 Locking Screw/ 4.5 Cortical Screw/ 6.5 Cannulaetd Locking Screw

Cable Positioning Pin (Stardrive)




| REF       | Holes |
|-----------|-------|
| 5100-4803 | 3.5mm |
| 5100-4804 | 5.0mm |

Titanium Cable



| REF       | Specification |
|-----------|---------------|
| 5100-4803 | 1.3mm         |
| 5100-4804 | 1.7mm         |

FNS Locking Plate (130° CCD angle) (Use 5.0 Locking Screw)



| REF       | Specification |
|-----------|---------------|
| 3300-0101 | 1 hole        |
| 3300-0102 | 2 holes       |

FNS Bolts+Antirotation Screws



| REF       | Specification | REF       | Specification |
|-----------|---------------|-----------|---------------|
| 3300-0201 | 75mm          | 3300-0206 | 100mm         |
| 3300-0202 | 80mm          | 3300-0207 | 105mm         |
| 3300-0203 | 85mm          | 3300-0208 | 110mm         |
| 3300-0204 | 90mm          | 3300-0209 | 115mm         |
| 3300-0205 | 95mm          |           |               |



Calcaneus Locking Plate-I, Curved



| REF       | Holes        | Thickness | Width  | Length |
|-----------|--------------|-----------|--------|--------|
| 5100-3801 | Small Right  | 2mm       | 34mm   | 60mm   |
| 5100-3802 | Small Left   | 2mm       | 34mm   | 60mm   |
| 5100-3803 | Medium Right | 2mm       | 34.5mm | 67mm   |
| 5100-3804 | Medium Left  | 2mm       | 34.5mm | 67mm   |
| 5100-3805 | Large Right  | 2mm       | 35mm   | 73mm   |
| 5100-3806 | Large Left   | 2mm       | 35mm   | 73mm   |

Use 3.5 Locking Screw

Calcaneus Locking Plate



| REF       | Holes               | Thickness | Width | Length |
|-----------|---------------------|-----------|-------|--------|
| 5100-4905 | Short 3 holes Left  | 1.5mm     | 8.5mm | 76.5mm |
| 5100-4906 | Short 3 holes Right | 1.5mm     | 8.5mm | 76.5mm |
| 5100-4907 | Long 3 holes Left   | 1.5mm     | 8.5mm | 76.5mm |
| 5100-4908 | Long 3 holes Right  | 1.5mm     | 8.5mm | 76.5mm |
| 5100-4909 | Short 4 holes Left  | 1.5mm     | 8.5mm | 76.5mm |
| 5100-4910 | Short 4 holes Right | 1.5mm     | 8.5mm | 76.5mm |
| 5100-4911 | Long 4 holes Left   | 1.5mm     | 8.5mm | 76.5mm |
| 5100-4912 | Long 4 holes Right  | 1.5mm     | 8.5mm | 76.5mm |

Use 3.5 Locking Screw

MD AP and PT Calcaneal Locking Plate



| REF        | Holes        | Thickness | Width | Length |
|------------|--------------|-----------|-------|--------|
| 01.4070-01 | Small Right  | 1.8mm     | 10mm  | 56.3mm |
| 01.4070-02 | Small Left   | 1.8mm     | 10mm  | 58.9mm |
| 01.4070-03 | Medium Right | 1.8mm     | 10mm  | 67.8mm |
| 01.4070-11 | Medium Left  | 1.8mm     | 10mm  | 56.3mm |
| 01.4070-12 | Large Right  | 1.8mm     | 10mm  | 58.9mm |
| 01.4070-13 | Large Left   | 1.8mm     | 10mm  | 67.8mm |

Use 3.5 Locking Screw

MD Posterior Tuberosity Calcaneal Locking Plate



| REF        | Holes       | Thickness | Width | Length |
|------------|-------------|-----------|-------|--------|
| 01.4071-01 | Small Right | 1.8mm     | 10mm  | 39.8mm |
| 01.4071-02 | Small Left  | 1.8mm     | 10mm  | 49.3mm |
| 01.4071-11 | Large Right | 1.8mm     | 10mm  | 39.8mm |
| 01.4071-12 | Large Left  | 1.8mm     | 10mm  | 49.3mm |

Use 3.5 Locking Screw

MD Anterior Process Calcaneal Locking Plate



| REF        | Holes       | Thickness | Width | Length |
|------------|-------------|-----------|-------|--------|
| 01.4072-01 | Small Right | 1.8mm     | 10mm  | 46mm   |
| 01.4072-02 | Small Left  | 1.8mm     | 10mm  | 50.8mm |
| 01.4072-11 | Large Right | 1.8mm     | 10mm  | 46mm   |
| 01.4072-12 | Large Left  | 1.8mm     | 10mm  | 50.8mm |

Use 3.5 Locking Screw

MT Fifth Metatarsal Locking Plate



| REF        | Holes     | Thickness | Width | Length |
|------------|-----------|-----------|-------|--------|
| 01.4073-01 | 6 holes L | 1.5mm     | 7mm   | 48mm   |
| 01.4073-11 | 6 holes R | 1.5mm     | 7mm   | 48mm   |

Use 2.7 Locking Screw/ 2.7 Cortical Screw

MT First Metatarsocuneiform Fusion Locking Plate



| REF        | Holes   | Thickness | Width | Length |
|------------|---------|-----------|-------|--------|
| 01.4074-01 | 8 holes | 1.5mm     | 7.7mm | 60mm   |

Use 2.7 Locking Screw/ 2.7 Cortical Screw

MT L-shape Metatarsal Locking Plate



| REF        | Holes     | Thickness | Width | Length |
|------------|-----------|-----------|-------|--------|
| 01.4075-01 | 4 holes L | 1.5mm     | 8.5mm | 45.1mm |
| 01.4075-02 | 5 holes L | 1.5mm     | 8.5mm | 51.6mm |
| 01.4075-11 | 4 holes R | 1.5mm     | 8.5mm | 45.1mm |
| 01.4075-12 | 5 holes R | 1.5mm     | 8.5mm | 51.6mm |

Use 3.5 Locking Screw/ 3.5 Cortical Screw

MT T-shape Locking Plate



| REF        | Holes   | Thickness | Width | Length |
|------------|---------|-----------|-------|--------|
| 01.4076-01 | 2 holes | 1.5mm     | 7.1mm | 27.4mm |
| 01.4076-02 | 3 holes | 1.5mm     | 7.1mm | 33.9mm |
| 01.4076-03 | 4 holes | 1.5mm     | 7.1mm | 40.4mm |
| 01.4076-04 | 5 holes | 1.5mm     | 7.1mm | 47mm   |
| 01.4076-05 | 6 holes | 1.5mm     | 7.1mm | 53.7mm |

Use 2.7 Locking Screw/ 2.7 Cortical Screw

MT Straight Fusion Locking Plate



| REF        | Holes   | Thickness | Width | Length |
|------------|---------|-----------|-------|--------|
| 01.4077-01 | 2 holes | 1.5mm     | 7.5mm | 18.5mm |
| 01.4077-02 | 4 holes | 1.5mm     | 7.5mm | 35.6mm |
| 01.4077-03 | 3 holes | 1.5mm     | 7.5mm | 30mm   |

Use 2.7 Locking Screw/ 2.7 Cortical Screw



### MT First Metatarsophalangeal Fusion Locking Plate



| REF        | Holes       | Thickness | Width  | Length |
|------------|-------------|-----------|--------|--------|
| 01.4078-01 | Small Right | 1.5mm     | 13.7mm | 55.3mm |
| 01.4078-02 | Small Left  | 1.5mm     | 13.7mm | 55.3mm |
| 01.4078-11 | Large Right | 1.5mm     | 13.7mm | 60mm   |
| 01.4078-12 | Large Left  | 1.5mm     | 13.7mm | 60mm   |

Use 2.7 Locking Screw/ 2.7 Cortical Screw

### MT X-Locking Plate



| REF        | Holes  | Thickness | Width | Length |
|------------|--------|-----------|-------|--------|
| 01.4079-01 | 24mm   | 1.5mm     | 19mm  | 24mm   |
| 01.4079-02 | 28.5mm | 1.5mm     | 19mm  | 28.5mm |

Use 2.7 Locking Screw

### MT X-Locking Plate



| REF       | Holes | Thickness | Width | Length |
|-----------|-------|-----------|-------|--------|
| 5100-9901 | L     | 1.5mm     | 15mm  | 28mm   |
| 5100-9902 | R     | 1.5mm     | 15mm  | 28mm   |

Use 2.7 Locking Screw

### MT Cuboid Locking Plate



| REF        | Holes     | Thickness | Width  | Length |
|------------|-----------|-----------|--------|--------|
| 01.4080-01 | 9 holes L | 1.5mm     | 27.9mm | 29.5mm |
| 01.4080-11 | 9 holes R | 1.5mm     | 27.9mm | 29.5mm |

Use 2.7 Locking Screw

### MT Astragalus Locking Plate



| REF        | Holes   | Thickness | Width | Length |
|------------|---------|-----------|-------|--------|
| 01.4081-01 | 5 holes | 1.5mm     | 6mm   | 43mm   |

Use 2.7 Locking Screw

### MT Navicular Locking Plate



| REF        | Holes     | Thickness | Width | Length |
|------------|-----------|-----------|-------|--------|
| 01.4082-01 | 9 holes L | 1.5mm     | 6.3mm | 46mm   |
| 01.4082-11 | 9 holes R | 1.5mm     | 6.3mm | 46mm   |

Use 2.7 Locking Screw

### MT Calcaneal H-Locking Plate



| REF       | Holes       | Thickness | Width | Length |
|-----------|-------------|-----------|-------|--------|
| 5100-9701 | 6mm Height  | 1.8mm     | 22mm  | 25mm   |
| 5100-9702 | 8mm Height  | 1.8mm     | 22mm  | 25mm   |
| 5100-9703 | 10mm Height | 1.8mm     | 22mm  | 25mm   |

Use 3.5 Locking Screw

### MT Calcaneal H-Locking Plate



| REF       | Holes       | Thickness | Width | Length |
|-----------|-------------|-----------|-------|--------|
| 5100-9801 | 0mm Height  | 2mm       | 17mm  | 29mm   |
| 5100-9802 | 8mm Height  | 2mm       | 17mm  | 29mm   |
| 5100-9803 | 10mm Height | 2mm       | 17mm  | 29mm   |
| 5100-9804 | 12mm Height | 2mm       | 17mm  | 29mm   |

Use 3.5 Locking Screw

### MT Butterfly-shaped Locking Plate



| REF       | Holes        | Thickness | Width | Length |
|-----------|--------------|-----------|-------|--------|
| 5100-9601 | 0mm Height L | 1.4mm     | 15mm  | 25.3mm |
| 5100-9602 | 3mm Height L | 1.4mm     | 15mm  | 25.3mm |
| 5100-9603 | 4mm Height L | 1.4mm     | 15mm  | 25.3mm |
| 5100-9604 | 5mm Height L | 1.4mm     | 15mm  | 25.3mm |
| 5100-9605 | 6mm Height L | 1.4mm     | 15mm  | 25.3mm |
| 5100-9606 | 7mm Height L | 1.4mm     | 15mm  | 25.3mm |
| 5100-9607 | 0mm Height R | 1.4mm     | 15mm  | 25.3mm |
| 5100-9608 | 3mm Height R | 1.4mm     | 15mm  | 25.3mm |
| 5100-9609 | 4mm Height R | 1.4mm     | 15mm  | 25.3mm |
| 5100-9610 | 5mm Height R | 1.4mm     | 15mm  | 25.3mm |
| 5100-9611 | 6mm Height R | 1.4mm     | 15mm  | 25.3mm |
| 5100-9612 | 7mm Height R | 1.4mm     | 15mm  | 25.3mm |

Use 2.7 Locking Screw



Subtalar Screw(Tarsal Tunnel Approach)



| REF        | Specification |
|------------|---------------|
| 3105355021 | Φ5.5×21       |
| 3105365021 | Φ6.5×21       |
| 3105373021 | Φ7.3×21       |
| 3105385021 | Φ8.5×21       |
| 3105395021 | Φ9.5×21       |

Subtalar Screw(Tarsal Sinus Approach)



| REF        | Specification |
|------------|---------------|
| 3152385015 | Φ8.5*15       |
| 3152310015 | Φ10*15        |
| 3152311015 | Φ11*15        |
| 3152311020 | Φ11*20        |
| 3152312015 | Φ12*15        |
| 3152312020 | Φ12*20        |

2.4mm Cortex Screw, T8 Stardrive, Self-tapping



| REF       | Specification | REF       | Specification |
|-----------|---------------|-----------|---------------|
| 030350006 | 2.4*6         | 030350020 | 2.4*20        |
| 030350008 | 2.4*8         | 030350022 | 2.4*22        |
| 030350010 | 2.4*10        | 030350024 | 2.4*24        |
| 030350012 | 2.4*12        | 030350026 | 2.4*26        |
| 030350014 | 2.4*14        | 030350028 | ∴             |
| 030350016 | 2.4*16        | 030350030 | ∴             |
| 030350018 | 2.4*18        |           |               |

2.4mm Locking Screw, T8 Stardrive, Self-tapping



| REF       | Specification | REF       | Spe |
|-----------|---------------|-----------|-----|
| 040240006 | 2.4*6         | 040240020 | ∴   |
| 040240008 | 2.4*8         | 040240022 | ∴   |
| 040240010 | 2.4*10        | 040240024 | ∴   |
| 040240012 | 2.4*12        | 040240026 | ∴   |
| 040240014 | 2.4*14        | 040240028 | ∴   |
| 040240016 | 2.4*16        | 040240030 | ∴   |
| 040240018 | 2.4*18        |           |     |

2.7mm Cortex Screw, T8 Stardrive Self-tapping



| REF       | Specification | REF       | Specification |
|-----------|---------------|-----------|---------------|
| 030390010 | 2.7*10        | 030390022 | 2.7*22        |
| 030390012 | 2.7*12        | 030390024 | 2.7*24        |
| 030390014 | 2.7*14        | 030390026 | 2.7*26        |
| 030390016 | 2.7*16        | 030390028 | 2.7*28        |
| 030390018 | 2.7*18        | 030390030 | 2.7*30        |
| 030390020 | 2.7*20        |           |               |

2.7mm Locking Screw, T8 Stardrive Self-tapping



| REF       | Specification | REF       | Specification |
|-----------|---------------|-----------|---------------|
| 030340010 | 2.7*10        | 030340026 | 2.7*26        |
| 030340012 | 2.7*12        | 030340028 | 2.7*28        |
| 030340014 | 2.7*14        | 030340030 | 2.7*30        |
| 030340016 | 2.7*16        | 030340032 | 2.7*32        |
| 030340018 | 2.7*18        | 030340034 | 2.7*34        |
| 030340020 | 2.7*20        | 030340036 | 2.7*36        |
| 030340022 | 2.7*22        | 030340038 | 2.7*38        |
| 030340024 | 2.7*24        | 030340040 | 2.7*40        |

2.4mm Universal Locking Screw, T8 Stardrive



| REF       | Specification | REF       | Specification |
|-----------|---------------|-----------|---------------|
| 5100-6701 | 2.4*8         | 5100-6707 | 2.4*20        |
| 5100-6702 | 2.4*10        | 5100-6708 | 2.4*22        |
| 5100-6703 | 2.4*12        | 5100-6709 | 2.4*24        |
| 5100-6704 | 2.4*14        | 5100-6710 | 2.4*26        |
| 5100-6705 | 2.4*16        | 5100-6711 | 2.4*28        |
| 5100-6706 | 2.4*18        | 5100-6712 | 2.4*30        |

2.7mm Universal Locking Screw, T8 Stardrive



| REF       | Specification | REF       | Specification |
|-----------|---------------|-----------|---------------|
| 5100-6801 | 2.7*10        | 5100-6808 | 2.7*24        |
| 5100-6802 | 2.7*12        | 5100-6809 | 2.7*26        |
| 5100-6803 | 2.7*14        | 5100-6810 | 2.7*28        |
| 5100-6804 | 2.7*16        | 5100-6811 | 2.7*30        |
| 5100-6805 | 2.7*18        | 5100-6812 | 2.7*32        |
| 5100-6806 | 2.7*20        | 5100-6813 | 2.7*34        |
| 5100-6807 | 2.7*22        | 5100-6814 | 2.7*36        |

3.5mm VA Locking Screw



| REF       | Specification | REF       | Specification |
|-----------|---------------|-----------|---------------|
| 5200-0801 | 3.5*12        | 5200-0817 | 3.5*44        |
| 5200-0802 | 3.5*14        | 5200-0818 | 3.5*46        |
| 5200-0803 | 3.5*16        | 5200-0819 | 3.5*48        |
| 5200-0804 | 3.5*18        | 5200-0820 | 3.5*50        |
| 5200-0805 | 3.5*20        | 5200-0821 | 3.5*52        |
| 5200-0806 | 3.5*22        | 5200-0822 | 3.5*54        |
| 5200-0807 | 3.5*24        | 5200-0823 | 3.5*56        |
| 5200-0808 | 3.5*26        | 5200-0824 | 3.5*58        |
| 5200-0809 | 3.5*28        | 5200-0825 | 3.5*60        |
| 5200-0810 | 3.5*30        | 5200-0826 | 3.5*65        |
| 5200-0811 | 3.5*32        | 5200-0827 | 3.5*70        |
| 5200-0812 | 3.5*34        | 5200-0828 | 3.5*75        |
| 5200-0813 | 3.5*36        | 5200-0829 | 3.5*80        |
| 5200-0814 | 3.5*38        | 5200-0830 | 3.5*85        |
| 5200-0815 | 3.5*40        | 5200-0831 | 3.5*90        |
| 5200-0816 | 3.5*42        | 5200-0832 | 3.5*95        |



3.5mm Locking Screw (Stardrive)



| REF       | Specification |
|-----------|---------------|
| 5100-3901 | 3.5*12        |
| 5100-3902 | 3.5*14        |
| 5100-3903 | 3.5*16        |
| 5100-3904 | 3.5*18        |
| 5100-3905 | 3.5*20        |
| 5100-3906 | 3.5*22        |
| 5100-3907 | 3.5*24        |
| 5100-3908 | 3.5*26        |
| 5100-3909 | 3.5*28        |
| 5100-3910 | 3.5*30        |
| 5100-3911 | 3.5*32        |
| 5100-3912 | 3.5*34        |
| 5100-3913 | 3.5*36        |
| 5100-3914 | 3.5*38        |
| 5100-3915 | 3.5*40        |
| 5100-3916 | 3.5*42        |

| REF       | Specification |
|-----------|---------------|
| 5100-3917 | 3.5*44        |
| 5100-3918 | 3.5*46        |
| 5100-3919 | 3.5*48        |
| 5100-3920 | 3.5*50        |
| 5100-3921 | 3.5*52        |
| 5100-3922 | 3.5*54        |
| 5100-3923 | 3.5*56        |
| 5100-3924 | 3.5*58        |
| 5100-3925 | 3.5*60        |
| 5100-3926 | 3.5*65        |
| 5100-3927 | 3.5*70        |
| 5100-3928 | 3.5*75        |
| 5100-3929 | 3.5*80        |
| 5100-3930 | 3.5*85        |
| 5100-3931 | 3.5*90        |
| 5100-3932 | 3.5*95        |

5.0mm Locking Screw (Stardrive)



| REF       | Specification |
|-----------|---------------|
| 5100-4001 | 5.0*22        |
| 5100-4002 | 5.0*24        |
| 5100-4003 | 5.0*26        |
| 5100-4004 | 5.0*28        |
| 5100-4005 | 5.0*30        |
| 5100-4006 | 5.0*32        |
| 5100-4007 | 5.0*34        |
| 5100-4008 | 5.0*36        |
| 5100-4009 | 5.0*38        |
| 5100-4010 | 5.0*40        |
| 5100-4011 | 5.0*42        |
| 5100-4012 | 5.0*44        |
| 5100-4013 | 5.0*46        |
| 5100-4014 | 5.0*48        |

| REF       | Specification |
|-----------|---------------|
| 5100-4015 | 5.0*50        |
| 5100-4016 | 5.0*52        |
| 5100-4017 | 5.0*54        |
| 5100-4018 | 5.0*56        |
| 5100-4019 | 5.0*58        |
| 5100-4020 | 5.0*60        |
| 5100-4021 | 5.0*65        |
| 5100-4022 | 5.0*70        |
| 5100-4023 | 5.0*75        |
| 5100-4024 | 5.0*80        |
| 5100-4025 | 5.0*85        |
| 5100-4026 | 5.0*90        |
| 5100-4027 | 5.0*95        |

3.5mm Cortical Screw (Stardrive)



| REF       | Specification |
|-----------|---------------|
| 5100-4101 | 3.5*12        |
| 5100-4102 | 3.5*14        |
| 5100-4103 | 3.5*16        |
| 5100-4104 | 3.5*18        |
| 5100-4105 | 3.5*20        |
| 5100-4106 | 3.5*22        |
| 5100-4107 | 3.5*24        |
| 5100-4108 | 3.5*26        |
| 5100-4109 | 3.5*28        |
| 5100-4110 | 3.5*30        |
| 5100-4111 | 3.5*32        |

| REF       | Specification |
|-----------|---------------|
| 5100-4112 | 3.5*34        |
| 5100-4113 | 3.5*36        |
| 5100-4114 | 3.5*38        |
| 5100-4115 | 3.5*40        |
| 5100-4116 | 3.5*42        |
| 5100-4117 | 3.5*44        |
| 5100-4118 | 3.5*46        |
| 5100-4119 | 3.5*48        |
| 5100-4120 | 3.5*50        |
| 5100-4121 | 3.5*55        |
| 5100-4122 | 3.5*60        |

4.5mm Cortical Screw (Stardrive)



| REF       | Specification |
|-----------|---------------|
| 5100-4201 | 4.5*22        |
| 5100-4202 | 4.5*24        |
| 5100-4203 | 4.5*26        |
| 5100-4204 | 4.5*28        |
| 5100-4205 | 4.5*30        |
| 5100-4206 | 4.5*32        |
| 5100-4207 | 4.5*34        |
| 5100-4208 | 4.5*36        |
| 5100-4209 | 4.5*38        |
| 5100-4210 | 4.5*40        |

| REF       | Specification |
|-----------|---------------|
| 5100-4211 | 4.5*42        |
| 5100-4212 | 4.5*44        |
| 5100-4213 | 4.5*46        |
| 5100-4214 | 4.5*48        |
| 5100-4215 | 4.5*50        |
| 5100-4216 | 4.5*52        |
| 5100-4217 | 4.5*54        |
| 5100-4218 | 4.5*56        |
| 5100-4219 | 4.5*58        |
| 5100-4220 | 4.5*60        |

4.0mm Cancellous Screw (Stardrive)



| REF       | Specification |
|-----------|---------------|
| 5100-4301 | 4.0*12        |
| 5100-4302 | 4.0*14        |
| 5100-4303 | 4.0*16        |
| 5100-4304 | 4.0*18        |
| 5100-4305 | 4.0*20        |
| 5100-4306 | 4.0*22        |
| 5100-4307 | 4.0*24        |
| 5100-4308 | 4.0*26        |
| 5100-4309 | 4.0*28        |
| 5100-4310 | 4.0*30        |

| REF       | Specification |
|-----------|---------------|
| 5100-4311 | 4.0*32        |
| 5100-4312 | 4.0*34        |
| 5100-4313 | 4.0*36        |
| 5100-4314 | 4.0*38        |
| 5100-4315 | 4.0*40        |
| 5100-4316 | 4.0*42        |
| 5100-4317 | 4.0*44        |
| 5100-4318 | 4.0*46        |
| 5100-4319 | 4.0*48        |
| 5100-4320 | 4.0*50        |

6.5 Cannulated Full-threaded Locking Screw (Stardrive)



| REF       | Specification |
|-----------|---------------|
| 5100-4401 | 6.5*50        |
| 5100-4402 | 6.5*55        |
| 5100-4403 | 6.5*60        |
| 5100-4404 | 6.5*65        |
| 5100-4405 | 6.5*70        |
| 5100-4406 | 6.5*75        |
| 5100-4407 | 6.5*80        |

| REF       | Specification |
|-----------|---------------|
| 5100-4408 | 6.5*85        |
| 5100-4409 | 6.5*90        |
| 5100-4410 | 6.5*95        |
| 5100-4411 | 6.5*100       |
| 5100-4412 | 6.5*105       |
| 5100-4413 | 6.5*110       |

Rib Reconstruction Locking Plate



| REF       | Specification |
|-----------|---------------|
| 5100-2001 | 15 holes L    |
| 5100-2002 | 18 holes L    |
| 5100-2003 | 15 holes R    |
| 5100-2004 | 18 holes R    |

Rib Plate Locking Screw(Self-tapping)



| REF       | Specification |
|-----------|---------------|
| 5100-2101 | 2.9*8mm       |
| 5100-2102 | 2.9*10mm      |
| 5100-2103 | 2.9*12mm      |
| 5100-2104 | 2.9*14mm      |

| REF       | Specification |
|-----------|---------------|
| 5100-2105 | 2.9*16mm      |
| 5100-2106 | 2.9*18mm      |
| 5100-2107 | 2.9*20mm      |



Sternum Locking Plate-I



| REF       | Holes   | Thickness | Length |
|-----------|---------|-----------|--------|
| 5100-4501 | 4 holes | 1.7mm     | 20mm   |

Sternum Locking Plate-II



| REF       | Holes   | Thickness | Length |
|-----------|---------|-----------|--------|
| 5100-4502 | 6 holes | 1.7mm     | 29mm   |

Sternum Locking Plate-III



| REF       | Holes   | Thickness | Length |
|-----------|---------|-----------|--------|
| 5100-4503 | 8 holes | 1.7mm     | 20mm   |

Sternum Locking Plate-IV



| REF       | Holes           | Thickness | Length |
|-----------|-----------------|-----------|--------|
| 5100-4504 | 10 holes Narrow | 1.7mm     | 24mm   |

Sternum Locking Plate-V



| REF       | Holes         | Thickness | Length |
|-----------|---------------|-----------|--------|
| 5100-4505 | 10 holes Wild | 1.7mm     | 33mm   |

Straight Sternum Locking Plate-I



| REF       | Holes    | Thickness | Length |
|-----------|----------|-----------|--------|
| 5100-4506 | 12 holes | 1.7mm     | 105mm  |



2.7mm Universal Locking Screw, T8 Stardrive



3.5mm Locking Screw (Stardrive)



5.0mm Locking Screw (Stardrive)



3.5mm Cortical Screw (Stardrive)



| REF       | Specification |
|-----------|---------------|
| 5100-6801 | 2.7*10        |
| 5100-6802 | 2.7*12        |
| 5100-6803 | 2.7*14        |
| 5100-6804 | 2.7*16        |
| 5100-6805 | 2.7*18        |
| 5100-6806 | 2.7*20        |
| 5100-6807 | 2.7*22        |

| REF       | Specification |
|-----------|---------------|
| 5100-6808 | 2.7*24        |
| 5100-6809 | 2.7*26        |
| 5100-6810 | 2.7*28        |
| 5100-6811 | 2.7*30        |
| 5100-6812 | 2.7*32        |
| 5100-6813 | 2.7*34        |
| 5100-6814 | 2.7*36        |

| REF       | Specification |
|-----------|---------------|
| 5100-3901 | 3.5*12        |
| 5100-3902 | 3.5*14        |
| 5100-3903 | 3.5*16        |
| 5100-3904 | 3.5*18        |
| 5100-3905 | 3.5*20        |
| 5100-3906 | 3.5*22        |
| 5100-3907 | 3.5*24        |
| 5100-3908 | 3.5*26        |
| 5100-3909 | 3.5*28        |
| 5100-3910 | 3.5*30        |
| 5100-3911 | 3.5*32        |
| 5100-3912 | 3.5*34        |
| 5100-3913 | 3.5*36        |
| 5100-3914 | 3.5*38        |
| 5100-3915 | 3.5*40        |
| 5100-3916 | 3.5*42        |

| REF       | Specification |
|-----------|---------------|
| 5100-3917 | 3.5*44        |
| 5100-3918 | 3.5*46        |
| 5100-3919 | 3.5*48        |
| 5100-3920 | 3.5*50        |
| 5100-3921 | 3.5*52        |
| 5100-3922 | 3.5*54        |
| 5100-3923 | 3.5*56        |
| 5100-3924 | 3.5*58        |
| 5100-3925 | 3.5*60        |
| 5100-3926 | 3.5*65        |
| 5100-3927 | 3.5*70        |
| 5100-3928 | 3.5*75        |
| 5100-3929 | 3.5*80        |
| 5100-3930 | 3.5*85        |
| 5100-3931 | 3.5*90        |
| 5100-3932 | 3.5*95        |

| REF       | Specification |
|-----------|---------------|
| 5100-4001 | 5.0*22        |
| 5100-4002 | 5.0*24        |
| 5100-4003 | 5.0*26        |
| 5100-4004 | 5.0*28        |
| 5100-4005 | 5.0*30        |
| 5100-4006 | 5.0*32        |
| 5100-4007 | 5.0*34        |
| 5100-4008 | 5.0*36        |
| 5100-4009 | 5.0*38        |
| 5100-4010 | 5.0*40        |
| 5100-4011 | 5.0*42        |
| 5100-4012 | 5.0*44        |
| 5100-4013 | 5.0*46        |
| 5100-4014 | 5.0*48        |

| REF       | Specification |
|-----------|---------------|
| 5100-4015 | 5.0*50        |
| 5100-4016 | 5.0*52        |
| 5100-4017 | 5.0*54        |
| 5100-4018 | 5.0*56        |
| 5100-4019 | 5.0*58        |
| 5100-4020 | 5.0*60        |
| 5100-4021 | 5.0*65        |
| 5100-4022 | 5.0*70        |
| 5100-4023 | 5.0*75        |
| 5100-4024 | 5.0*80        |
| 5100-4025 | 5.0*85        |
| 5100-4026 | 5.0*90        |
| 5100-4027 | 5.0*95        |

| REF       | Specification |
|-----------|---------------|
| 5100-4101 | 3.5*12        |
| 5100-4102 | 3.5*14        |
| 5100-4103 | 3.5*16        |
| 5100-4104 | 3.5*18        |
| 5100-4105 | 3.5*20        |
| 5100-4106 | 3.5*22        |
| 5100-4107 | 3.5*24        |
| 5100-4108 | 3.5*26        |
| 5100-4109 | 3.5*28        |
| 5100-4110 | 3.5*30        |
| 5100-4111 | 3.5*32        |

| REF       | Specification |
|-----------|---------------|
| 5100-4112 | 3.5*34        |
| 5100-4113 | 3.5*36        |
| 5100-4114 | 3.5*38        |
| 5100-4115 | 3.5*40        |
| 5100-4116 | 3.5*42        |
| 5100-4117 | 3.5*44        |
| 5100-4118 | 3.5*46        |
| 5100-4119 | 3.5*48        |
| 5100-4120 | 3.5*50        |
| 5100-4121 | 3.5*55        |
| 5100-4122 | 3.5*60        |



#### 4.5mm Cortical Screw (Stardrive)



| REF       | Specification | REF       | Specification |
|-----------|---------------|-----------|---------------|
| 5100-4201 | 4.5*22        | 5100-4211 | 4.5*42        |
| 5100-4202 | 4.5*24        | 5100-4212 | 4.5*44        |
| 5100-4203 | 4.5*26        | 5100-4213 | 4.5*46        |
| 5100-4204 | 4.5*28        | 5100-4214 | 4.5*48        |
| 5100-4205 | 4.5*30        | 5100-4215 | 4.5*50        |
| 5100-4206 | 4.5*32        | 5100-4216 | 4.5*52        |
| 5100-4207 | 4.5*34        | 5100-4217 | 4.5*54        |
| 5100-4208 | 4.5*36        | 5100-4218 | 4.5*56        |
| 5100-4209 | 4.5*38        | 5100-4219 | 4.5*58        |
| 5100-4210 | 4.5*40        | 5100-4220 | 4.5*60        |

#### 4.0mm Cancellous Screw (Stardrive)



| REF       | Specification | REF       | Specification |
|-----------|---------------|-----------|---------------|
| 5100-4301 | 4.0*12        | 5100-4311 | 4.0*32        |
| 5100-4302 | 4.0*14        | 5100-4312 | 4.0*34        |
| 5100-4303 | 4.0*16        | 5100-4313 | 4.0*36        |
| 5100-4304 | 4.0*18        | 5100-4314 | 4.0*38        |
| 5100-4305 | 4.0*20        | 5100-4315 | 4.0*40        |
| 5100-4306 | 4.0*22        | 5100-4316 | 4.0*42        |
| 5100-4307 | 4.0*24        | 5100-4317 | 4.0*44        |
| 5100-4308 | 4.0*26        | 5100-4318 | 4.0*46        |
| 5100-4309 | 4.0*28        | 5100-4319 | 4.0*48        |
| 5100-4310 | 4.0*30        | 5100-4320 | 4.0*50        |

#### 6.5 Cannulated Full-threaded Locking Screw (Stardrive)



| REF       | Specification | REF       | Specification |
|-----------|---------------|-----------|---------------|
| 5100-4401 | 6.5*50        | 5100-4408 | 6.5*85        |
| 5100-4402 | 6.5*55        | 5100-4409 | 6.5*90        |
| 5100-4403 | 6.5*60        | 5100-4410 | 6.5*95        |
| 5100-4404 | 6.5*65        | 5100-4411 | 6.5*100       |
| 5100-4405 | 6.5*70        | 5100-4412 | 6.5*105       |
| 5100-4406 | 6.5*75        | 5100-4413 | 6.5*110       |
| 5100-4407 | 6.5*80        |           |               |

#### Rib Reconstruction Locking Plate



| REF       | Specification |
|-----------|---------------|
| 5100-2001 | 15 holes L    |
| 5100-2002 | 18 holes L    |
| 5100-2003 | 15 holes R    |
| 5100-2004 | 18 holes R    |

#### Rib Plate Locking Screw(Self-tapping)



| REF       | Specification | REF       | Specification |
|-----------|---------------|-----------|---------------|
| 5100-2101 | 2.9*8mm       | 5100-2105 | 2.9*16mm      |
| 5100-2102 | 2.9*10mm      | 5100-2106 | 2.9*18mm      |
| 5100-2103 | 2.9*12mm      | 5100-2107 | 2.9*20mm      |
| 5100-2104 | 2.9*14mm      |           |               |

#### Straight Sternum Locking Plate-II



| REF       | Holes    | Thickness | Length |
|-----------|----------|-----------|--------|
| 5100-4507 | 20 holes | 1.7mm     | 169mm  |

#### Sternum Locking Plate with Clip-I



| REF       | Thickness | Length |
|-----------|-----------|--------|
| 5100-4508 | 2.3mm     | 20mm   |

#### Sternum Locking Plate with Clip-II



| REF       | Thickness | Length |
|-----------|-----------|--------|
| 5100-4509 | 2.3mm     | 29mm   |

#### Sternum Locking Plate with Clip-III



| REF       | Thickness | Length |
|-----------|-----------|--------|
| 5100-4510 | 2.3mm     | 20mm   |

#### Sternum Locking Plate with Clip-IV



| REF       | Thickness | Length |
|-----------|-----------|--------|
| 5100-4511 | 2.3mm     | 24mm   |

#### Sternum Locking Plate with Clip-V



| REF       | Thickness | Length |
|-----------|-----------|--------|
| 5100-4512 | 2.3mm     | 33mm   |



Straight Sternum Locking Plate with Clip-I



| REF       | Thickness | Length |
|-----------|-----------|--------|
| 5100-4513 | 2.3mm     | 105mm  |

Straight Sternum Locking Plate with Clip-II



| REF       | Thickness | Length |
|-----------|-----------|--------|
| 5100-4514 | 2.3mm     | 169mm  |

Sternum Locking Screw



| REF       | Specification |
|-----------|---------------|
| 5100-4515 | 3.0*6mm       |
| 5100-4516 | 3.0*8mm       |
| 5100-4517 | 3.0*10mm      |
| 5100-4518 | 3.0*12mm      |
| 5100-4519 | 3.0*14mm      |
| 5100-4520 | 3.0*16mm      |

Pectus Excavatum Plate



| REF        | Specification | REF        | Specification |
|------------|---------------|------------|---------------|
| T4100-4101 | 12*170mm      | T4100-4106 | 12*230mm      |
| T4100-4102 | 12*180mm      | T4100-4107 | 12*240mm      |
| T4100-4103 | 12*190mm      | T4100-4108 | 12*255mm      |
| T4100-4104 | 12*205mm      | T4100-4109 | 12*265mm      |
| T4100-4105 | 12*215mm      |            |               |

Fixing Plate



| REF              | Specification |
|------------------|---------------|
| 17.10.151.150501 | 15*50mm       |

Plate Bender



| REF          | Specification     |
|--------------|-------------------|
| 50.36.001.00 | Large 300*65*32mm |

Gripper





| REF          | Specification        |
|--------------|----------------------|
| 50.36.002.00 | Straight 250*50*30mm |

Guider



| REF          | Specification |
|--------------|---------------|
| 50.36.003.00 | 480*60*35mm   |



Pediatric Hinge Locking Plate  
Instrument Set 3200-18



| NO. | REF       | Description                        | Qty. |
|-----|-----------|------------------------------------|------|
| 1   | 3200-1801 | Screwdriver SW2.5                  | 1    |
| 2   | 3200-1802 | Cannulated Screwdriver SW2.5       | 1    |
| 3   | 3200-1803 | K-wire 1.5*150                     | 1    |
| 4   | 3200-1804 | K-wire 1.5*150                     | 1    |
| 5   | 3200-1805 | K-wire 1.5*150                     | 1    |
| 6   | 3200-1806 | K-wire 1.5*150                     | 1    |
| 7   | 3200-1807 | Target Mould Houlder               | 1    |
| 8   | 3200-1808 | Target Mould Houlder               | 1    |
| 9   | 3200-1809 | K-wire Depth Gague 0-60mm          | 1    |
| 10  | 3200-1810 | Screw Holding Forcep               | 1    |
| 11  | 3200-1811 | Hign Plate Holding Forcep          | 1    |
| 1.2 | 3200-1812 | Drill Sleeve 4.5mm                 | 1    |
| 13  | 3200-1813 | Drill Sleeve 4.5mm                 | 1    |
| 14  | 3200-1814 | Hex Wrench                         | 1    |
| 15  | 3200-1815 | Plate & Screw Box                  | 1    |
| 16  | 3200-1816 | Target Mould 2 holes*10mm          | 1    |
| 17  | 3200-1817 | Target Mould 2 holes*12mm          | 1    |
| 18  | 3200-1818 | Target Mould 2 holes*18mm          | 1    |
| 19  | 3200-1819 | Target Mould 2 holes*20mm          | 1    |
| 20  | 3200-1820 | Target Mould 2 holes*25mm          | 1    |
| 21  | 3200-1821 | Target Mould 4 holes*16mm          | 1    |
| 22  | 3200-1822 | Target Mould 4 holes*20mm          | 2    |
| 23  | 3200-1823 | Target Mould 4 holes*12mm          | 1    |
| 24  | 3200-1824 | 1.5mm Locking Sleeve               | 1    |
| 25  | 3200-1825 | 1.5mm Locking Sleeve               | 1    |
| 26  | 3200-1826 | 2.5mm Drill Sleeve                 | 1    |
| 27  | 3200-1827 | 2.5mm Drill Sleeve                 | 1    |
| 28  | 3200-1828 | Sleeve                             | 1    |
| 29  | 3200-1829 | Sleeve                             | 1    |
| 30  | 3200-1830 | 10mm Target Mould Holder           | 1    |
| 31  | 3200-1831 | 10mm Target Mould Holder           | 1    |
| 32  | 3200-1832 | 1.5mm Wire Sleeve                  | 1    |
| 33  | 3200-1833 | 1.5mm Wire Sleeve                  | 1    |
| 34  | 3200-1834 | 3.2mm Drill Sleeve                 | 1    |
| 35  | 3200-1835 | 1.2mm Wire Sleeve                  | 1    |
| 36  | 3200-1836 | 1.2mm Wire Sleeve                  | 1    |
| 37  | 3200-1837 | 3.2mm Drill Sleeve                 | 1    |
| 38  | 3200-1838 | Cannulated Drill Bit 3.2*200       | 1    |
| 39  | 3200-1839 | Cannulated Drill Bit 3.2*200 Blcok | 1    |
| 40  | 3200-1840 | Cannulated Drill Bit 2.5*200       | 1    |
| 41  | 3200-1841 | Cannulated Drill Bit 2.5*200 Blcok | 1    |
| 42  | 3200-1842 | K-wire 1.2*150                     | 1    |
| 43  | 3200-1843 | K-wire 1.2*150                     | 1    |
| 44  | 3200-1844 | Aluminium Box                      | 1    |





2.7/3.5/5.0 Pediatric Osteotomy Locking  
Plate Instrument Set 3200-22





| NO. | REF         | Description                                | Qty. |
|-----|-------------|--------------------------------------------|------|
| 1   | 3200-2208   | T15 Star Screwdriver                       | 2    |
| 2   | 3200-2209   | T20 Star Screwdriver                       | 2    |
| 3   | 3200-2210   | T8 Star Screwdriver                        | 2    |
| 4   | 3200-2221   | 2.7mm Locking Sleeve                       | 3    |
| 5   | 3200-2205   | Tap HA2.7*96                               | 1    |
| 6   | 3200-2223   | 2.7mm Sleeve                               | 2    |
| 7   | 3200-2212   | AO Drill Bit 2.0*200 With Limitator        | 2    |
| 8   | 3200-2212-2 | φ2.0 Limitator                             | 2    |
| 9   | 3200-2213   | AO Drill Bit 2.5*200                       | 2    |
| 10  | 3200-2214   | K-wire Depth Gague 0-100mm                 | 1    |
| 11  | 3200-2220   | AO Drill Bit 2.9*250 With Limitator        | 2    |
| 12  | 3200-2220-2 | φ2.9 Limitator                             | 2    |
| 13  | 3200-2219   | AO Drill Bit 4.2*300 With Limitator        | 2    |
| 14  | 3200-2219-2 | φ4.2 Limitator                             | 2    |
| 15  | 3200-2215   | K-wire 2.0*230                             | 4    |
| 16  | 3200-2217   | Threaded K-wire 2.0*230                    | 4    |
| 17  | 3200-2216   | K-wire 2.5*250                             | 4    |
| 18  | 3200-2218   | Threaded K-wire 2.5*300                    | 4    |
| 19  | 3200-2232   | T-handle                                   | 1    |
| 20  | 3200-2229   | Triangle Location Plate 90°/50°/40°        | 1    |
| 21  | 3200-2230   | Triangle Location Plate 80°/70°/30°        | 1    |
| 22  | 3200-2231   | Triangle Location Plate 100°/60°/20°       | 1    |
| 23  | 3200-2225   | 2.7 /3.5Locating Device                    | 2    |
| 24  | 3200-2226   | 3.5/5.0 Locating Device                    | 2    |
| 25  | 3200-2207   | Tap HA3.5*140                              | 1    |
| 26  | 3200-2206   | Tap HA4.5*150                              | 1    |
| 27  | 3200-2211   | AO Drill Bit 3.2*150                       | 2    |
| 28  | 3200-2201   | 4.2mm Locking Sleeve                       | 3    |
| 29  | 3200-2203   | 3.0mm Locking Sleeve                       | 3    |
| 30  | 3200-2233   | Guider 2.0/2.7                             | 1    |
| 31  | 3200-2234   | Guider 2.5/3.5                             | 1    |
| 32  | 3200-2235   | Guider 3.2/4.5                             | 1    |
| 33  | 3200-2227   | 3.5/5.0 Angle Device 80-160°               | 1    |
| 34  | 3200-2228   | 2.7 Angle Device 80-160°                   | 1    |
| 35  | 3200-2238   | 2.7 Osteotome Locator                      | 1    |
| 36  | 3200-2239   | 3.5/5.0 Osteotome Locator                  | 1    |
| 37  | 4200-0214   | Self-centering Bone Holding Forcep (270mm) | 1    |
| 38  | 4200-0117   | Self-centering Bone Holding Forcep (190mm) | 1    |
| 39  | 4200-0210   | Obilique Reduction Forcecp (230mm)         | 1    |
| 40  | 4200-0119   | Obilique Reduction Forcecp (170mm)         | 1    |
| 41  | 3200-2236   | Straight Handle                            | 1    |
| 42  | 3200-2240   | 2.7 Target Device                          | 1    |
| 43  | 3200-2241   | 3.5 Target Device                          | 1    |
| 44  | 3200-2242   | 5.0 Target Device                          | 1    |
| 45  | 3200-2246   | Torque Wrench 4.0 N.m                      | 1    |
| 46  | 3200-2237   | Torque Wrench 1.5 N.m                      | 1    |
| 47  | 3200-2247   | Torque Wrench 0.8 N.m                      | 1    |
| 48  | 3200-2243   | 2.7 Depth Gague 0-50mm                     | 1    |
| 49  | 3200-2244   | 3.5 Depth Gague 10-60mm                    | 1    |
| 50  | 3200-2245   | 5.0 Depth Gague 10-100mm                   | 1    |
| 51  | 3200-2248   | Aluminium Box                              | 1    |

Mini Fragment Locking Plate  
Instrument Set 3200-10



| NO. | REF       | Description                    | Qty |
|-----|-----------|--------------------------------|-----|
| 1   | 3200-1001 | 1.5mm PlateTray                | 1   |
| 2   | 3200-1002 | 2.0mm Plate Tray               | 1   |
| 3   | 3200-1003 | 2.4mm Plate Tray               | 1   |
| 4   | 3200-1004 | 1.5mm Plate Holder             | 2   |
| 5   | 3200-1005 | 1.1mm AO Drill Bit             | 1   |
| 6   | 3200-1006 | T5 Stardriver                  | 1   |
|     | 3200-1042 | T5 Screw Holder                | 1   |
| 7   | 3200-1007 | 1.5mm AO Tap                   | 1   |
| 8   | 3200-1008 | 3.0mm Countersink              | 1   |
| 9   | 3200-1009 | 1.1/1.5 Guider                 | 1   |
| 10  | 3200-1010 | 1.1mm Locking Sleeve           | 1   |
| 11  | 3200-1011 | 1.1mm Locking Sleeve           | 1   |
| 12  | 3200-1012 | 2.0mm Plate Holder             | 2   |
| 13  | 3200-1013 | 1.5mm AO Drill Bit             | 1   |
| 14  | 3200-1014 | T6 Stardriver                  | 1   |
|     | 3200-1043 | T6 Screw Holder                | 1   |
| 15  | 3200-1015 | 2.0mm AO Tap                   | 1   |
| 16  | 3200-1016 | 4.0mm Countersink              | 1   |
| 17  | 3200-1017 | 1.5/2.0 Guider                 | 1   |
| 18  | 3200-1018 | 1.5mm Locking Sleeve           | 1   |
| 19  | 3200-1019 | 1.5mm Locking Sleeve           | 1   |
| 20  | 3200-1020 | 2.4mm Plate Holder             | 2   |
| 21  | 3200-1021 | 2.0mm AO Drill Bit             | 1   |
| 22  | 3200-1022 | T8 Stardriver                  | 1   |
|     | 3200-1043 | T8 Screw Holder                | 1   |
| 23  | 3200-1023 | 4.5mm Countersink              | 1   |
| 24  | 3200-1024 | 2.4mm AO Tap                   | 1   |
| 25  | 3200-1025 | 2.0/2.4 Guider                 | 1   |
| 26  | 3200-1026 | 2.0mm Locking Sleeve           | 1   |
| 27  | 3200-1027 | 2.0mm Locking Sleeve           | 1   |
| 28  | 3200-1028 | Bending Forcecp                | 1   |
| 29  | 3200-1029 | Plate Cutter                   | 1   |
| 30  | 3200-1030 | 0.4 N.m Torque Wrench          | 1   |
| 31  | 3200-1031 | 0.8 N.m Torque Wrench          | 1   |
| 32  | 3200-1032 | Self-locking Reduction Forcecp | 1   |
| 33  | 3200-1033 | Screwholding Forcecp           | 1   |
| 34  | 3200-1034 | AO Quick Coupling Handle       | 1   |
| 35  | 3200-1035 | Tip Reduction Forcecp          | 1   |
| 36  | 3200-1036 | Depth Gague 0-30mm             | 1   |
| 37  | 3200-1037 | Depth Gague 0-40mm             | 1   |
| 38  | 3200-1038 | Periosteal Elevator Flat       | 1   |
| 39  | 3200-1039 | Periosteal Elevator Round      | 1   |
| 40  | 3200-1040 | Aluminium Box                  | 1   |




2.7mm Mini Fragment Locking Plate  
Instrument Set 3200-11



| NO. | REF       | Description          | Qty |
|-----|-----------|----------------------|-----|
| 1   | 3200-1101 | 2.7mm PlateTray      | 1   |
| 2   | 3200-1102 | 2.0/2.7 Guider       | 1   |
| 3   | 3200-1103 | 4.5mm Countersink    | 1   |
| 4   | 3200-1104 | 2.7mm AO Tap         | 1   |
| 5   | 3200-1105 | T8 Stardriver        | 1   |
|     | 3200-1111 | T8 Screw Holder      | 1   |
| 6   | 3200-1106 | 2.0mm AO Drill Bit   | 1   |
| 7   | 3200-1107 | 2.7mm Plate Holder   | 2   |
| 8   | 3200-1108 | 2.7mm Locking Sleeve | 1   |
| 9   | 3200-1109 | 2.7mm Locking Sleeve | 1   |






Small Fragment Locking Plate  
Instrument Set (AO) 3200-13



| NO.  | REF         | Description                                | Qty. |
|------|-------------|--------------------------------------------|------|
| 1    | 3200-1301   | Hex Key SW2.5                              | 1    |
| 2    | 3200-1302   | Locking Sleeve 3.0mm for Thin Plates       | 1    |
| 3    | 3200-1302   | Locking Sleeve 3.0mm for Thin Plates       | 1    |
| 4    | 3200-1304   | Wire Sleeve 1.2mm                          | 1    |
| 5    | 3200-1305   | Locking Sleeve 3.0mm for Thick Plates      | 1    |
| 6    | 3200-1305   | Locking Sleeve 3.0mm for Thick Plates      | 1    |
| 7    | 3200-1305   | Locking Sleeve 3.0mm for Thick Plates      | 1    |
| 8    | 3200-1304   | Wire Sleeve 1.2mm                          | 1    |
| 9    | 3200-1309   | Threaded Guider Wire $\phi$ 1.2*150mm      | 1    |
| 10   | 3200-1310   | Threaded Guider Wire $\phi$ 1.2*150mm      | 1    |
| 11   | 3200-1311   | K Wire $\phi$ 1.2*150mm                    | 1    |
| 12   | 3200-1311   | K Wire $\phi$ 1.2*150mm                    | 1    |
| 13   | 3200-1313   | AO Drill Bit $\phi$ 2.5*200                | 2    |
| 14   | 3200-1346   | AO Drill Bit $\phi$ 1.9*200                | 2    |
| 15   | 3200-1315   | AO Drill Bit $\phi$ 2.9*250 with Block     | 1    |
| 15.1 | 3200-1315-1 | Block                                      | /    |
| 16   | 3200-1316   | AO Drill Bit $\phi$ 2.9*250                | 1    |
| 17   | 3200-1317   | Countersink $\phi$ 6.0                     | 1    |
| 18   | 3200-1318   | Hollow $\phi$ 6.0                          | 1    |
| 19   | 3200-1319   | Retractor 8.5mm (0-60mm)                   | 1    |
| 20   | 3200-1320   | Retractor 15.5mm                           | 1    |
| 21   | 3200-1321   | Depth Gauge                                | 1    |
| 22   | 3200-1322   | T-handle Quick Coupling AO                 | 1    |
| 23   | 3200-1323   | Periosteal Dissector 12mm                  | 1    |
| 24   | 3200-1324   | Periosteal Dissector 8mm                   | 1    |
| 25   | 3200-1325   | Drill & Tap Sleeve combined 2.5/3.5        | 1    |
| 26   | 3200-1326   | Guider 2.0/4.0                             | 1    |
| 27   | 3200-1327   | Transition Interface                       | 1    |
| 28   | 3200-1328   | Screwholding Forcep                        | 1    |
| 29   | 3200-1329   | Extraction Screw                           | 1    |
| 30   | 3200-1330   | Sleeve Key $\phi$ 4.0                      | 1    |
| 31   | 3200-1331   | Star Screwdriver T15                       | 2    |
| 32   | 3200-1345   | Star Screw-Holder (Open Type)              | 1    |
| 33   | 3200-1333   | Tap for Cancellous Screw 4.0mm             | 1    |
| 34   | 3200-1334   | Tap for Locking Screw HC3.5*125mm          | 1    |
|      | 3200-1335   | Tap for Cortical Screw HA3.5*140mm         | 1    |
| 35   | 3200-1336   | Staight Handle Qucik Coupling AO           | 1    |
| 36   | 3200-1337   | Sharp Reduction Forcep (190mm)             | 1    |
| 37   | 3200-1338   | Torque Handle 1.5N.m                       | 1    |
| 38   | 3200-1339   | Self-centering Bone Holding Forcep (190mm) | 2    |
| 39   | 3200-1340   | Obilique Reduction Forcep (170mm)          | 1    |
| 40   | 3200-1341   | Bending Iron                               | 1    |
| 40.1 | 3200-1341-1 | Bending Iron L                             | /    |
| 40.2 | 3200-1341-2 | Bending Iron R                             | /    |
| 41   | 3200-1342-N | Aluminium Box                              | 1    |


Large Fragment Locking Plate  
Instrument Set(AO)3200-14



| NO.  | REF         | Description                                                         | Qty. |
|------|-------------|---------------------------------------------------------------------|------|
| 1    | 3200-1401   | Drill & Tap Sleeve combined 3.2/4.5mm                               | 1    |
| 2    | 3200-1402   | Block Hex Key Small SW2.5                                           | 1    |
|      | 3200-1403   | Block Hex Key Large SW3.0                                           | 1    |
| 3    | 3200-1404   | Sleeve Key $\phi$ 4.0                                               | 1    |
| 4    | 3200-1405   | Extraction Screw                                                    | 1    |
| 5    | 3200-1406   | Star Screwdriver T20 (Open Type)                                    | 1    |
| 6    | 3200-1407   | Star Screwdriver T20                                                | 2    |
| 7    | 3200-1408   | AO Drill Bit $\phi$ 3.2                                             | 1    |
| 8    | 3200-1409   | AO Drill Bit $\phi$ 3.2                                             | 1    |
| 9    | 3200-1410   | Tap for Cortical Screw 4.5mm                                        | 1    |
| 10   | 3200-1411   | Tap for Locking Screw 5.0mm                                         | 1    |
| 11   | 3200-1412   | AO Drill Bit $\phi$ 4.2/AO Drill Bit $\phi$ 4.2 with Limiator Block | 1+1  |
| 11.1 | 3200-1412-1 | AO Drill Bit $\phi$ 4.2*290                                         | /    |
| 11.2 | 3200-1412-2 | $\phi$ 4.2 with Limiator Block                                      | /    |
| 12   | 3200-1413   | 5.0 Cannulated Drill Bit with Block                                 | 1    |
| 12.1 | 3200-1413-1 | $\phi$ 5.0 *280 Cannulated Drill                                    | /    |
| 12.2 | 3200-1413-2 | $\phi$ 5.0 with Limiator Block                                      | /    |
| 13   | 3200-1414   | Threaded K Wire $\phi$ 2.5*300mm                                    | 2    |
| 14   | 3200-1415   | K Wire/Threaded Guider Wire $\phi$ 2.0*250mm                        | 2+2  |
| 14.1 | 3200-1415-1 | K Wire $\phi$ 2.0*250mm                                             | /    |
| 14.2 | 3200-1415-2 | Threaded Guider Wire $\phi$ 2.0*250mm                               | /    |
| 15   | 3200-1416   | Locking Sleeve 5.0mm                                                | 1    |
| 16   | 3200-1417   | Wire Sleeve 2.5*5.0mm                                               | 1    |
| 17   | 3200-1418   | Countersink $\phi$ 8.0                                              | 1    |
| 18   | 3200-1419   | Hollow $\phi$ 8.0                                                   | 1    |
| 19   | 3200-1420   | Sharp Reduction Forcep (200mm)                                      | 1    |
| 20   | 3200-1421   | 6.5 Cannulated Locking Screwdriver T27                              | 1    |
| 21   | 3200-1422   | Periosteal Dissector 8mm                                            | 1    |
| 22   | 3200-1423   | Periosteal Dissector 15mm                                           | 1    |
| 23   | 3200-1424   | Screwholding Forcep                                                 | 1    |
| 24   | 3200-1425   | T-handle Quick Coupling AO                                          | 1    |
| 25   | 3200-1426   | Transition Interface                                                | 1    |
| 26   | 3200-1427   | Locking Sleeve 4.2mm for Thin Plates                                | 2    |
| 27   | 3200-1428   | Locking Sleeve 4.2mm for Thick Plates                               | 1    |
| 28   | 3200-1428   | Locking Sleeve 4.2mm for Thick Plates                               | 1    |
| 29   | 3200-1428   | Locking Sleeve 4.2mm for Thick Plates                               | 1    |
| 30   | 3200-1431   | Wire Sleeve 2.0*4.2mm                                               | 1    |
| 31   | 3200-1431   | Wire Sleeve 2.0*4.2mm                                               | 1    |
| 32   | 3200-1433   | Depth Gauge 40-90mm                                                 | 1    |
| 33   | 3200-1434   | Bending Iron                                                        | 1    |
| 33.1 | 3200-1434-1 | Bending Iron L                                                      | /    |
| 33.2 | 3200-1434-2 | Bending Iron R                                                      | /    |
| 34   | 3200-1435   | Obilique Reduction Forcep (230mm)                                   | 1    |
| 35   | 3200-1436   | Torque Handle 4.0N.m                                                | 1    |
| 36   | 3200-1437   | Self-centering Bone Hold Forcep(270mm)                              | 2    |
| 37   | 3200-1438   | Dissector                                                           | 2    |
| 37.1 | 3200-1438-1 | Dissector 16mm                                                      | /    |
| 37.2 | 3200-1438-2 | Dissector 44mm                                                      | /    |
| 38   | 3200-1439   | Staight Handle Qucik Coupling                                       | 1    |
| 39   | 3200-1440-N | Aluminium Box                                                       | 1    |




2.4/2.7 Multi-axial Locking Plates  
Instrument Set 3200-16



| NO. | REF         | Description                               | Qty |
|-----|-------------|-------------------------------------------|-----|
| 1   | 3200-1601   | Alumium Box & Screw Rack                  | 1   |
| 2   | 3200-1602   | Multi-axial Cone Guider Ø2.0              | 1   |
| 3   | 3200-1603   | Multi-axial Cone Guider Ø1.8              | 1   |
| 4   | 3200-1604   | Depth Gauge (0-40mm)                      | 1   |
| 5   | 3200-1605   | Torx Head Screwdriver T8                  | 1   |
| 6   | 3200-1606   | Multi-axial Straight Guider Ø2.0          | 1   |
| 7   | 3200-1607   | Multi-axial Straight Guider Ø1.8          | 1   |
| 8   | 3200-1608   | Screwdriver Handle                        | 1   |
| 9   | 3200-1609   | Drill Sleeve Ø1.8                         | 1   |
| 10  | 3200-1610   | Drill Sleeve Ø1.8                         | 1   |
| 11  | 3200-1611   | Drill Sleeve Ø2.0                         | 1   |
| 12  | 3200-1612   | Drill Sleeve Ø2.0                         | 1   |
| 13  | 3200-1613   | Plate Bender 2.4/2.7mm                    | 1   |
| 14  | 3200-1614   | Plate Bender 2.4/2.7mm                    | 1   |
| 15  | 3200-1615   | Drill Guider Ø2.0/2.7                     | 1   |
| 16  | 3200-1616   | Drill Guider Ø1.8/2.4                     | 1   |
| 17  | 3200-1617   | AO Drill Bit Ø1.8mm                       | 1   |
| 18  | 3200-1618   | AO Drill Bit Ø1.8mm                       | 1   |
| 19  | 3200-1619   | AO Drill Bit Ø2.0mm                       | 1   |
| 20  | 3200-1620   | AO Drill Bit Ø2.0mm                       | 1   |
| 21  | 3200-1621   | Drill Guider Ø2.0/2.7                     | 1   |
| 22  | 3200-1622   | Drill Guider Ø1.8/2.4                     | 1   |
| 23  | 3200-1623   | Torque Wrench 0.8N.m.                     | 1   |
| 24  | 3200-1624   | Screw Holding Sleeve                      | 1   |
| 25  | 3200-1625   | Countersink Drill Ø4.5                    | 1   |
| 26  | 3200-1626   | Countersink Drill Ø4.5                    | 1   |
| 27  | 3200-1627   | Torx Head Screwdriver T8 (Quick-coupling) | 1   |
| 28  | 3200-1628   | Quick Coupling Tap HC2.4                  | 1   |
| 29  | 3200-1629   | Quick Coupling Tap HA2.4                  | 1   |
| 30  | 3200-1630   | Quick Coupling Tap HA2.7                  | 1   |
| 31  | 3200-1631   | Quick Coupling Tap HC2.7                  | 1   |
| 32  | 3200-1632   | Guider Pin 1.2*150                        | 1   |
| 33  | 3200-1633   | Guider Pin 1.2*150                        | 1   |
| 34  | 3200-1634-N | Aluminium Box                             | 1   |




Osteotomy Locking Plate Instrument  
Set 3200-08



| NO. | REF         | Description                                     | Qty |
|-----|-------------|-------------------------------------------------|-----|
| 1   | 3200-0812   | Depth Gauge 0-120mm                             | 1   |
| 2   | 3200-0809   | Drill Bit Kirschner Wire Ø2.5*300               | 2   |
| 3   | 3200-0810   | Pointed Kirschner Wire Ø2.5*280                 | 2   |
| 4   | 3200-0811   | Pointed Kirschner Wire Ø2.0*280                 | 2   |
| 5   | 3200-0805   | Osteotome 25mm                                  | 1   |
| 6   | 3200-0806   | Osteotome 20mm                                  | 1   |
| 7   | 3200-0807   | Osteotome 15mm                                  | 1   |
| 8   | 3200-0808   | Osteotome 10mm                                  | 1   |
| 9   | 3200-0801   | Qucik Coupling Drill Bit Ø4.3*220               | 2   |
| 9.1 | 3200-0801-1 | Qucik Coupling Drill Bit Ø4.3*220 With Measures | /   |
| 9.2 | 3200-0801-2 | Ø4.3 Limitator                                  | /   |
| 10  | 3200-0802   | Qucik Coupling Tap HC 5.0*220                   | 1   |
| 11  | 3200-0803   | Qucik Coupling Drill Bit Ø3.2*195               | 1   |
| 12  | 3200-0803   | Qucik Coupling Drill Bit Ø3.2*195               | 1   |
| 13  | 3200-0815   | Torque Handle 4.0N.m                            | 1   |
| 14  | 3200-0818   | Retractor Forcep                                | 1   |
| 15  | 3200-0813   | Angle Measurer                                  | 1   |
| 16  | 3200-0816   | Spring Drill Ø3.2/4.3                           | 1   |
| 17  | 3200-0817   | Straight Quick Coupling Handle                  | 1   |
| 18  | 3200-0819   | Quick coupling Hex Screwdriver T4.3*100         | 1   |
| 19  | 3200-0820   | Tap Guide Drill Sleeve Ø5.0*100                 | 1   |
| 20  | 3200-0823   | Drill-oriented Drill Sleeve Ø4.3*100            | 1   |
| 21  | 3200-0824   | Drill-oriented Drill Sleeve Ø4.3*100            | 1   |
| 22  | 3200-0827   | Measurer                                        | 1   |
| 23  | 3200-0821   | Guide Pin Duide Ø2.0                            | 1   |
| 24  | 3200-0821   | Guide Pin Duide Ø2.0                            | 1   |
| 25  | 3200-0814   | Section Height Monitor                          | 1   |
| 26  | 3200-0828   | Section Height Monitor Screwdriver T20/130      | 1   |
| 27  | 3200-0825-N | Aluminium Box                                   | 1   |





Ankle Locking Plate Instrument Set  
3200-27





| NO. | REF         | Description                                         | Qty |
|-----|-------------|-----------------------------------------------------|-----|
| 1   | 3200-2712   | AO Drill Bit $\Phi$ 2.9x150                         | 2   |
| 2   | 3200-2713   | AO Drill Bit $\Phi$ 2.5x150                         | 2   |
| 3   | 3200-2714   | AO Drill Bit $\Phi$ 2.0x150                         | 2   |
| 4   | 3200-2715   | AO Drill Bit $\Phi$ 1.9x150                         | 2   |
| 5   | 3200-2702   | Threaded Wire $\Phi$ 2.5*150                        | 5   |
| 6   | 3200-2701   | Threaded Wire $\Phi$ 1.5*150                        | 5   |
| 7   | 3200-2737   | Bone Hook Double Head 180x9x22                      | 2   |
| 8   | 3200-2721   | Locking Sleeve $\Phi$ 3.0*50                        | 1   |
| 9   | 3200-2721   | Locking Sleeve $\Phi$ 3.0*50                        | 1   |
| 10  | 3200-2716   | AO T-handle                                         | 1   |
| 11  | 3200-2720   | Locking Sleeve $\Phi$ 2.0*50                        | 1   |
| 12  | 3200-2720   | Locking Sleeve $\Phi$ 2.0*50                        | 1   |
| 13  | 3200-2725   | Neutral and Load Drill Guide $\Phi$ 2.5/ $\Phi$ 3.5 | 1   |
| 14  | 3200-2726   | Neutral and Load Drill Guide $\Phi$ 2.0/ $\Phi$ 2.7 | 1   |
| 15  | 3200-2727   | Screw Holding Forcep                                | 1   |
| 16  | 3200-2719   | Wire Sleeve $\Phi$ 1.5x70                           | 1   |
| 17  | 3200-2719   | Wire Sleeve $\Phi$ 1.5x70                           | 1   |
| 18  | 3200-2710   | Depth Gauge 0-30mm                                  | 1   |
| 19  | 3200-2711   | Depth Gauge 0-60mm                                  | 1   |
| 20  | 3200-2736   | Periosteal Elevator 8mm                             | 1   |
| 21  | 3200-2739   | Bone Elevator 160x5mm                               | 1   |
| 22  | 3200-2740   | Bone Elevator 160x8mm                               | 1   |
| 23  | 3200-2741   | Bone Elevator 160x13mm                              | 1   |
| 24  | 3200-2738   | Skin Retractor 2 Teeth 170x5                        | 2   |
| 25  | 3200-2734   | Dissector Double Head 3mm                           | 1   |
| 26  | 3200-2735   | Dissector Double Head 5mm                           | 1   |
| 27  | 3200-2722   | HC3.5 Multi-angle Guider 30mm                       | 1   |
| 28  | 3200-2704   | Torque Wrench 0.8NM                                 | 1   |
| 29  | 3200-2718   | Small Straight Handle AO                            | 1   |
| 30  | 3200-2703   | Torque Wrench 1.5NM                                 | 1   |
| 31  | 3200-2723   | Torx Screwdriver T8*106                             | 2   |
| 32  | 3200-2724   | Torx Screwdriver T15*110mm                          | 2   |
| 33  | 3200-2742   | HC3.5 Holding Device 110mm                          | 2   |
| 34  | 3200-2706   | AO Tap 2.7mm Cortical                               | 1   |
| 35  | 3200-2707   | AO Tap 2.7mm Locking                                | 1   |
| 36  | 3200-2708   | AO Tap 3.5mm Cortical                               | 1   |
| 37  | 3200-2705   | AO Tap 3.5mm Locking                                | 1   |
| 38  | 3200-2709   | AO Tap 4.0mm Cancellous                             | 1   |
| 39  | 3200-2717   | Straight Handle AO                                  | 1   |
| 40  | 3200-2728   | Sharp Reduction Forcep                              | 1   |
| 41  | 3200-2731   | Plate Cutter                                        | 1   |
| 42  | 3200-2730   | Plate Bender                                        | 1   |
| 43  | 3200-2733   | K-wire Distractor 200x $\Phi$ 1.6x $\Phi$ 2.6       | 1   |
|     | 3200-2732   | Knee Distractor 190x9.5x25                          | 1   |
| 44  | 3200-2730   | Plate Bender                                        | 1   |
| 45  | 3200-2743   | Screw Rack                                          | 1   |
| 46  | 3200-2744-N | Aluminium Box                                       | 1   |

Rib Reconstruction Locking Plate  
3200-09



| NO. | REF         | Description                | Qty |
|-----|-------------|----------------------------|-----|
| 1   | 3200-0901   | Locking Sleeve             | 1   |
| 2   | 3200-0902   | Locking Sleeve             | 1   |
| 3   | 3200-0903   | Screwdriver                | 1   |
| 4   | 3200-0904   | Screwdriver                | 1   |
| 5   | 3200-0905   | Limited Drill Bit 2.2*8mm  | 1   |
| 6   | 3200-0906   | Limited Drill Bit 2.2*10mm | 1   |
| 7   | 3200-0907   | Limited Drill Bit 2.2*12mm | 1   |
| 8   | 3200-0908   | Limited Drill Bit 2.2*14mm | 1   |
| 9   | 3200-0909   | Limited Drill Bit 2.2*16mm | 1   |
| 10  | 3200-0910   | Limited Drill Bit 2.2*18mm | 1   |
| 11  | 3200-0911   | Limited Drill Bit 2.2*20mm | 1   |
| 12  | 3200-0912   | Depth Gauge                | 1   |
| 13  | 3200-0913   | Straight Handle            | 1   |
| 14  | 3200-0913   | Straight Handle            | 1   |
| 15  | 3200-0915   | Plate Holding Forcep       | 1   |
| 16  | 3200-0915   | Plate Holding Forcep       | 1   |
| 17  | 3200-0917   | Reduction Forcep           | 1   |
| 18  | 3200-0917   | Reduction Forcep           | 1   |
| 19  | 3200-0919   | Plate Bender               | 1   |
| 20  | 3200-0920   | Plate Bender L             | 1   |
| 21  | 3200-0921   | Plate Bender R             | 1   |
| 22  | 3200-0922   | Drill Guider               | 1   |
| 23  | 3200-0923-N | Aluminium Box              | 1   |


Sternum Locking Plate Instrument Set  
3200-17



| NO. | REF       | Description                              | Qty |
|-----|-----------|------------------------------------------|-----|
| 1   | 3200-1701 | Screw Candy                              | 1   |
| 2   | 3200-1702 | Temporary holder HC3.0                   | 1   |
| 3   | 3200-1703 | Guide drill 2.3                          | 2   |
| 4   | 3200-1704 | Temporary holder HC3.0                   | 2   |
| 5   | 3200-1705 | Quadrangle screwdriver SW1.7             | 1   |
| 6   | 3200-1706 | 2.3 Limited Drill Bit 6mm                | 1   |
| 7   | 3200-1707 | 2.3 Limited Drill Bit 8mm                | 1   |
| 8   | 3200-1708 | 2.3 Limited Drill Bit 10mm               | 1   |
| 9   | 3200-1709 | 2.3 Limited Drill Bit 12mm               | 1   |
| 10  | 3200-1710 | 2.3 Limited Drill Bit 14mm               | 1   |
| 11  | 3200-1711 | 2.3 Limited Drill Bit 16mm               | 1   |
| 12  | 3200-1712 | Plate bender L=200mm                     | 2   |
| 13  | 3200-1713 | Plate shears L=190mm                     | 2   |
| 14  | 3200-1714 | Plate Holding Compression Spreader HC3.0 | 1   |
| 15  | 3200-1715 | Percutaneous Guiding Drill 3.0           | 1   |
| 16  | 3200-1716 | Depth Gauge 0-30mm                       | 1   |
| 17  | 3200-1717 | Straight Handle                          | 1   |
| 18  | 3200-1718 | Aluminium Box                            | 1   |




### FNS Instrument Set 3300-04



| NO. | REF       | Description                                              | Qty |
|-----|-----------|----------------------------------------------------------|-----|
| 1   | 3300-0401 | Torque Limiting Attachment, 4 Nm, for AO Reaming Coupler | 1   |
| 2   | 3300-0402 | Large Cannulated Handle, Quick Coupling, 12 mm Hexagonal | 1   |
| 3   | 3300-0403 | Transition Interface                                     | 1   |
| 4   | 3300-0404 | Protection Sleeve                                        | 1   |
| 5   | 3300-0405 | T25 StarDrive Screwdriver Shaft                          | 1   |
| 6   | 3300-0405 | T25 StarDrive Screwdriver Shaft                          | 1   |
| 7   | 3300-0407 | T20 StarDrive Screwdriver Shaft                          | 1   |
| 8   | 3300-0408 | Multifunction Rod, for Insertion Instruments             | 1   |
| 9   | 3300-0409 | Depth Gauge up to 100 mm                                 | 1   |
| 10  | 3300-0410 | Direct Measuring Device, for 3.2 mm Guide Wires          | 1   |
| 11  | 3300-0411 | Insertion Handle                                         | 1   |
| 12  | 3300-0412 | 130 Degree Angled Guide, for 3.2 mm Guide Wires          | 1   |
| 13  | 3300-0413 | Insert for Insertion Handle                              | 1   |
| 14  | 3300-0414 | Cylinder, for Insertion Instruments                      | 1   |
| 15  | 3300-0415 | Correction Guide, for 3.2 mm Guide Wires                 | 1   |
| 16  | 3300-0416 | Complete Opening Drill Bit/Reamer                        | 1   |
| 17  | 3300-0417 | 3.2 mm Guide Wire, 400 mm                                | 4   |
| 18  | 3300-0418 | 4.3 mm Drill Bi with stopper, 413 mm Length              | 1   |
| 19  | 3300-0419 | Aluminium Box                                            | 1   |

### Subtalar Screw Instrument Set 3200-40




| NO. | REF       | Description                  | Qty |
|-----|-----------|------------------------------|-----|
| 1   | 3200-4017 | Peeler                       | 1   |
| 2   | 3200-4015 | Cannulated Screwdriver SW4.0 | 1   |
| 3   | 3200-4016 | Cannulated Screwdriver SW4.5 | 1   |
| 4   | 3200-4010 | Trail $\Phi$ 5.5 $\times$ 21 | 1   |
| 5   | 3200-4011 | Trail $\Phi$ 6.5 $\times$ 21 | 1   |
| 6   | 3200-4012 | Trail $\Phi$ 7.3 $\times$ 21 | 1   |
| 7   | 3200-4013 | Trail $\Phi$ 8.5 $\times$ 21 | 1   |
| 8   | 3200-4014 | Trail $\Phi$ 9.5 $\times$ 21 | 1   |
| 9   | 3200-4006 | Trail $\Phi$ 8.5 $\times$ 15 | 1   |
| 10  | 3200-4007 | Trail $\Phi$ 10 $\times$ 15  | 1   |
| 11  | 3200-4008 | Trail $\Phi$ 11 $\times$ 15  | 1   |
| 12  | 3200-4009 | Trail $\Phi$ 12 $\times$ 15  | 1   |
| 13  | 3200-4005 | Bone Hook                    | 1   |
| 14  | 3200-4004 | Screw Holding Forcep         | 1   |
| 15  | 3200-4003 | Bone Cutter                  | 1   |
| 16  | 3200-4001 | Rechat Handle                | 1   |
| 17  | 3200-4002 | K-wire 2 $\times$ 250mm      | 2   |
| 18  | 3200-4019 | Alumium Box                  | 1   |

### Small Fragment Locking Screw Rack 1000-0107




Size:40\*14\*11cm  
 3.5mm Locking Screw Length:12-95mm  
 3.5mm Cortex Screw Length:12-60mm  
 4.0mm Cancellous Screw Length:12-50mm

### Large Fragment Locking Screw Rack 1000-0108



Size:35\*13\*12cm  
 5.0mm Locking Screw Length:22-95mm  
 4.5mm Cortex Screw Length:22-60mm  
 6.5mm Cannulated Locking Screw Length:50-110mm

### Radius Plate Box 1000-0128



For 3-7 holes VA Radius Locking Plate and Locking Screw



# CONTÁCTANOS



+56 9 9851 8862



SURGINEST\_GLOBAL



SURGINEST GLOBAL



**SURGINEST  
GLOBAL**  
INSTRUMENTAL Y EQUIPOS MÉDICOS