

SPINE

Instrumental de Acero Quirúrgico

CATÁLOGO



**SURGINEST
GLOBAL**

INSTRUMENTAL Y EQUIPOS MÉDICOS

PROFESSIONAL OPTICAL ENDOSCOPE

SPINE CATALOG



- | TRANSFORAMINAL ENDOSCOPE
- | INTERLAMINAR ENDOSCOPE
- | UBE ENDOSCOPE
- | ENDOSCOPE INSTRUMENT



**SURGINEST
GLOBAL**
INSTRUMENTAL Y EQUIPOS MÉDICOS

Table of Content:

Transforaminal & Interlaminar Endoscopic Surgical System:

3.3/5.3/5.9/6.3/6.9/7.3/8.4/10mm Transforaminal & Interlaminar Endoscope

Transforaminal & Interlaminar Endoscopic Instrument:

Rongeur Forceps

Nucleus Pulposus Forceps

Punch Forceps

Grasping Forceps

Flexible Forceps

Guide Bar Dilator

Cannula(Endoscopic Working Cannula)

Cannula(Endoscopic Working Cannula With Sealing Cap)

Metal Puncture Needle

Guide Needle

Trephine(Endoscope Use)

Safe Trephine/Safe Trephine Cannula

Visual Trephine

Safe Bone Drill

Visual Reamer

Nerve Hook/Nerve Dissector/Osteotome/Nerve probe/Operating knife

Bone Hammer

Large Channel Interlaminar Endoscope

Large Channel Rongeur Forceps

Large Channel Nucleus Pulposus Forceps

Delta Cannula

Delta Trephine

Delta Guide Needle/Dilator/Operating Knife

UBE Endoscope And Instruments

Disposable RF & Plasma Surgical Electrode

Endoscope And Instruments Box

Transforaminal / Interlaminar endoscope



Model number	Outerside diameter	Working length	Angle	Field angle	Working channel	Notes
CY3309000	3.3mm	90mm	0°	60°	2.1mm	HD autoclave
CY5313030	5.3mm	130mm	30°	80°	2.7mm	HD autoclave
CY5318130	5.3mm	181mm	30°	80°	2.7mm	HD autoclave
CY5913630	5.9mm	136mm	30°	80°	3.3mm	HD autoclave
CY5918130	5.9mm	181mm	30°	80°	3.3mm	HD autoclave
CY6313015	6.3mm	130mm	15°	80°	3.7mm	HD autoclave
★ CY6313030	6.3mm	130mm	30°	80°	3.7mm	HD autoclave
CY6318100	6.3mm	181mm	0°	80°	3.7mm	HD autoclave
CY6318115	6.3mm	181mm	15°	80°	3.7mm	HD autoclave
★ CY6318130	6.3mm	181mm	30°	80°	3.7mm	HD autoclave
CY6319330	6.3mm	193mm	30°	80°	3.7mm	HD autoclave
CY6913015	6.9mm	130mm	15°	80°	4.3mm	HD autoclave
★ CY6913030	6.9mm	130mm	30°	80°	4.3mm	HD autoclave
CY6918100	6.9mm	181mm	0°	80°	4.3mm	HD autoclave
CY6918115	6.9mm	181mm	15°	80°	4.3mm	HD autoclave
★ CY6918130	6.9mm	181mm	30°	80°	4.3mm	HD autoclave
CY6919330	6.9mm	193mm	30°	80°	4.3mm	HD autoclave
CY7313030	7.3mm	130mm	30°	80°	4.7mm	HD autoclave
CY7318130	7.3mm	181mm	30°	80°	4.7mm	HD autoclave
★ CY8412015	8.4mm	120mm	15°	80°	5.7mm	HD autoclave
CY1013630	10mm	136mm	30°	80°	7.1mm	HD autoclave
CY1013615	10mm	136mm	15°	80°	7.1mm	HD autoclave
★ CY1012515	10mm	125mm	15°	80°	7.3mm	HD autoclave




Rongeur forceps

Model number	Outerside diameter	Working length	Angle
CR2233045	2.2mm	330mm	45°
CR2833045	2.8mm	330mm	45°
CR3533045	3.5mm	330mm	45°
CR3833045	3.8mm	330mm	45°
CR4033045	4.0mm	330mm	45°
CR2233090	2.2mm	330mm	90°
CR2833090	2.8mm	330mm	90°
CR3533090	3.5mm	330mm	90°
CR3833090	3.8mm	330mm	90°
CR4033090	4.0mm	330mm	90°
CR3533045B(Head 1.7mm)	3.5mm	330mm	45°
CR3526045B(Head 1.7mm)	3.5mm	260mm	45°
CR3533090B(Head 1.7mm)	3.5mm	330mm	90°
CR3526090B(Head 1.7mm)	3.5mm	260mm	90°
CR2226045	2.2mm	260mm	45°
CR2826045	2.8mm	260mm	45°
CR3526045	3.5mm	260mm	45°
CR3826045	3.8mm	260mm	45°
CR4026045	4.0mm	260mm	45°
CR2226090	2.2mm	260mm	90°
CR2826090	2.8mm	260mm	90°
CR3526090	3.5mm	260mm	90°
CR3826090	3.8mm	260mm	90°
CR4026090	4.0mm	260mm	90°



Transforaminal / Interlaminar endoscopic surgical system

Nucleus pulposus forceps/Punch forceps/Grasping forceps/Flexible forceps



Product name	Model number	Diameter	Working length	Angle	Functional characteristics
Nucleus pulposus forceps	R253300BS	2.5mm	330mm	0°	Overload protection
Nucleus pulposus forceps	R303300BS	3.0mm	330mm	0°	Overload protection
Nucleus pulposus forceps	R353300BS	3.5mm	330mm	0°	Overload protection
Nucleus pulposus forceps	R403300BS	4.0mm	330mm	0°	Overload protection
Nucleus pulposus forceps	R253345BS	2.5mm	330mm	45°	Overload protection
Nucleus pulposus forceps	R303345BS	3.0mm	330mm	45°	Overload protection
Nucleus pulposus forceps	R253300CS	2.5mm	330mm	0°	Overload protection
Nucleus pulposus forceps	R303300CS	3.0mm	330mm	0°	Overload protection
Nucleus pulposus forceps	R353300CS	3.5mm	330mm	0°	Overload protection
Nucleus pulposus forceps	R253315CS	2.5mm	330mm	15°	Overload protection
Punch forceps	R253300BC	2.5mm	330mm	0°	Overload protection
Punch forceps	R303300BC	3.0mm	330mm	0°	Overload protection
Punch forceps	R403300BC	4.0mm	330mm	0°	Overload protection
Punch forceps	R253315BC	2.5mm	330mm	15°	Overload protection
Punch forceps	R303315BC	3.0mm	330mm	15°	Overload protection
Punch forceps	R253345BC	2.5mm	330mm	45°	Overload protection
Punch forceps	R353345BC	3.5mm	330mm	45°	Overload protection
Grasping forceps, pawl	R253300FB	2.5mm	330mm	0°	Overload protection
Grasping forceps, pawl	R303300FB	3.0mm	330mm	0°	Overload protection
Grasping forceps, pawl	R353300FB	3.5mm	330mm	0°	Overload protection
Grasping forceps, pawl	R303300FT	3.0mm	330mm	0°	Overload protection
Suture cutter forceps	R303300HS	3.0mm	330mm	0°	Overload protection
Flexible forceps	F253300ST	2.5mm	330mm	0°	Overload protection
Flexible forceps	F303300ST	3.0mm	330mm	0°	Overload protection
Flexible forceps	F253345BS	2.5mm	330mm	45°	Overload protection
Flexible forceps	F303345BS	3.0mm	330mm	45°	Overload protection
Flexible forceps	F253300BC	2.5mm	330mm	0°	Overload protection
Flexible forceps	F303300BC	3.0mm	330mm	0°	Overload protection
Flexible forceps	F253300CS	2.5mm	330mm	0°	Overload protection
Flexible forceps	F303300CS	3.0mm	330mm	0°	Overload protection



With overload and pressure relief device
protection - improve the service life of the instrument



BC
bite cut 0°
Flexible forceps



BS
bowl shape 45°
Flexible forceps



CS
cup shape 0°
Flexible forceps



ST
saw-tooth 0°
Flexible forceps



Product name	Model number	Diameter	Working length	Angle	Functional characteristics
Nucleus pulposus forceps	R252600BS	2.5mm	260mm	0°	Overload protection
Nucleus pulposus forceps	R302600BS	3.0mm	260mm	0°	Overload protection
Nucleus pulposus forceps	R352600BS	3.5mm	260mm	0°	Overload protection
Nucleus pulposus forceps	R402600BS	4.0mm	260mm	0°	Overload protection
Nucleus pulposus forceps	R252645BS	2.5mm	260mm	45°	Overload protection
Nucleus pulposus forceps	R302645BS	3.0mm	260mm	45°	Overload protection
Nucleus pulposus forceps	R252600CS	2.5mm	260mm	0°	Overload protection
Nucleus pulposus forceps	R302600CS	3.0mm	260mm	0°	Overload protection
Nucleus pulposus forceps	R352600CS	3.5mm	260mm	0°	Overload protection
Nucleus pulposus forceps	R252615CS	2.5mm	260mm	15°	Overload protection
Punch forceps	R252600BC	2.5mm	260mm	0°	Overload protection
Punch forceps	R302600BC	3.0mm	260mm	0°	Overload protection
Punch forceps	R402600BC	4.0mm	260mm	0°	Overload protection
Punch forceps	R252615BC	2.5mm	260mm	15°	Overload protection
Punch forceps	R302615BC	3.0mm	260mm	15°	Overload protection
Punch forceps	R252645BC	2.5mm	260mm	45°	Overload protection
Punch forceps	R352645BC	3.5mm	260mm	45°	Overload protection
Grasping forceps, pawl	R252600FB	2.5mm	260mm	0°	Overload protection
Grasping forceps, pawl	R302600FB	3.0mm	260mm	0°	Overload protection
Grasping forceps, pawl	R352600FB	3.5mm	260mm	0°	Overload protection
Grasping forceps, pawl	R302600FT	3.0mm	260mm	0°	Overload protection
Suture cutter forceps	R302600HS	3.0mm	260mm	0°	Overload protection
Flexible forceps	F252600ST	2.5mm	260mm	0°	Overload protection
Flexible forceps	F302600ST	3.0mm	260mm	0°	Overload protection
Flexible forceps	F252645BS	2.5mm	260mm	45°	Overload protection
Flexible forceps	F302645BS	3.0mm	260mm	45°	Overload protection
Flexible forceps	F252600BC	2.5mm	260mm	0°	Overload protection
Flexible forceps	F302600BC	3.0mm	260mm	0°	Overload protection
Flexible forceps	F252600CS	2.5mm	260mm	0°	Overload protection
Flexible forceps	F302600CS	3.0mm	260mm	0°	Overload protection



**SURGINEST
GLOBAL**
INSTRUMENTAL Y EQUIPOS MÉDICOS

Guide bar dilator



Dilator(4 grades for 6.3mm endoscope)

Product name	Model number	Outside diameter	Working length	Inside diameter	Functional characteristics
Guide bar (Dilator)	CY4631	7.5mm	175mm	6.5mm	With anti slip design
Guide bar (Dilator)	CY4632	6.4mm	190mm	5.2mm	With anti slip design
Guide bar (Dilator)	CY4633	5.0mm	212mm	3.4mm	With anti slip design
Guide bar (Dilator)	CY4634	3.2mm	235mm	1.5mm	With anti slip design

Dilator(4 grades for 6.9mm endoscope)

Product name	Model number	Outside diameter	Working length	Inside diameter	Functional characteristics
Guide bar (Dilator)	CY4691	8.3mm	180mm	7.3mm	With anti slip design
Guide bar (Dilator)	CY4692	7.0mm	200mm	5.2mm	With anti slip design
Guide bar (Dilator)	CY4693	5.0mm	220mm	3.4mm	With anti slip design
Guide bar (Dilator)	CY4694	3.2mm	240mm	1.6mm	With anti slip design



Dilator(2 grades for 6.3mm endoscope)

Product name	Model number	Outside diameter	Working length	Inside diameter	Functional characteristics
Guide bar (Dilator)	CY2631	6.3mm	210mm	3.1mm	With anti slip design
Guide bar (Dilator)	CY2632	3.0mm	230mm	1.5mm	With anti slip design

Dilator(2 grades for 6.9mm endoscope)

Product name	Model number	Outside diameter	Working length	Inside diameter	Functional characteristics
Guide bar (Dilator)	CY2691	7.0mm	210mm	3.1mm	With anti slip design
Guide bar (Dilator)	CY2692	3.0mm	230mm	1.5mm	With anti slip design



Dilator(3 grades for for cervical vertebra)

Product name	Model number	Outside diameter	Working length	Inside diameter	Functional characteristics
Guide bar (Dilator)	CY6952	6.9mm	150mm	5.2mm	With anti slip design
Guide bar (Dilator)	CY5034	5.0mm	160mm	3.4mm	With anti slip design
Guide bar (Dilator)	CY3210	3.2mm	170mm	1.0mm	With anti slip design

Guide bar dilator



Guide bar(Single channel)

Product name	Model number	Diameter	Working length	Functional characteristics
Guide bar(Single channel)	CY6523	6.5mm	230mm	With anti slip design
Guide bar(Single channel)	CY7023	7.0mm	230mm	With anti slip design
Guide bar(Single channel)	CY7523	7.5mm	230mm	With anti slip design
Guide bar(Single channel)	CY8023	8.0mm	230mm	With anti slip design
Guide bar(Single channel)	CY8523	8.5mm	230mm	With anti slip design
Guide bar(Single channel)	CY9523	9.5mm	230mm	With anti slip design



Guide bar(Double channel)

Product name	Model number	Diameter	Working length	Functional characteristics
Guide bar(Single channel)	CY6423D	6.4mm	230mm	With anti slip design
Guide bar(Single channel)	CY6923D	6.9mm	230mm	With anti slip design
Guide bar(Single channel)	CY7523D	7.5mm	230mm	With anti slip design



Eccentric rod

Product name	Model number	Diameter	Working length	Functional characteristics
Guide bar(Single channel)	CY6423E	6.4mm	230mm	With anti slip design
Guide bar(Single channel)	CY7023E	7.0mm	230mm	With anti slip design



Cannula



Cannula (Endoscopic working cannula)

Product name	Model number	Outside diameter	Inside diameter	Working length	Angle	Functional characteristics
Cannula	C75651653	7.5mm	6.5mm	165mm	30°	For 6.3mm endoscope
Cannula	C75651203	7.5mm	6.5mm	120mm	30°	For 6.3mm endoscope
Cannula	C75651654	7.5mm	6.5mm	165mm	45°	For 6.3mm endoscope
Cannula	C75651204	7.5mm	6.5mm	120mm	45°	For 6.3mm endoscope
Cannula	C75651657	7.5mm	6.5mm	165mm	70°	For 6.3mm endoscope
Cannula	C75651207	7.5mm	6.5mm	120mm	70°	For 6.3mm endoscope
Cannula	C7565165D	7.5mm	6.5mm	165mm	Duckbilled	For 6.3mm endoscope
Cannula	C7565120D	7.5mm	6.5mm	120mm	Duckbilled	For 6.3mm endoscope
Cannula	C82711653	8.2mm	7.1mm	165mm	30°	For 6.9mm endoscope
Cannula	C82711203	8.2mm	7.1mm	120mm	30°	For 6.9mm endoscope
Cannula	C82711654	8.2mm	7.1mm	165mm	45°	For 6.9mm endoscope
Cannula	C82711204	8.2mm	7.1mm	120mm	45°	For 6.9mm endoscope
Cannula	C82711657	8.2mm	7.1mm	165mm	70°	For 6.9mm endoscope
Cannula	C82711207	8.2mm	7.1mm	120mm	70°	For 6.9mm endoscope
Cannula	C8271165D	8.2mm	7.1mm	165mm	Duckbilled	For 6.9mm endoscope
Cannula	C8271120D	8.2mm	7.1mm	120mm	Duckbilled	For 6.9mm endoscope

Metal puncture needle



Product name	Model number	Specifications	Working length
Puncture needle	P16150	16G	150mm
Puncture needle	P16200	16G	200mm
Puncture needle	P16300	16G	300mm
Puncture needle	P18150	18G	150mm
Puncture needle	P18200	18G	200mm
Puncture needle	P18300	18G	300mm



Cannula



Cannula (Endoscopic working cannula with sealing cap)

Product name	Model number	Outside diameter	Inside diameter	Working length	Angle	Functional characteristics
Cannula	C75651653T	7.5mm	6.5mm	165mm	30°	For 6.3mm endoscope
Cannula	C75651203T	7.5mm	6.5mm	120mm	30°	For 6.3mm endoscope
Cannula	C75651654T	7.5mm	6.5mm	165mm	45°	For 6.3mm endoscope
Cannula	C75651204T	7.5mm	6.5mm	120mm	45°	For 6.3mm endoscope
Cannula	C75651657T	7.5mm	6.5mm	165mm	70°	For 6.3mm endoscope
Cannula	C75651207T	7.5mm	6.5mm	120mm	70°	For 6.3mm endoscope
Cannula	C7565165DT	7.5mm	6.5mm	165mm	Duckbilled	For 6.3mm endoscope
Cannula	C7565120DT	7.5mm	6.5mm	120mm	Duckbilled	For 6.3mm endoscope
Cannula	C82711653T	8.2mm	7.1mm	165mm	30°	For 6.9mm endoscope
Cannula	C82711203T	8.2mm	7.1mm	120mm	30°	For 6.9mm endoscope
Cannula	C82711654T	8.2mm	7.1mm	165mm	45°	For 6.9mm endoscope
Cannula	C82711204T	8.2mm	7.1mm	120mm	45°	For 6.9mm endoscope
Cannula	C82711657T	8.2mm	7.1mm	165mm	70°	For 6.9mm endoscope
Cannula	C82711207T	8.2mm	7.1mm	120mm	70°	For 6.9mm endoscope
Cannula	C8271165DT	8.2mm	7.1mm	165mm	Duckbilled	For 6.9mm endoscope
Cannula	C8271120DT	8.2mm	7.1mm	120mm	Duckbilled	For 6.9mm endoscope

Guide needle

Product name	Model number	Specifications	Working length
Guide needle	G08450	0.8mm	450mm
Guide needle	G10450	1.0mm	450mm
Guide needle	G15450	1.5mm	450mm



**SURGINEST
GLOBAL**
INSTRUMENTAL Y EQUIPOS MEDICOS


Trephine (For endoscope)



Product name	Model number	Diameter	Working length
Trephine (Integrated handle)	T3533	3.5mm	330mm
Trephine (Integrated handle)	T4033	4.0mm	330mm



Product name	Model number	Diameter	Working length
Trephine (Trilateral)	T3533T	3.5mm	330mm
Trephine (Trilateral)	T4033T	4.0mm	330mm
Trephine handle (Trilateral)	T000B		



Product name	Model number	Diameter	Working length
Trephine (Pentagonal)	T3533P	3.5mm	330mm
Trephine (Pentagonal)	T4033P	4.0mm	330mm
Trephine handle (Pentagonal)	T000C		

Safe trephine / Safe trephine cannula



Product name	Model number	Diameter	Working length
Trephine handle integrated	T5023	5.0mm	230mm
Trephine handle integrated	T6523	6.5mm	230mm
Trephine handle integrated	T7523	7.5mm	230mm
Trephine handle integrated	T8223	8.2mm	230mm
Trephine handle integrated	T8523	8.5mm	230mm
Trephine handle split(Trilateral)	T5021T	5.0mm	210mm
Trephine handle split(Trilateral)	T6521T	6.5mm	210mm
Trephine handle split(Trilateral)	T7521T	7.5mm	210mm
Trephine handle split(Trilateral)	T8221T	8.2mm	210mm
Trephine handle split(Trilateral)	T8521T	8.5mm	210mm
Trephine handle split(Pentagonal)	T5021P	5.0mm	210mm
Trephine handle split(Pentagonal)	T6521P	6.5mm	210mm
Trephine handle split(Pentagonal)	T7521P	7.5mm	210mm
Trephine handle split(Pentagonal)	T8221P	8.2mm	210mm
Trephine handle split(Pentagonal)	T8521P	8.5mm	210mm
Trephine cannula	T6017	6.0mm	170mm
Trephine cannula	T7517	7.5mm	170mm
Trephine cannula	T8517	8.5mm	170mm
Trephine cannula	T9017	9.0mm	170mm
Trephine cannula	T9517	9.5mm	170mm




Visual trephine(For 5.3mm endoscope)




Visual trephine



Handle



Product name	Model number	Outside Diameter	Inside diameter	Working length	Functional characteristics
Visual trephine	T531204	6.5mm	5.5mm	120mm	For 130mm endoscope
Visual trephine	T531704	6.5mm	5.5mm	170mm	For 181mm endoscope



Product name	Model number	Outside Diameter	Inside diameter	Working length	Functional characteristics
Straight trephine	T531202	6.5mm	5.5mm	120mm	For 130mm endoscope
Straight trephine	T531702	6.5mm	5.5mm	170mm	For 181mm endoscope
Trilateral trephine	T531203	6.5mm	5.5mm	120mm	For 130mm endoscope
Trilateral trephine	T531703	6.5mm	5.5mm	170mm	For 181mm endoscope



Product name	Model number	Outside Diameter	Inside diameter	Working length	Functional characteristics
Trephine cannula	C7767102D	7.7mm	6.7mm	102mm	For 120mm trephine
Trephine cannula	C7767153D	7.7mm	6.7mm	153mm	For 170mm trephine



Product name	Model number	Outside Diameter	Inside diameter	Working length	Functional characteristics
Endoscope cannula	C6555120D	6.5mm	5.5mm	120mm	For 130mm endoscope
Endoscope cannula	C6555170D	6.5mm	5.5mm	170mm	For 181mm endoscope


Visual trephine(For 5.9mm endoscope)




Visual trephine



Handle



Product name	Model number	Outside Diameter	Inside diameter	Working length	Functional characteristics
Visual trephine	T591204	7.1mm	6.1mm	120mm	For 130mm endoscope
Visual trephine	T591704	7.1mm	6.1mm	170mm	For 181mm endoscope



Product name	Model number	Outside Diameter	Inside diameter	Working length	Functional characteristics
Straight trephine	T591202	7.1mm	6.1mm	120mm	For 130mm endoscope
Straight trephine	T591702	7.1mm	6.1mm	170mm	For 181mm endoscope
Trilateral trephine	T591203	7.1mm	6.1mm	120mm	For 130mm endoscope
Trilateral trephine	T591703	7.1mm	6.1mm	170mm	For 181mm endoscope



Product name	Model number	Outside Diameter	Inside diameter	Working length	Functional characteristics
Trephine cannula	C8373102D	8.3mm	7.3mm	102mm	For 120mm trephine
Trephine cannula	C8373153D	8.3mm	7.3mm	153mm	For 170mm trephine



Product name	Model number	Outside Diameter	Inside diameter	Working length	Functional characteristics
Endoscope cannula	C7161120D	7.1mm	6.1mm	120mm	For 130mm endoscope
Endoscope cannula	C7161170D	7.1mm	6.1mm	170mm	For 181mm endoscope




Visual trephine(For 6.3mm endoscope)




Visual trephine



Handle



Product name	Model number	Outside Diameter	Inside diameter	Working length	Functional characteristics
Visual trephine	T631204	7.5mm	6.5mm	120mm	For 130mm endoscope
Visual trephine	T631704	7.5mm	6.5mm	170mm	For 181mm endoscope



Product name	Model number	Outside Diameter	Inside diameter	Working length	Functional characteristics
Straight trephine	T631202	7.5mm	6.5mm	120mm	For 130mm endoscope
Straight trephine	T631702	7.5mm	6.5mm	170mm	For 181mm endoscope
Trilateral trephine	T631203	7.5mm	6.5mm	120mm	For 130mm endoscope
Trilateral trephine	T631703	7.5mm	6.5mm	170mm	For 181mm endoscope



Product name	Model number	Outside Diameter	Inside diameter	Working length	Functional characteristics
Trephine cannula	C8777102D	8.7mm	7.7mm	102mm	For 120mm trephine
Trephine cannula	C8777153D	8.7mm	7.7mm	153mm	For 170mm trephine



Product name	Model number	Outside Diameter	Inside diameter	Working length	Functional characteristics
Endoscope cannula	C7565120D	7.5mm	6.5mm	120mm	For 130mm endoscope
Endoscope cannula	C7565170D	7.5mm	6.5mm	170mm	For 181mm endoscope




Visual trephine(For 6.9mm endoscope)




Visual trephine



Handle



Product name	Model number	Outside Diameter	Inside diameter	Working length	Functional characteristics
Visual trephine	T691204	8.1mm	7.1mm	120mm	For 130mm endoscope
Visual trephine	T691704	8.1mm	7.1mm	170mm	For 181mm endoscope



Product name	Model number	Outside Diameter	Inside diameter	Working length	Functional characteristics
Straight trephine	T691202	8.1mm	7.1mm	120mm	For 130mm endoscope
Straight trephine	T691702	8.1mm	7.1mm	170mm	For 181mm endoscope
Trilateral trephine	T691203	8.1mm	7.1mm	120mm	For 130mm endoscope
Trilateral trephine	T691703	8.1mm	7.1mm	170mm	For 181mm endoscope



Product name	Model number	Outside Diameter	Inside diameter	Working length	Functional characteristics
Trephine cannula	C9383102D	9.3mm	8.3mm	102mm	For 120mm trephine
Trephine cannula	C9383153D	9.3mm	8.3mm	153mm	For 170mm trephine



Product name	Model number	Outside Diameter	Inside diameter	Working length	Functional characteristics
Endoscope cannula	C8171120D	8.1mm	7.1mm	120mm	For 130mm endoscope
Endoscope cannula	C8171170D	8.1mm	7.1mm	170mm	For 181mm endoscope




Visual trephine(For 7.3mm endoscope)




Visual trephine



Handle




Product name	Model number	Outside Diameter	Inside diameter	Working length	Functional characteristics
Visual trephine	T731204	8.5mm	7.5mm	120mm	For 130mm endoscope
Visual trephine	T731704	8.5mm	7.5mm	170mm	For 181mm endoscope



Product name	Model number	Outside Diameter	Inside diameter	Working length	Functional characteristics
Straight trephine	T731202	8.5mm	7.5mm	120mm	For 130mm endoscope
Straight trephine	T731702	8.5mm	7.5mm	170mm	For 181mm endoscope
Trilateral trephine	T731203	8.5mm	7.5mm	120mm	For 130mm endoscope
Trilateral trephine	T731703	8.5mm	7.5mm	170mm	For 181mm endoscope




Product name	Model number	Outside Diameter	Inside diameter	Working length	Functional characteristics
Trephine cannula	C9787102D	9.7mm	8.7mm	102mm	For 120mm trephine
Trephine cannula	C9787153D	9.7mm	8.7mm	153mm	For 170mm trephine



Product name	Model number	Outside Diameter	Inside diameter	Working length	Functional characteristics
Endoscope cannula	C8575120D	8.5mm	7.5mm	120mm	For 130mm endoscope
Endoscope cannula	C8575170D	8.5mm	7.5mm	170mm	For 181mm endoscope

Safe bone drill



Product name	Model number	Diameter	Working length
Handle of bone drill	D00001		
Bone drill	D40200	4.0mm	200mm
Bone drill	D50200	5.0mm	200mm
Bone drill	D60200	6.0mm	200mm
Bone drill	D70200	7.0mm	200mm
Bone drill	D75200	7.5mm	200mm
Bone drill	D80200	8.0mm	200mm
Bone drill	D85200	8.5mm	200mm
Bone drill	D90200	9.0mm	200mm
TOM needle(Sharp head)	D21185S	2.1mm	185mm
TOM needle(Blunt head)	D21185B	2.1mm	185mm
TOM needle(Drill head)	D21185D	2.1mm	185mm



Safe bone drill



Visual reamer

A variety of visual reamers can safely and accurately handle intervertebral discs and endplates.

Product name	Model number	Specifications
Eccentric reamer	ER636910	For 6.3/6.9/10mm endoscope
Disc reamer	DR636910	For 6.3/6.9/10mm endoscope
Endplate reamer	ER636910S	For 6.3/6.9/10mm endoscope



Nerve hook/Nerve dissector/Osteotome/Nerve probe/Operating knife



Product name	Model number	Working length	Diameter
Nerve hook	NH3323	330mm	2.5/3.0mm
Nerve hook	NH2623	260mm	2.5/3.0mm
Nerve dissector	ND3323	330mm	2.5/3.0mm
Nerve dissector	ND2623	260mm	2.5/3.0mm
Osteotome	OS3323	330mm	2.5/3.0/3.5mm
Osteotome	OS2623	260mm	2.5/3.0/3.5mm
Nerve probe	NP3330	330mm	3.0mm
Nerve probe	NP2630	260mm	3.0mm
Operating knife	OK3330	330mm	3.0mm
Operating knife	OK2630	260mm	3.0mm



Product name	Model number	Working length	Diameter
Telescopic probe	TP3330	330mm	3.0mm

Bone hammer



BH0001 Bone hammer



**SURGINEST
GLOBAL**
INSTRUMENTAL Y EQUIPOS MÉDICOS

Large channel stenosis endoscope




Product name	Model number	Diameter	Working length	Angle	Instrument channel	Functional
Interlaminar endoscope	CY1012515	10mm	125mm	15°	7.1mm	HD autoclave
Interlaminar endoscope	CY1012530	10mm	125mm	30°	7.1mm	HD autoclave

Large channel stenosis rongeur forceps



Product name	Model number	Diameter	Working length	Angle	Optional
Handle of forceps	CR00001				
Rongeur forceps	CR5226040	5.5/2.5mm	260mm	40°	Black-gold coated rubber 360° quick installation
Rongeur forceps	CR5326040	5.5/3.5mm	260mm	40°	
Rongeur forceps	CR5226090	5.5/2.5mm	260mm	90°	
Rongeur forceps	CR5326090	5.5/3.5mm	260mm	90°	

Large channel stenosis nucleus pulposus forceps



Product name	Model number	Diameter	Working length	Angle
Punch forceps	R402600BC	4.0mm	260mm	0°
Punch forceps	R402615BC	4.0mm	260mm	15°
Punch forceps	R352640BC	3.5mm	260mm	40°
Nucleus pulposus forceps	R352600BS	3.5mm	260mm	0°
Nucleus pulposus forceps	R402600BS	4.0mm	260mm	0°
Nucleus pulposus forceps	R302645BS	3.0mm	260mm	45°
Nucleus pulposus forceps	R352645BS	3.5mm	260mm	45°
Nucleus pulposus forceps	R352600CS	3.5mm	260mm	0°
Nucleus pulposus forceps	R402600CS	4.0mm	260mm	0°
Grasping forceps, pawl	R352600FB	3.5mm	260mm	0°
Grasping forceps	R402600FT	4.0mm	260mm	0°
Flexible forceps	F252600ST	2.5mm	260mm	0°
Flexible forceps	F302600ST	3.0mm	260mm	0°
Flexible forceps	F252600BC	2.5mm	260mm	0°
Flexible forceps	F302600BC	3.0mm	260mm	0°

Delta cannula(For 10mm stenosis endoscope)




Product name	Model number	Diameter	Working length	Angle	Functional characteristics
Endoscope cannula	C1412D1	14mm	120mm	Duck mouth	Perspective development with thread and sealing cap to continuously maintain the water pressure under the endoscope.
Endoscope cannula	C1412F1	14mm	120mm	Flat mouth	
Endoscope cannula	C1412D2	14mm	120mm	Duck mouth	With screw thread and sealing cap to continuously maintain the water pressure under the endoscope.
Endoscope cannula	C1412F2	14mm	120mm	Flat mouth	
Endoscope cannula	C11212	11.2mm	120mm	Inclined mouth	Non threaded but with sealing cap to continuously maintain the water pressure under the endoscope.

Delta cannula(For 10mm stenosis endoscope)



Product name	Model number	Diameter	Working length	Angle	Functional
Endoscope cannula	C6195F	6mm	195mm	Flat mouth	Inside endoscope
Endoscope cannula	C6231D	6mm	231mm	Duck mouth	Inside endoscope

Delta trephine(For 10mm stenosis endoscope)



Product name	Model number	Diameter	Inside diameter	Working length	Angle	Functional
Visible trephine	T111015F	11.2mm	10.2mm	115mm	Flat mouth	Outside endoscope
Trephine cannula	T121195D	12.5mm	11.5mm	95mm	Flat mouth	Outside endoscope
Visible trephine	T655526F	6.5mm	5.5mm	260mm	Flat mouth	Inside endoscope

Large channel interlaminar endoscope




Product name	Model number	Diameter	Working length	Angle	Instrument channel	Functional
Interlaminar endoscope	CY8412015	8.4mm	120mm	15°	5.7mm	HD autoclave
Interlaminar endoscope	CY8412030	8.4mm	120mm	30°	5.7mm	HD autoclave

Large channel interlaminar rongeur forceps



Product name	Model number	Diameter	Working length	Angle	Optional
Handle of forceps	CR00001				Black-gold coated rubber 360° quick installation
Rongeur forceps	CR4226040	4.8/2.5mm	260mm	40°	
Rongeur forceps	CR4326040	4.8/3.5mm	260mm	40°	
Rongeur forceps	CR4226090	4.8/2.5mm	260mm	90°	
Rongeur forceps	CR4326090	4.8/3.5mm	260mm	90°	

Large channel interlaminar nucleus pulposus forceps



Product name	Model number	Diameter	Working length	Angle
Punch forceps	R402600BC	4.0mm	260mm	0°
Punch forceps	R402615BC	4.0mm	260mm	15°
Punch forceps	R352640BC	3.5mm	260mm	40°
Nucleus pulposus forceps	R352600BS	3.5mm	260mm	0°
Nucleus pulposus forceps	R402600BS	4.0mm	260mm	0°
Nucleus pulposus forceps	R302645BS	3.0mm	260mm	45°
Nucleus pulposus forceps	R352645BS	3.5mm	260mm	45°
Nucleus pulposus forceps	R352600CS	3.5mm	260mm	0°
Nucleus pulposus forceps	R402600CS	4.0mm	260mm	0°
Grasping forceps, pawl	R352600FB	3.5mm	260mm	0°
Grasping forceps	R402600FT	4.0mm	260mm	0°
Flexible forceps	F252600ST	2.5mm	260mm	0°
Flexible forceps	F302600ST	3.0mm	260mm	0°
Flexible forceps	F252600BC	2.5mm	260mm	0°
Flexible forceps	F302600BC	3.0mm	260mm	0°

Delta cannula(For 8.4mm interlaminar endoscope)




Product name	Model number	Diameter	Working length	Angle	Functional characteristics
Endoscope cannula	C1211D1	12mm	115mm	Duck mouth	Perspective development with thread and sealing cap to continuously maintain the water pressure under the endoscope.
Endoscope cannula	C1211F1	12mm	115mm	Flat mouth	
Endoscope cannula	C1211D2	12mm	115mm	Duck mouth	With screw thread and sealing cap to continuously maintain the water pressure under the endoscope.
Endoscope cannula	C1211F2	12mm	115mm	Flat mouth	
Endoscope cannula	C96115	9.6mm	115mm	Inclined mouth	Non threaded but with sealing cap to continuously maintain the water pressure under the endoscope.

Delta cannula(For 8.4mm interlaminar endoscope)



Product name	Model number	Diameter	Working length	Angle	Functional
Endoscope cannula	C5190F	5mm	190mm	Flat mouth	Inside endoscope
Endoscope cannula	C5226D	5mm	226mm	Duck mouth	Inside endoscope

Delta trephine(For 8.4mm interlaminar endoscope)



Product name	Model number	Diameter	Inside diameter	Working length	Angle	Functional
Visible trephine	T968611F	9.6mm	8.6mm	110mm	Flat mouth	Outside endoscope
Trephine cannula	T109890F	10.8mm	9.8mm	90mm	Flat mouth	Outside endoscope
Visible trephine	T514126F	5.1mm	4.1mm	260mm	Flat mouth	Inside endoscope



**SURGINEST
GLOBAL**
INSTRUMENTAL Y EQUIPOS MEDICOS

Delta guide needle

Product name	Model number	Diameter	Working length
Guide needle	G15250	1.5mm	250mm

Delta dilator




Product name	Model number	Diameter	Inside diameter	Working length
Dilator(1 grade)	D031723	3.0mm	1.7mm	230mm
Dilator(2 grade)	D053221	5.0mm	3.2mm	210mm
Dilator(3 grade)	D085219	8.0mm	5.2mm	190mm
Dilator(4 grade)	D108217	10.0mm	8.2mm	170mm
Dilator(5 grade)	D111015	11.3mm	10.2mm	150mm
Dilator(single one)	D101218	10.0mm	1.2mm	180mm

Operating knife




	Product name	Model number	Diameter	Working length
	Nerve hook	NH3026	3.0mm	260mm
	Nerve dissector	ND3026	3.0mm	260mm
	Curette	CU4026	4.0mm	260mm
	Osteotome	OS4026	4.0mm	260mm

Multi-function surgical serpentine pneumatic arm



Performance	
Bearing capacity	1.5kg
Power source meal	Central air supply, compressed gas cylinder, compressor
Air pressure	0.4-0.8MPa air, nitrogen, carbon dioxide

UBE (Unilateral biportal endoscopy)



Product name	Model number	Diameter	Working length	Angle
UBE endoscope	4170.00U	4.0mm	170mm	0°
UBE endoscope	4170.30U	4.0mm	170mm	30°
UBE endoscope	4170.70U	4.0mm	170mm	70°


Sheath & Obturator



Product name	Model number	Diameter	Working length	Angle
UBE sheath	U3001	5.5mm	135mm	0°/30°/70°

Forceps

Rongeur forceps/Punch forceps/Nucleus pulposus forceps



Model number	Outside diameter	Working length	Angle
CR0001	Laminar rongeur handle	360° rotating	Detachable
CR2023040	2.0mm	230mm	40°
CR3023040	3.0mm	230mm	40°
CR4023040	4.0mm	230mm	40°
CR3023040U	3.0mm	230mm	40° (Upward)



Product name	Model number	Diameter	Working length	Angle
Punch forceps	R4018000BC	4.0mm	180mm	0°
Punch forceps	R4018045BC	4.0mm	180mm	45°
Nucleus pulposus forceps	R4018000WG	4.0mm	180mm	0°
Nucleus pulposus forceps	R4018030WG	4.0mm	180mm	30°
Nucleus pulposus forceps	R4018000ST	4.0mm	180mm	0°
Nucleus pulposus forceps	R3018045BS	3.0mm	180mm	45°



Product name	Model number	Specifications
Adjustable reamer	R0814A	Minimun 8mm-Maximun 14mm



SURGINEST
GLOBAL
INSTRUMENTAL Y EQUIPOS MÉDICOS

UBE instrument



Product name	Model number	Diameter	Working length
Guide needle	G4210	4.0mm	210mm



Product name	Model number	Specifications
Muscle detacher	CY1001	Opener



Product name	Model number	Specifications
Dilator	CY0005	5mm
Dilator	CY0007	7mm
Dilator	CY0009	9mm
Dilator	CY0011	11mm
Dilator	CY0013	13mm
Dilator	CY0015	15mm



UBE instrument



Product name	Model number	Specifications
Cage guider	CY1002	



Product name	Model number	Specifications
Root retractor	CY4108	4mm x 108mm
Root retractor	CY1108	10mm x 108mm



Product name	Model number	Specifications
Reamer	CY8180	8mm x 180mm




Product name	Model number	Specifications
Half casing tube	CY1060	10mm x 60mm



UBE instrument



	Product name	Model number	Specifications
	Biportal operate knife	CY2005B	5mm left bend
	Biportal operate knife	CY2005R	5mm left bend
	Biportal operate knife	CY2005S	5mm left bend
	Biportal operate knife	CY2005LB	5mm left bend
	Biportal operate knife	CY2005T	5mm left bend



	Product name	Model number	Specifications
	Biportal curette	CY2003RP	3mm Reverse port
	Biportal curette	CY2003LB	3mm Large bend
	Biportal curette	CY2003FB	3mm Forward bend
	Biportal curette	CY2005FB	5mm Forward bend
	Biportal curette	CY2006FB	6mm Forward bend



UBE instrument



Product name	Model number	Specifications
Double ended tetacher	CY0515	3mm x (5° ~ 15°)
Double ended tetacher	CY2535	3mm x (25° ~ 35°)



Product name	Model number	Specifications
Nerve hook	CY2004	4mm



Product name	Model number	Specifications
Bone grafting funnel	CY7565	7.5mm x 165mm



Product name	Model number	Specifications
Bone grafting push rod	CY2014	



Product name	Model number	Specifications
Bone hammer	BH0002	



Disposable plasma & radiofrequency surgical electrode




Specifications

Surgical electrode(For transforaminal & interlaminar)

Product name	Model number	Probe Diameter	Working length	Plug Type
Surgical electrode	E20400A	2.0mm	400mm	Type A-Rectangle
Surgical electrode	E20340A	2.0mm	340mm	Type A-Rectangle
Surgical electrode	E20240A	2.0mm	240mm	Type A-Rectangle
Surgical electrode	E20150A	2.0mm	150mm	Type A-Rectangle
Surgical electrode	E20400B	2.0mm	400mm	Type B-Round
Surgical electrode	E20340B	2.0mm	340mm	Type B-Round
Surgical electrode	E20240B	2.0mm	240mm	Type B-Round
Surgical electrode	E20150B	2.0mm	150mm	Type B-Round
Surgical electrode	E20400C	2.0mm	400mm	Type C-Split
Surgical electrode	E20340C	2.0mm	340mm	Type C-Split
Surgical electrode	E20240C	2.0mm	240mm	Type C-Split
Surgical electrode	E20150C	2.0mm	150mm	Type C-Split




Disposable plasma surgical electrode



A

Without suction tube



B

With suction tube

Surgical electrode(For sports medicine)

Product name	Model number	Diameter	Working length	Angle	Suction type	
Surgical electrode	E2907030A	2.9mm	70mm	30°	A	
Surgical electrode	E2909030A	2.9mm	90mm	30°	A	
Surgical electrode	E2912030A	2.9mm	120mm	30°	A	
Surgical electrode	E3412030A	3.4mm	120mm	30°	A	
Surgical electrode	E3414030A	3.4mm	140mm	30°	A	
Surgical electrode	E3416030A	3.4mm	160mm	30°	A	
Surgical electrode	E3812090B	3.8mm	120mm	90°	B	
Surgical electrode	E3814090B	3.8mm	140mm	90°	B	
Surgical electrode	E3816090B	3.8mm	160mm	90°	B	
Surgical electrode	E4012090A	4.0mm	120mm	90°	A	
Surgical electrode	E4014090A	4.0mm	140mm	90°	A	
Surgical electrode	E4016090A	4.0mm	160mm	90°	A	
Surgical electrode	E3912040B	3.9mm	120mm	40°	B	
Surgical electrode	E3915040B	3.9mm	150mm	40°	B	
Surgical electrode	E3918040B	3.9mm	180mm	40°	B	
Surgical electrode	E1208010A	1.2mm	80mm	10°	A	
Surgical electrode	E1211010A	1.2mm	110mm	10°	A	
Surgical electrode	E1214010A	1.2mm	140mm	10°	A	
Surgical electrode	E1217010A	1.2mm	170mm	10°	A	
Surgical electrode	E2912030A	2.9mm	120mm	30°	A	
Surgical electrode	E2915030A	2.9mm	150mm	30°	A	
Surgical electrode	E2918030A	2.9mm	180mm	30°	A	




**SURGINEST
GLOBAL**
INSTRUMENTAL Y EQUIPOS MEDICOS

Instruments sterilization box




S0001 Instruments sterilization box
560x260x130mm



Endoscope sterilization box



S0002 Endoscope sterilization box
430x145x75mm



Transforaminal / Interlaminar endoscope



Model number	Outerside diameter	Working length	Angle	Field angle	Working channel	Notes
CY3309000	3.3mm	90mm	0°	60°	2.1mm	HD autoclave
CY5313030	5.3mm	130mm	30°	80°	2.7mm	HD autoclave
CY5318130	5.3mm	181mm	30°	80°	2.7mm	HD autoclave
CY5913630	5.9mm	136mm	30°	80°	3.3mm	HD autoclave
CY5918130	5.9mm	181mm	30°	80°	3.3mm	HD autoclave
CY6313015	6.3mm	130mm	15°	80°	3.7mm	HD autoclave
★ CY6313030	6.3mm	130mm	30°	80°	3.7mm	HD autoclave
CY6318100	6.3mm	181mm	0°	80°	3.7mm	HD autoclave
CY6318115	6.3mm	181mm	15°	80°	3.7mm	HD autoclave
★ CY6318130	6.3mm	181mm	30°	80°	3.7mm	HD autoclave
CY6319330	6.3mm	193mm	30°	80°	3.7mm	HD autoclave
CY6913015	6.9mm	130mm	15°	80°	4.3mm	HD autoclave
★ CY6913030	6.9mm	130mm	30°	80°	4.3mm	HD autoclave
CY6918100	6.9mm	181mm	0°	80°	4.3mm	HD autoclave
CY6918115	6.9mm	181mm	15°	80°	4.3mm	HD autoclave
★ CY6918130	6.9mm	181mm	30°	80°	4.3mm	HD autoclave
CY6919330	6.9mm	193mm	30°	80°	4.3mm	HD autoclave
CY7313030	7.3mm	130mm	30°	80°	4.7mm	HD autoclave
CY7318130	7.3mm	181mm	30°	80°	4.7mm	HD autoclave
★ CY8412015	8.4mm	120mm	15°	80°	5.7mm	HD autoclave
CY1013630	10mm	136mm	30°	80°	7.1mm	HD autoclave
CY1013615	10mm	136mm	15°	80°	7.1mm	HD autoclave
★ CY1012515	10mm	125mm	15°	80°	7.3mm	HD autoclave

CONTÁCTANOS



+56 9 9851 8862



SURGINEST_GLOBAL



SURGINEST GLOBAL



SURGINEST
GLOBAL
INSTRUMENTAL Y EQUIPOS MÉDICOS